



APE SOFTWARE
TOOLS FOR MANAGEMENT SYSTEMS

Calibration Control Manual

Calibration Management Software

Table of Contents

Sample Calibration Databases.....	15
Explore Calibration Control with Sample Data	15
Selecting the MS Access Sample Database on First Start.....	15
Selecting Sample Database After Program Start	16
Reset Configuration	17
Sample SQL Database	18
Remove Sample Data	18
Remove Sample Data to Begin Entering Your Own Records	18
Double Click (Automatic Method).....	18
Move Application Database	20
Move Application Database to a New Location.....	20
STEP 1: Close Application.....	20
STEP 2: Find Database File	20
STEP 3: Move Database File.....	20
STEP 4: Point Application to New Location	20
Locating Database File (MS Access).....	22
Find the Database and Configuration Files	22
Method 1 - Version 5.5 and Higher	22
Method 2 - Version 8.0 and Higher	22
Method 3 - General Config File	23
Data Imported To MS Access	23
How to Deploy Your Imported Data	23
Install Current Version.....	23
Deploy Database (Scenario A).....	23
Remove Old Database (Scenario B).....	24
Additional Help.....	25
Creating a Multi-User Environment.....	25
Setting up Multiple Users.....	25
Move Your Database (MS Access Database Only)	26
Set Up a Common Files Folder (Access & SQL)	26
Add Additional Users	26
Add Ape Database to SQL Server	27
Install Calibration Control on a SQL Server.....	27
Restore Back-Up Database	27
Restore Failed	29
User Security	30
Connect to SQL Server Database.....	31
Configure Calibration Control to Connect to a SQL Server	31

Start Calibration Control	31
Restarting the Process	33
Connection Troubleshooting	33
Calibration Control Folders	33
Settings Folder	35
Grid Features	35
Features Common to All Data Grids	35
Highlight Colors of the Equipment Grid	39
Understanding the Highlighted Colors of Equipment Records	39
Showing/Hiding Colors	39
Default Colors	40
Highlight Past Due	41
Highlight Standards	42
Highlight Received Out of Tolerance	43
Received for Calibration	44
Filtering Records Using the Filter Row	44
How to Use the Filter Row in the Data Grids	44
Display Filter Row	44
Create New Filter	45
Clear Filter(s)	45
Filter Methods	46
Custom Filters	46
Hiding Equipment Records by Status	47
Use Equipment Status to Manage Visibility	47
Status Codes Visibility	47
Show/Hide All Hidden Records	48
Manage Out-Of-Tolerance (OOT) Equipment	48
Manage and Investigate OOT Equipment	48
Calibration Event	48
Red Highlighted Records	49
Calibration History	50
Calibration Event Quick Reference and Filter	50
Top N Records	50
Filtering and Sorting	51
Context Menu	51
Equipment Dialog	52
Edit Equipment Records	52
Equipment Dialog Fields	52
Calibrations Tab of Equipment Record	53
Maintenance	54

Attachments Documents Notes Parts Tab.....	55
Misc (Miscellaneous) Tab	55
Other Tabs in the Equipment Dialog	56
Calibration Events.....	57
Create a Simple Calibration Event	57
Calibration Event Fields	58
Attachments Documents Notes Parts Tab	58
Standards Misc. Tab	59
Calibration Due for Certificate.....	60
Measurement Group Fields	61
Standards Used in Measurement Group.....	63
Additional Tabs in the Edit Calibration Event Dialog	64
Calibration Frequencies.....	64
How do Calibration Frequencies Work?	64
Calibration Fields	65
“Month of” and “Week of” Frequencies	67
Referencing Calibration Standards.....	68
How to Link Standards to Calibrated Equipment	68
ID a Calibration Standard.....	68
Grid Highlight.....	69
Link Calibration Standards	69
Available Standards	70
Results.....	70
Calibration Templates.....	70
How to Use Calibration Templates	70
Default Calibration Templates.....	70
Create User-Defined Calibration Templates	71
Calibration Template Structure and Hierarchy.....	72
Link Templates (Hierarchy).....	72
Attachments.....	73
Add Attachments to Equipment and Calibration Records	73
Add File, Link to File, Folder or URL	74
Edit Attachment	74
Attachments Grid.....	75
Image in Equipment Dialog.....	76
Companies	77
Manage Companies within Calibration Control	77
Company Dialog	78
Email, Phones, People, and Types Dialog.....	78
Notes and Attachments Tab	79

Other Tabs	80
Company Types.....	81
Edit Company Type Codes	81
Company Types Dialog.....	82
Link to Company Records.....	82
Other Tabs	83
Equipment Systems	83
Groups of Equipment Records.....	83
View Systems Grid.....	83
Equipment System Dialog.....	84
Equipment Types.....	90
Classify Equipment and Link to Calibration Templates	90
Equipment Types Grid	90
Equipment Types Dialog	91
Other Tabs	92
Size and Range in Calibration Control.....	92
Edit the Size and Range Combo Values of Equipment	92
Menu.....	92
Dialog.....	93
Model Numbers	94
Keep Track of Equipment Model Numbers.....	94
Model Number Dialog.....	94
Uncertainty Tab.....	95
Overrides Tab.....	95
Custom Tab	96
Meta Tab.....	96
Departments	96
Edit Department Dialog	97
Other Tabs	98
Part Categories.....	98
Create Part Category Codes to Assign to Parts	98
Other Tabs	99
Assign a Part Category to Part Number Record	99
Part Numbers	99
Manage Part Numbers Related to Equipment.....	99
Edit Part Number Dialog	100
Part Number Sources.....	100
Equipment & Calibrations.....	101
Work Orders.....	102
Other Tabs	102

People	102
Person Dialog	103
Update Equipment Locations.....	104
Grids Dialog	104
Other Tabs	105
Documents Dialog	105
Procedures Are Now Known as Documents	105
Edit Document Dialog.....	106
Link Documents to Equipment Records	107
Email Addresses	108
Email Address Dialog.....	109
Companies Tab.....	109
People Tab	110
Other Tabs	110
Phone Numbers.....	110
Phone Number Dialog.....	111
Companies Tab.....	111
People Tab	112
Other Tabs	112
Notes	112
Notes Grid	113
Notes Dialog	113
Other Tabs	115
Application Codes	115
Codes List.....	116
Adding a New Code	116
Custom Tab	117
Meta Tab.....	118
Program Options.....	118
Modify Application Settings	118
General Tab.....	118
Advanced Tab.....	120
Options - Calibration	122
Options - Check Out	123
Options - Email	124
Options - Folder & Path	125
Change Field Names.....	126
Change the Text of Field Labels	126
Edit Form Labels.....	126
Edit Fields and Labels.....	128

Reset Changes	129
Custom Fields	129
Options - Reference Records	130
Add Missing	131
Refresh All	132
Reports & Labels	132
Default Reports	132
Default Labels Tab	133
Options - Security	134
Feature Visibility Options	135
Make Features Invisible	135
Minimize Ribbon	138
Options - Spell Checker.....	138
Measurements in Metric or Imperial.....	139
Display Metric or Imperial as Default Measurement Units.....	139
Asset Transfer Dialog.....	140
Use the Asset Transfer Dialog	140
Asset Issue & Receive Dialogs	141
Asset Issue Dialog	141
3-Step Process to Issue Equipment.....	141
Asset Receive Dialog.....	142
3- Step Process to Receive Equipment.....	142
Status Change Dialog	143
Quickly Change Status of Multiple Equipment	143
Status Change in Ribbon	143
Using the Status Change Dialog	143
Equipment Check In/Out.....	145
Opening the Check Out Dialog	145
Check Out.....	146
Check In.....	147
Check Out by Equipment System Name.....	148
Check Out with Barcode Scanner	148
Using Scanner for Barcoded Equipment ID	148
Barcode Recognition Order	149
Equipment Use Count or Use Days.....	149
Default Increments in Check Out.....	149
Verify Check Out Status.....	151
Check Out Log.....	152
Recording and Researching Jobs	153
Record Equipment Used on Jobs for Traceability.....	153

Browse Jobs.....	153
Auto Notify Utility.....	154
Send Emails and Desktop Notifications Automatically	154
User Authentication	154
Auto Notify Grid	154
Auto Notification Dialog.....	155
SQL tab.....	155
Body Template tab	156
Email Activity Log.....	157
Troubleshooting.....	157
Terminal Mode.....	158
Minor Tasks in Ape Terminal Mode	158
Limited Use.....	158
Further Limitations.....	159
Projects.....	159
Projects Dialog.....	160
Project Dialog Fields	160
Equipment/People/Parts/Work Orders	160
Notes/Attachments/Documents.....	161
Other Tabs	161
Work Orders	162
Work Orders Dialog	162
Submit Tab	162
Service Tab.....	163
Documents/Notes/Parts/Attachments Tab	163
Other Tabs	164
Custom Barcodes	164
Use Unique Barcodes Already on Equipment	164
Custom Barcode Field	164
Barcode Scan Dialog	165
Date Masks in Calibration Control	165
Number Masks in Calibration Control	167
Number Masks	168
Reserved Characters	169
Leading Zeros Mask	169
Calibration Measurement Group Mask.....	169
Troubleshooting Number Masks	170
Procedural Steps and Groups	170
Step-by-step Procedures for Worksheets.....	170
Concept	171

Procedural Groups.....	171
Procedural Steps	171
Add Procedural Steps to a Procedural Group.....	172
Linking Procedural Groups to Equipment	173
Set the Cal Worksheet with Procedural Steps as the Global Default Report	175
Override Instructions at Equipment Level	176
Required Fields.....	176
Make Procedural Groups feature Visible	177
Create a New Procedural Group.....	177
Custom Highlight Colors	179
Get Ape Support	179
Publish a Custom Website.....	180
Publish Website of Equipment Grouped by Departments	180
One-Click Website Publishing (Create).....	180
Auto Publication with Auto Notify Utility.....	180
Change Location Where Site is Published	182
Measurement Uncertainty Studies.....	182
Perform Uncertainty Analysis of a Measurement.....	182
Uncertainty Study Dialog	183
Uncertainty Detail Dialog.....	184
Learning about Measurement Uncertainty	185
References	185
THUM Temperature and Humidity Device	186
Automatically Record Temperature and Humidity.....	186
Setup	186
Troubleshooting.....	187
Best Label Printers for Ape Software.....	188
Which label printers will work with Ape's Software?	188
Brother P-touch Label Printers	188
TZe Label Tape.....	188
Extra Adhesive Label Tape.....	188
Free Label Editing Software and Drivers	189
Free Brother Printer Drivers and Label Editing Software.....	189
Printer Drivers.....	189
P-touch Editor - Label Editing Software	189
Labels Quick Start Instructions	189
Sample Labels	189
Important Points	189
Access All the Label Fields.....	190
Understanding Labels Fields in Calibration Control.....	190

Sample Due Cal Labels.....	191
Print Future Calibration Labels	193
Print Calibration Labels to Apply at a Future Date.....	193
Printing Future Labels	193
Sample Asset Labels.....	194
Calibration Not Required Labels	196
Chain Printing Labels	198
How to Chain Print Labels to save Label Tape	198
Individually Printed Labels	198
Chain Printed Labels.....	199
Labels & Reports Overrides.....	200
Override Global Defaults for Reports and Labels	200
Override Settings.....	200
Label Fields	201
Report Fields	201
Print Reports.....	202
Printing the Standard Reports is Easy.....	202
Calibration Due Report	203
Quick Access to the Due Cal Report.....	203
Due Cal: Print Preview.....	203
Calibration Worksheets.....	204
Record Calibration Data While Away from a Computer.....	204
Create a Custom Report.....	207
Create Custom Reports Using Report Designer.....	207
Create a Custom Report	207
Define Report Data Source	208
Connection String - Choosing the Provider	209
Connection String - Setting the Connection	210
Writing a SQL Select Statement	210
Adding Fields, Labels, and Report Info	211
Menu Settings.....	211
Preview and save the Report Layout	211
Report Parameters	211
Pass Data to SQL Statements at Run Time	211
Combobox Parameter Type.....	213
Report Settings.....	213
Find and Change Settings of Custom Reports	213
Open Settings	214
Page Setup	214
Printer Settings	214

Styles	215
Global	216
Update Default Reports	216
Quickly Update Default Report and Label Templates.....	216
Manually Update Default Reports Using a 2-Step Process	217
Steps to Update Default Reports	217
SQL SELECT Statement.....	217
Useful with Calibration Management Software Reports.....	217
A Primer on Tables.....	218
The Equipment Master View.....	218
The SQL SELECT Statement.....	219
Dashboard Charts.....	220
Edit Dashboard Charts.....	220
Displaying Charts.....	224
Use Pie, Column, and Bar Charts to Visualize Your Data	224
Chart Data Tables	224
Saving and Printing Charts	225
3D Charts.....	226
‘Min Others Pie Slice’ (for pie charts only)	227
Due Cal Calendars.....	228
Visualize Due Cal Equipment in Calendars	228
Printing Calendars	230
Security Methods - How They Work	231
Database Security	231
Activating User-Authentication (Sign-In Mode)	231
What an Administrator Can Do	231
Administrators and Passwords	231
Pre-Defined Roles.....	231
Hierarchy of Roles	231
Permissions	232
Change Log.....	232
Users.....	233
Creating a User.....	233
User Tab	233
Advanced tab	234
Notes Tab	235
Sign-In Mode (Enabling Authentication)	235
Set up and Turn on User Authentication	235
Creating an Admin User.....	235
Turning On User Authentication.....	235

Password Security	237
Using the Password Security Dialog	237
Maximum Settings	238
Minimum Characters	238
Check Settings.....	238
Reset Login by Admin.....	239
Database Password Encryption	239
Encrypting an MS Access Database with a Password	239
Permissions	240
Change the Roles Required for Any Permission	240
Reset Options	241
Change Log.....	241
Audit Changes Made in Calibration Control	241
Columns of the Change Log Grid.....	242
FDA 21 CFR Part 11 Compliance	242
FDA Software Compliance	242
Authentication at User Level	242
File Integrity	243
Change Management (Audit Trails)	243
Electronic Signatures	243
Password Strength.....	243
Inactivity Settings.....	243
Sessions.....	243
See Who is Using Calibration Control	243
Sessions Grid & Dialog	243
View All or Current Sessions	244
Machines	245
Control Computer Access in the Machines Grid	245
Licensing and Authentication	247
Machine Authentication (default)	247
User Authentication	247
Simple Troubleshooting	248
Fix Most Setup or Configuration Problems.....	248
Problem Steps Recorder.....	250
Crash Reports.....	250
Tell Ape About This Crash	250
Description	250
Followup.....	250
Find Lost Records in Equipment Grid	251
Troubleshooting Missing Records	251

Clearing the Filter Row	251
Showing Hidden Records	252
Lost or Invalid Product Key / Lost All Records	253
Troubleshoot an Invalid Product Key or Lost Records	253
How It Happens	253
How to Fix It.....	253
Related Help Topics.....	253
Troubleshoot Calibration Control Not Working.....	254
Why Did Calibration Control Stop Working?	254
Missing Database Drivers	254
Missing .NET Drivers	254
Missing Read/Write Privileges	254
Cannot Connect to an Unstable Database.....	254
Server is Waiting for a Response.....	254
Get Ape Support	254
Troubleshoot Database Location	255
Why Am I Seeing the "Choose Database Type" Dialog?	255
Last Known Location of Database	255
Troubleshooting Label Printing	256
Label Print Failed	256
Solution 1	257
Solution 2.....	257
Get Ape Support	257
No Label File Error.....	257
Label File Could Not Be Loaded	257
Label Printer Drivers	257
P-touch Editor Software	257
Troubleshoot SQL Server Connection.....	258
For Use with Ape Software Databases.....	258
Universal Data Link (UDL) File.....	258
Create a UDL File	258
Data Link Properties	259
SMTP Test Procedure.....	262
Introduction	262
Procedure.....	262
Verify Domain:	262
Adding TELNET:	263
Connect to SMTP server:.....	264
UDL File for Database Connection.....	265
Universal Data Link for Database Connection	265

Create Universal Data Link (UDL) File	265
Create Shortcut.....	268
Copy Shortcut to Clients	269
Startup INI	269
Single Connection Configuration for All Users	269
Before Using the INI Utility.....	269
Start the INI Utility.....	270
Startup INI Editor	270
Updating Access DB	271
Problem Description	271
Automatic Backups	271
Local Update	271
Current Access Database Drivers	272
Ape Software Free Update Service	272
Update a SQL Server Database.....	272
Update an Ape Software SQL Server Database	272
Before Updating	272
Updating Database	273
Get Ape Support	276
Backup and Restore SQL Server Database.....	276
Create a Full Database Backup	276
Restore a SQL Server Database Backup.....	279
Downgrading a Database	281
Return a Database to the Previous Version	281
Restore an Access Database.....	281
Restore a SQL Server Database.....	281


Sample Calibration Databases

Explore Calibration Control with Sample Data


The Sample Databases are most useful during the trial period but can be accessed at any time.

Selecting the MS Access Sample Database on First Start

When Calibration Control starts for the first time there are two database formats to choose from, MS Access or SQLServer. Click the [MS Access] button to continue.




In the next dialog (Choose Database Type), check the "CREATE New SAMPLE Database" checkbox and click the [Continue] button.




The Trial Period Notification dialog pops up each time Calibration Control starts without a current Product Key. If a Product

Key has been purchased click the [Yes] button to use it in the Sample Database. Otherwise, click the [No] button to continue.



When Calibration Control starts, about 300 sample records will be visible.




The screenshot shows the software interface with the "Calibration" tab selected in the ribbon. The "Equipment" table is visible, showing 262 records (58 past due). The table has columns for Equipment ID, Model, Description, Serial, and Manufacturer. The first seven rows are highlighted in alternating colors (yellow, pink, green, grey, pink, grey, pink).


Equipment ID	Model	Description	Serial	Man
SAMPLE-001	BMS10-83	Type III	112	
SAMPLE-002	IDR-8111	Indicator	2213	East
SAMPLE-003	GBS0-5	Gage Block Set		
SAMPLE-004	CRG-101	Cylindrical Ring G	8411	
SAMPLE-005	CRG-101	Cylindrical Ring G	199	
SAMPLE-006		10-32 UNF-3B	1008	
SAMPLE-007	CRG-101	Cylindrical Ring G	195	

Selecting Sample Database After Program Start

Open the Options dialog (from either the File or Utilities tab of the ribbon menu) and navigate to the Advanced tab. Click the **Switch to Sample Database** button.



Then click the [Yes] button in the following 'Show Sample Database' dialog.



Reset Configuration


If you selected MS Access and now you'd like to try using a SQL database or vice versa, you'll first need to reset the configuration. You can use the keyboard shortcut by holding CTRL and launch Calibration Control OR you can manually reset the configuration by navigating to this path on your computer:

C:\Users\Public\Documents\Ape Software\Calibration Control

Then, delete the general.config file and restart Calibration Control. You should now see the Choose Database Format dialog.

Sample SQL Database

To get started using our Sample SQL Server database, choose SQL Server from the Choose Database Format dialog, then select SAMPLE DB in Cloud. If you'd like a blank database to try, first ensure you can connect to our Sample DB in Cloud, then email us with the subject "Blank SQL Database Request" and include your full name, company name and company address/general location. We'll set up a blank database just for you and send you the connection instructions.



Remove Sample Data

Remove Sample Data to Begin Entering Your OwnRecords

If, during the trial period of Calibration Control (our **Calibration Management Software**), the sample database is used to evaluate the software, the sample data should be removed before entering live equipment records.

Double Click (Automatic Method)


The easiest way to remove the sample data is to double-click on the graphic in the bottom-right of the Equipment grid . . .

r	Yearly	0	04/20/2019	04/20/20
	Yearly	0	12/26/2018	12/26/20
	Yearly	0	05/24/2018	05/24/20
	Yearly	0	05/24/2018	05/24/20
	Yearly	0	05/24/2018	05/24/20
r	Yearly	0		20
	Yearly	0		20
n	Months	8		20


SAMPLE DATA
DOUBLE-CLICK HERE
FOR BLANK DATABASE

General Configuration

... and then click the [Yes] button.



After the application restarts, select the 'CREATE New Blank Database' or 'BROWSE to Existing Database' checkbox, then click the [Continue] button.



Note: If the CREATE New Blank Database option is selected, then the new database will be created at this path C:\Users\Public\Documents\Ape Software\Calibration Control on the computer being used. Therefore, if the database needs to be moved to a network location, follow the instructions in the move database help topic.

If you need any help, contact Ape Software and we can do a remote support screen share to do this together.

Move Application Database

Move Application Database to a New Location

Watch this short video on how to move an application database file (apecal.mdb) in Calibration Control (our **Calibration Management Software**) from its initial location, to a new location to make the database easier to backup or to place it in amore central location to share with additional users.

STEP 1: Close Application

Before moving the database file, ensure the application is closed.

STEP 2: Find Database File

Before moving the 'apecal.mdb' file, find its location. If you need help, see the help page on finding your database filebefore moving to the next step.

STEP 3: Move Database File

Choose the new location either on a shared network or the same computer and move the 'apecal.mdb' file from its originallocation to the new location.


Note 1: If moving the database file to a network location, grant appropriate user permissions to that new locationdepending on who needs access.

Note 2: When moving a database, rename or backup to a different location the database from the first location. The goalis to ensure that the program does not mistakenly connect to the old database sometime in the future thereby making the users assume that data is lost (old db has fewer records).


STEP 4: Point Application to New Location

Point the Ape Application to the new location by starting the application. When the application discovers the 'apecal.mdb' file is missing from its known location, it will show the following 'Choose Database Type' dialog:

From within this dialog, select 'BROWSE to Existing Database' and click the [Continue] button.



Browse to and select the new database location and click the [Open] button.



The Ape Database Software should start as normal and will remember this new location.

Locating Database File (MS Access)


Find the Database and Configuration Files

By default, the apecal.mdb database is stored in the Settings Folder unless the database was moved to another location, like a network folder. The default Data & Settings Folders is different for different versions of Calibration Control and different operation systems. Refer to the Folders Help Topic for an overview of the ones used by Ape Database Software.

Method 1 - Version 5.5 and Higher

This is the best method to find the database if the application **can be opened**.

The easiest way to find the Data Folder in versions 5.5 and higher is to select the Open Data Folder Option from the Filedropdown menu.



Otherwise, the default Data Folder is located at:

C:\Users\Public\Documents\Ape Software\Calibration Control

Method 2 - Version 8.0 and Higher

This is the best method to find the database if the application **cannot be opened**.

Versions 8.0.4 and higher automatically create a Connection Log file (**connectionlog____.csv**) found in the Settings Folder(see section above). The Connection Log can be read with a spreadsheet program, like MS Excel. This file logs each of the connections to the database (e.g., path). A blank path statement in the connection log means that the database was located in the default location (i.e., Settings Folder).

Method 3 - General Config File

If the apecal.mdb file is not located in the folders described above, the path to the database can be found inside the general.config file located in the Settings folder. After finding the general.config file, open it either by double-clicking or using a common text editor like Notepad.



```
general.config X
<?xml version="1.0" encoding="utf-8"?>
<configuration>
  <Global>
    <tblGlobal key="ConfigName" value="Access_Sample_Tester" />
    <tblGlobal key="ConfigColor" value="-41378" />
    <tblGlobal key="ConnectionType" value="OleDb" />
    <tblGlobal key="Server" value="" />
    <tblGlobal key="CustomDbName" value="apecal-sample_db_2016-06-20.mdb" />
    <tblGlobal key="ConnectionString" value="" />
    <tblGlobal key="Ace" value="True" />
    <tblGlobal key="OleDbServices" value="False" />
    <tblGlobal key="AutoSqlLogin" value="False" />
    <tblGlobal key="WindowsAuth" value="False" />
    <tblGlobal key="SqlUserName" value="" />
    <tblGlobal key="SqlPassword" value="" />
    <tblGlobal key="ApeUserName" value="" />
    <tblGlobal key="ApePassword" value="" />
    <tblGlobal key="FilesPath" value="" />
    <tblGlobal key="LCID" value="en-US" />
    <tblGlobal key="Email" value="" />
    <tblGlobal key="CrashEmail" value="False" />
    <tblGlobal key="CrashNotifyApe" value="False" />
    <tblGlobal key="DatabasePath" value="c:\users\public\documents\ape software\calibration control\" />
    <tblGlobal key="AccessPassword" value="" />
    <tblGlobal key="DefaultLabelPrinter" value="" />
    <tblGlobal key="CalibrationOrderBy" value="" />
    <tblGlobal key="CalibrationWhere" value="" />
    <tblGlobal key="TopCalRecordCount" value="" />
    <tblGlobal key="ViewStatusBar" value="True" />
  </Global>
</configuration>
```

Data Imported To MS Access

How to Deploy Your Imported Data

This help topic only applies to the [MS Access] database format. If using SQL Server, refer to the Add Ape Database to SQL Server help topic. If Ape Software performed your data import, you will receive your imported database in a zip file through either a download link or an email attachment. To deploy the database yourself, follow the instructions below. Otherwise, contact us and we'll assist with the process through a screen share.


Install Current Version

Before proceeding, download and install the most current version of Calibration Control. Ensure the Minimum Requirements are present before installation.


Deploy Database (Scenario A)

Proceed to Scenario B (below) if the Application has already been used on this computer to access a database. Otherwise, continue with this section.

- 1. Download & Unzip:** Download, unzip, and place the imported database in the desired folder. This may either be on the local computer or a network location, if the database will be accessed by multiple users.
- 2. Choose Database Format:** Start the Application and select the [MS Access] button. If the Application does not start with the Choose Database Format dialog, proceed to the next section (below) **Remove Old Database (Scenario B)**



3. **Choose Database Type:** Select [BROWSE to Existing Database] from the Choose Database Type window.




4. **Select Database:** Navigate to the database you downloaded and unzipped (step 1), select, and click the [Open] button and the software should open with the imported data visible.

5. **Common Files Folder:** If the database was placed in a network location and/or multiple users will access the software, refer to the **Common Files Folder** section of Creating a Multi-User Environment to setup the Files Folder.

Remove Old Database (Scenario B)

If the Application was previously used with a Blank or Sample database, the files already present should be removed or renamed before proceeding.

1. **Open Data Folder:** After the Application is running, open the Data Folder by clicking [File] in the top-left of the screen and selecting "Open Data Folder". Confirm the existence of the database file (apecal.mdb or apecal_sample.mdb). This is the quickest and most accurate way of finding the current database. Leave this window open and proceed to the next step.



2. **Open Settings Folder:** In the File menu, open the Settings Folder. If the Settings Folder is the same folder location as the Data Folder, you can skip this step. Note that by default, the Data Folder and the Settings Folder are the same location unless the database is moved.
3. **Close Program:** Close the Application so the database files can be moved, renamed, and/or replaced.
4. **Rename or Move Files:** Although the data files can be renamed, it's easier to move ALL the files in the Settings and Data folders to a backup location (e.g., "Data Backup - Current Date"). Cleaning out these folders will create the conditions required for a clean first-run of the imported data and will prevent Calibration Control from automatically connecting to a database file in the Settings folder location.
5. **Scenario A:** Proceed with Deploy Database (Scenario A).

Additional Help

If you have any general questions about the import process itself, refer to the Database Import Process article. Other useful help topics:

- Creating a Multi-User Environment
- Move Calibration Management Database

Creating a Multi-User Environment

Setting up Multiple Users

Setting up access for multiple users to use Calibration Control (our **Calibration Management Software**) is as easy as deleting or renaming the general.config file, starting the application (or hold CTRL while starting Calibration Control), and linking to the new file. Therefore, if you are using an Access database file, you should be familiar with the help topics for finding your database and moving your database. If you are using a SQL Server, then take a look at the help topics on Adding an Ape database to SQL Server and Connecting to a SQL Server Database.


Move Your Database (MS Access Database Only)

Find your database and move your database to a network location where all users have read-write access. Most of the difficulties users experience with implementing a multi-user environment are related to network folder permissions. So if you experience any difficulties, you should first seek assistance from your network administrator.

Set Up a Common Files Folder (Access & SQL)

Calibration Control uses a Common Files Folder that contains four sub folders (i.e., Attachments, Emails, Labels, and Reports). Ensure each networked user is using the Common Files Folder so its sub folders are equally accessible to all users. Do this by:

1. **Common Mapping:** Ensure each user has the same network mapping for the Files Folder. For instance, if a mapped drive (e.g., "**M:**\\Server2\\Common Files") is used for one user, the same path will be used for all users.
2. **Point Ape Database Software:** Open the Advanced tab of the Options dialog and click the Folder & Path



button. Change the 'Location of Files' (i.e., Common Files Folder) to the common network location (step 1).


3. **Move Folders:** If custom files already exist in the old location(s) of the Files Folder, move those files to their new Common Files Folder locations (step 1). Specifically, move all attachments and custom templates (e.g., reports, labels, & emails). Otherwise, the default Reports and Labels templates will be copied to the new Files Folder Location.

Add Additional Users

When running the software on a client for the first time, it presents the option to create a database (blank or with

samplerrecords) or to navigate to an existing database. Therefore, after starting the software on a new computer, check the 'BROWSE to Existing Database' option and [Continue] if a database already exists in a network location.

Otherwise, if a database has already been created on a client that needs to be networked (i.e., linked to a database on the network), delete or rename the local 'apecal.mdb' file to prevent the program from automatically connecting to it. If already running Calibration Control, find the folder for this database file by selecting Open Data Folder from the Files dropdown menu. Otherwise, refer to the finding your database help topic.



Add Ape Database to SQL Server

Install Calibration Control on a SQL Server

Before beginning, ensure you are signed into a SQL Server with enough permissions to create and update databases and their objects (e.g., system admin).


Installing the apecal database on an existing SQL Server (2012 or higher) is a simple database restore from a provided backed up file. All SQL files needed to perform the install reside in the SQL Tools folder under the CC program folder located at Program Files (x86)\Ape Software\Calibration Control.

There are two backup files (.bak) located in the SQL Tools folder. One is the sample database named `apecal_sample_{structure/records versions}.bak` (i.e. `apecal_sample_s84_r86.bak`) and the other is a blank database named `apecal_{structure/records versions}.bak` (i.e. `apecal_s84_r86.bak`). Copy the backup file you wish to restore to your SQL Server default backup location.


Restore Back-Up Database

Follow along with our Backup and Restore SQL Server Database help topic for more detailed instructions. From within SQL Server Management Studio with the SQL Server connected, right-click on the Database node and select 'Restore

Database'.




From within the Restore Database dialog, type the name apecal in the 'Destination' field, select 'Device' as the Source, and then click on the [. . .] button to the right of the 'Device' field.




From within the Specify Backup dialog, click the [Add] button.

Locate the backup file by navigating to or pasting in the 'Selected path'. Paste or type the name of the file to restore from (e.g., apecal_s51_r52.bak). Then click [OK] to save and close the backup file and [OK] again to close the Specify Backup dialog.



When the Restore Database dialog is on the top again, click the checkbox in the Restore column of the backup set (file) you just designated and then click the [OK] button to begin the restore.


If successful, the following confirmation dialog should display:



After closing the confirmation dialog, refresh the list of databases to see the new apecal database.

Restore Failed

If the restore failed, you may need to repeat the above steps and specify different file names or path statements. Use the Options Page of the Restore Database dialog to modify the Database (e.g., Rows Data) or Log file names or the Path statement (i.e., Restore As).




User Security

Users can connect to the apecal SQL Server database using either Windows or SQL Server authentication. Each user will either need Database Roles or Explicit Permissions.

When using Database Roles, ensure each user has **db_datareader** and **db_datawriter** roles.

When using Explicit SQL Server Permissions, ensure each user is granted at least the Connect, Delete, Execute, Insert, Select, and Update permissions.



If you need help, contact Ape and we can setup a screen share and do this together.

Connect to SQL Server Database

Configure Calibration Control to Connect to a SQL Server


After creating the Calibration Control database on your SQL Server, follow these instructions to connect Calibration Control (our **Calibration Management Software**) to your SQL database.

Start Calibration Control


When running Calibration Control for the first time, select **MS SQL Server** as the database type.



To connect to your database, select **Custom DB**. Alternatively, we offer cloud hosted databases (sample data or blank) that you can easily connect to without further configuration.



If you are already connected to a SQL Server database or an MS Access database, you can select the 'SQL Connect' option from the Utilities tab in the ribbon menu.



Enter the fully qualified name into the Server Name field, confirm the database name, and select your authentication method. Click the **Test Connection** button to confirm that you have entered the correct connection information.

The Remember Login checkbox is automatically checked. If it is deselected, you will be prompted to verify your server authentication method the next time you start Calibration Control. Once your connection test is successful, click the **Submit** button.

The screenshot shows a dialog box titled "SQL Server Connection (CC Version 9.2.3.0)". It has the following fields and controls:

- Server Name:** A text input field.
- Database Name:** A text input field.
- ... of ...:** A small text label between the Database Name and Connection String fields.
- Connection String:** A large text input area.
- Authentication:** Two radio buttons: "Windows Authentication" (selected) and "SQL Server Authentication".
- User Name:** A text input field.
- Password:** A text input field with a "Remember" checkbox to its right.
- Buttons:** "Test Connection" (centered), "Cancel" (bottom left), and "Submit" (bottom right).

Restarting the Process

This entire process creates a configuration file in the Settings folder on the computer where the application is running. If a mistake is made during the process, the best way to restart from scratch is to delete the config file and start over. You can do this easily by holding down CTRL while starting Calibration Control (Version 9.1.2 and later).

You can also delete the configuration file manually by closing the application and opening the Settings folder, which you can find at this path on your computer: C:\Users\Public\Documents\Ape Software\Calibration Control


Then, delete the **general.config** file and restart Calibration Control.

Connection Troubleshooting

If needed, refer to the help topic on Troubleshooting a SQL Server Connection.

Calibration Control Folders

Each of the folders (Data, Files, Settings, & Web) are accessible from the File drop-down menu in Calibration Control (our **Calibration Management Software**).




What They Do

The **DATA** Folder is the location where the access database (when using MS Access) resides. Visit the Move Database Folder help topic for steps on how to move the Data Folder to a new location.

The **FILES** Folder stores the Attachments folder and the Labels, Reports, and Emails template folders.

The **SETTINGS** Folder stores the program settings that Ape Software requires to function normally. This is the only FileFolder that cannot be moved from its original location and every user must have read/write access.

The **WEB** Folder is the target location for publishing the reporting web pages.




Settings Folder

The Settings folder is the only folder that cannot be moved. The location of the Settings Folder is:

Windows 7 & Higher: C:\Users\Public\Documents\Ape Software\Calibration Control

Although all Folders share the same default location as the Settings folder when the software is first installed, the other three Folders can be moved to different locations on the computer or local network for sharing and backup purposes. The location of the Files folder can be configured in program Options > Advanced tab > Folder & Path.



Grid Features


Features Common to All Data Grids

Each grid in Calibration Control (our **Calibration Management Software**) shares the same interface and basic features. Navigate data with sorting options, filters, grouping, and more.

Context Menu: Find more options by right-clicking anywhere in a grid for the context menu. Selections may vary by grid.

- **Alternate Grid Views:** Most Data Grids have additional context menu items located at the bottom of the menu. Commonly, these additional menu items will be alternate views of that data grid. For example, the Equipment Data grid has 4 additional menu items: All Equipment, All Equipment - Show Hidden, Equipment Due Cal and

EquipmentMaintenance Due which also can be accessed by pushing the corresponding key (e.g. F7).



- ◆ **Export:** Right-click in the grid and in the context menu, select the option to Export to either Excel or PDF.

	Print Certificate(s)		CMM-421
	Open Equipment Record		
	Open Record	Ctrl+O	CRG-101
	New Record	Ctrl+N	IDR-8111
	Delete Record	Ctrl+D	IDCM3
	Select All Rows	Ctrl+R	GBS0-5
	Copy to Clipboard	Ctrl+C	RGS263e
	Select All in Cell	Ctrl+A	SJSWWTest
<input checked="" type="checkbox"/>	Filter Row		CRG-101
	Show Fields		CRG-101
	GroupBy Box		01-MAM
	Row Sizing		DTM
<input checked="" type="checkbox"/>	Show Top N		DRPI-57
	Refresh Grid	F5	ODM-d80b85
	Reset Grid to Default		DRPI-57
	Print Data Grid	Ctrl+P	
	Export		Excel (*.xls)
			Excel (*.xlsx)
			PDF
			XPS
AGE-021	5/07/19	1225	
AGE-027	5/07/19	1222	
AGE-009	1/30/19	1178	

- **Change Font Size:** Change the grid's font size by clicking on the 'Row Sizing' option and entering the size value in the Font Size box. The default font size shows 0, but is actually equal to 8, and numbers larger than 8 will yield larger font sizes.

Column Change: Drag any column header to change its placement on the grid. Alternatively, use the small down arrow next to each header to select which column you want changed to that space.

Context Sort: Click directly on a column header to sort records in ascending or descending order for that column. Hold the shift key while clicking on additional column headers to create sorts with multiple columns.

Refresh: Refresh the data by pressing [F5] or by right-clicking the grid and selecting Refresh Grid from the context menu.

Add New: Create a new record by double-clicking in the gray area of the grid, pressing [Ctrl+N], or by right-clicking the grid and selecting New Record from the context menu.

Open and Edit: Open a record by first selecting the record and then double-clicking, pressing [Ctrl+O], or by right-clicking a record and selecting Open Record from the context menu.

Delete: Delete a record by first selecting the record and then pressing [Ctrl+D] or by right-clicking the selected record and selecting Delete Record from the context menu.

Display Count: Displays how many records are listed. In the Equipment grid, the count for Equipment Past Due will also be displayed.

- **Records:** The number of records contained within the grid is indicated in the Grid Header.

Equipment					
[Equipment] - 262 Records					
Equipment ID	Model	Description	Serial	Manu	
SAMPLE-020	DRPI-57	Drop Indicator	2119		
SAMPI F-021	IDR-8111	indicator	305		

- **Filtered:** The number of filtered records is also indicated within the Grid Header.

Equipment					
[Equipment] - 126 FILTERED Records					
Equipment ID	Model	Description	Serial	Manu	
SAMPLE-017	CMM-421	CMM	186		
SAMPI F-019	ODM-d80h85	OD Micrometer	386		

Filter Row: The Filter Row is a helpful tool for searching and reorganizing data.

- When the Filter Row is enabled from the context menu, a blank row will appear at the top of each column.
- Filters can be added to multiple columns at the same time.
- Create custom (i.e., combination) filters within a single column by selecting 'Custom' from the drop-down menu.

Equipment ID	Model	Description	Serial	Manufacturer	
SAMPLE-001	BMS10-83	Type III	112		
SAMPLE-002	IDR-8111	Indicator	2213	East Coast Cal	
SAMPLE-003	GBS0-5	Gage Block Set			
SAMPLE-004	CRG-101	Cylindrical Ring G	8411		
SAMPLE-005	CRG-101	Cylindrical Ring G	199		
SAMPLE-006		10-32 UNF-3B	1008		
SAMPLE-007	CRG-101	Cylindrical Ring G	195		


Clear Filters: Clear a single column filter by placing the cursor in the filter field for that column and pressing the clear filter button (looks like a crossed out filter). Clear all filters by selecting the clear filter button at the far left of the filter row.

Pinned Sort Column: Click the Pinned Sort button (looks like a thumb tac) located at the top right corner of the column header and the selected column will be "pinned" to the left hand side of the grid. When horizontally scrolling through a grid, a pinned column will stay visible.

Grid Splitter: Use the Grid Splitter to view multiple sections of the grid at the same time by dragging the Grid Splitter bar to the desired height of the grid. And the grid can be split as many times as needed. Remove a Grid Splitter by dragging the bar all the way to the top or bottom of the grid.

Show Fields: Add and remove columns (i.e., fields) by right-clicking and selecting Show Fields from the context menu.

- Check the desired field columns to [Add] from available fields or [Remove] from visible fields.
- Also, change the column order by selecting a visible field and clicking the [Move Up] or [Move Down] buttons.



Grouping Data: The GroupBy Box view is great for on-the-fly organization with “mini-data sets”, giving you the ability to arrange records by their columns.

- When the GroupBy Box mode is enabled from the context menu, a gray 'Drag To' section (the GroupBy Box area) appears at the top of the grid.
- Select and drag column headers (fields) into the gray area to create custom groupings.
- Remove any fields from the header by simply dragging them up and out of the gray boxed area and releasing.
- And if desired, the GroupBy mode can be disabled with records still organized in custom groupings.


Highlight Colors of the Equipment Grid

Understanding the Highlighted Colors of Equipment Records

The Equipment Grid currently uses five default colors to visually identify 'past due' (pink), 'calibration/maintenance due soon' (yellow), 'calibration standard' (green), 'calibration failed on receipt' (red), and 'received for calibration' (gray). Colors can be customized.


Showing/Hiding Colors

Show or hide highlights by right-clicking for the context menu in the Equipment grid and toggling any of the five options in the 'Highlight' section.



Default Colors

Manage the default highlight colors by clicking the [Default Colors] button in the General tab of the Program Options dialog. Change the colors to ones that makes sense for your company, or click the [Reset] button to restore the default colors.



Highlight Past Due

Equipment records with a Next Cal date less than the current date are highlighted Pink. Records marked as Calibration /Test Standards (usually Green) are also marked as Pink when they are Past Due.

Equipment records not past due but “due soon” are highlighted Yellow. Due soon is defined by the Reminder Lead Daysfield in Options > Advanced tab > Calibrations. By default the Reminder Lead Days field is set to 30.

Calibration Options

Status on Calibration Pass: Active

Status on Calibration Fail: Inactive

Reminder Lead Days: 30

Received for Calibration Status: In Calibration

Manually Set Due Dates In Equipment Dialog

Ask to Print OOT Worksheet When Received OOT

Null Cal Received Date on Successful Cal

Use Last Calibration as Template if No Template Linked

Prevent Entry of Future Date for Calibration Events

Block Printing of Due Cal Labels for Failed Calibrations


Work Days

S M T W T F S

Cancel Save


Highlight Standards

Equipment marked as Calibration / Test Standard are highlighted Green unless they are also Past Due or soon to be Due Cal. Mark Equipment as a Standard by checking the 'Is Calibration / Test Standard' check box in the Calibrations tab of the Equipment dialog.




Highlight Received Out of Tolerance

If any of the past Calibrations of a given Equipment have an As Found status of Out of Tolerance (OOT), the record is highlighted red in the main Equipment grid.




Remove the Red highlight of an Equipment record by marking all of the Calibration Events with an 'As Found OOT' status as 'OOT Investigated'.



Received for Calibration

Equipment that has been received for calibration will have a gray highlight in the grid.



Filtering Records Using the Filter Row


How to Use the Filter Row in the Data Grids

The Filter Row is a helpful tool for searching and reorganizing your data in Calibration Control (our **CalibrationManagement Software**).

Display Filter Row

By default, the filter row is already displayed at the top of the grid; it looks a bit like a blank row. If the filter row is

notdisplaying, right-click the grid and select Filter row from the context menu.




Equipment ID	Model	Description	Serial	Manufacturer
SAMPLE-001	BMS10-83	Type III	112	
SAMPLE-002	IDR-8111	Indicator	2213	East Coast Cal
SAMPLE-003	GBS0-5	Gage Block Set		
SAMPLE-004	CRG-101	Cylindrical Ring G	8411	
SAMPLE-005	CRG-101	Cylindrical Ring G	199	
SAMPLE-006		10-32 UNF-3B	1008	
SAMPLE-007	CRG-101	Cylindrical Ring G	195	

- When the Filter Row is enabled from the context menu, a blank row will appear at the top of each column.
- Filters can be added to multiple columns at the same time.

Create New Filter

Choose which column to add a filter to and either use the drop-down menu to search for the filter, or type it into the blankfield. Results will immediately appear as you type.



Model	Description	Serial	Department
LT	Loctite	321	(Custom)
BMS10-83	Type III	113	(Blanks)
BMS10-83	Type II	110	(NonBlanks)
	GRAY 16515	2146	Engineering
BMS10-83	Type III	111	Inspection
	10-32 UNF-3B	1008	Machine Shop
			Manufacturing
			Production
			Manufacturing

Clear Filter(s)

Single Filter: Clear a filter within a single column by placing your cursor within that column's filter field (in the filter row) and clicking on the [clear filter] button on the right side of the field. The [clear filter] button is an icon of a filter with a linethrough it.

All Filters: Clear all filters by clicking the [clear filter] button on the *far left* side of the filter row.

Serial	Department	Manufacturer	Location
	St		
321	Stockroom		
110	Stockroom		Shelf
5060372002	Stockroom		Shelf
310	Stockroom		Shelf
124	Stockroom		Shelf
139	Stockroom		Shelf
363	Stockroom		Shelf

Filter Methods

Additional filter operators (i.e., filter method) are available by clicking on the blue square to the left of each filter field.

Equipment ID	Model	Description	St
SAMPLE-001	BMS10-		3
SAMPLE-002	IDR-811		13
SAMPLE-003	GBS0-5		
SAMPLE-004	CRG-10		11
SAMPLE-005	CRG-10		9
SAMPLE-006		10-32 UNF-3B	1008

Custom Filters

Create combination filters within single columns by clicking the drop-down button (downward arrow) within the filter field of that column and selecting (Custom).

Custom Filter ✕

Filter based on All of the following conditions:

+	Add	Serial	Starts with	M
-	Delete	Serial	Does not equal	[Description]


Hiding Equipment Records by Status

Use Equipment Status to Manage Visibility


It's a good practice to not delete equipment records in Calibration Control even after the equipment is no longer in use. Keeping out-of-use equipment visible can clutter up the Equipment grid and Due Cal reports making them more difficult to use. To help fix this problem, choose which status codes should be visible and which should be hidden.

Status Codes Visibility

Select the Codes grid from the Data Grids tab of the ribbon menu.




The status codes under the 'Equipment' grouping includes an [Is Hidden] checkbox that, when checked, hides the equipment records with this status code and they will no longer be visible in the Equipment grid or default Due Cal reports. Uncheck the checkbox to make the records with that status visible again.



Select a record's status by using the drop-down field in that record's main Equipment dialog. Adding a NEW status code is easy, just click on the shaded Status box and create a new status code.

Show/Hide All Hidden Records

From the Equipment grid, right-click for the context menu and toggle the [All Equipment - Show Hidden] option to show or hide all records with a hidden Equipment Status. The Equipment grid will show an [Equipment Including Hidden] status as the header to remind users that there are records being hidden.



Manage Out-Of-Tolerance (OOT) Equipment

Manage and Investigate OOT Equipment

When performing calibrations on measuring equipment in Calibration Control (our **Calibration Management Software**), it is proper to check the condition of the equipment (aka unit-under-test or UUT) before calibration (As Found) and after calibration (As Left). If measurement equipment is found out-of-tolerance (OOT) before performing calibration, an assessment of the effect of the OOT condition on processes must be performed. If it has been 12 months since the last calibration of the UUT, all the product and processes the UUT measured in the last 12 months are suspect. For instance, were products or processes falsely judged good when they were bad?

Calibration Event

In a Calibration record, if the As Found field is OOT, the OOT 'Investigated' checkbox and the [Print OOT Form] button become visible. Upon finding OOT equipment, conduct an investigation in accordance with your organization's corrective and preventive action (CAPA) or nonconforming material (NCM) processes. After the investigation is complete, check the 'OOT Investigated' checkbox.

If useful to your organization's CAPA or NCM processes, print the OOT Form from the Calibration record for use as an investigative tool.

The screenshot shows the 'Edit Calibration Event' window with the following details:

- Certificate:** 1186
- Temp:** 67.00 °F
- As Found:** Out of Tolerance
- Date:** 02/22/2019
- Humidity:** 46.00%
- As Left:** Pass
- Company:** [Empty]
- Technician:** David Wilcox
- Print OOT Form:** [Button]
- OOT Investigated:**

Remarks: [Empty text area]

Make Template [Button]

Group	Units	Plus Tolerance	Minus Toleranc	As Found	As Left	Measurements	Tolerance Type
2		0	0	Out of Tolerance	Pass	1	Tolerance (+/-)

Buttons: Cancel, Print Certificate, Submit

Red Highlighted Records

Calibration Events with pending OOT investigations are highlighted red in the Equipment grid when the 'Highlight - Received Out of Tolerance' option is enabled from the context menu. Clear the red highlighting by checking the 'OOT Investigated' checkbox for all Calibration Events (see above) that have OOT conditions for a given Equipment record.

Equipment				
[Equipment] - 262 Records				
Equipment ID	Location	Manufacturer	Units	Next Cal
SAMPLE-001	Shelf		2	11/29/2020
SAMPLE-002	Shelf	East Coast Cal	0	08/19/2019
SAMPLE-003	Shelf		0	05/03/2020
SAMPLE-004	Shelf		0	04/26/2020
SAMPLE-005	Shelf		0	04/26/2020
SAMPLE-006	Test Area		0	07/03/2019
SAMPLE-007	Shelf		5	05/31/2019
SAMPLE-008	Shelf		0	04/26/2020

Highlight	>	Due Cal
Check Out/In		Maintenance Due
Asset Transfer		Calibration Standards
Scan Barcode	→	<input checked="" type="checkbox"/> Received Out of Tolerance
Open Record	Ctrl+O	Received for Service

Calibration History


Calibration Event Quick Reference and Filter

Use the Calibration History Grid in Calibration Control (our **Calibration Management Software**) to quickly find Calibration Events without having to go through their Equipment record.

Top N Records

By default, Calibration Control limits the default display of Calibration records to the top 100 in descending Calibration Dateorder.

Change the number of records displayed by selecting from the drop down menu showing 'Select Top 100'. Available fields: Select Top 100, 500, 1000, 2000, 3000, 4000, 5000, 10,000 or Select All.



Filtering and Sorting

Enter Filters and Sorts just like any other data grid. Refer to the Features Common to All Data Grids help topic for additional instructions. The difference here is to remember that the Calibration Grid is only displaying the Top N records (see above). This means that the Filters and Sorts only apply to the collection of records currently displayed.

To apply a Filter or Sort the entire list of records, select "Select All" from the dropdown in the grid.

Context Menu

In the Calibration History grid, right-click a record to open the context menu and quickly Print Certificates or Open the Related Equipment record.

Certificate	Date	Equipment ID	Model
CC-0055	06/29/2020	SAMPLE-006	
CC-0057	06/29/2020	SAMPLE-006	
▶ CC-0104			LT
CC-0058			
CC-0310			DRPI-57
CC-0183			74c284-ODM
CC-0165			979e_IDR
CC-0196			BMS10-83
CC-0257			01-MAM
CC-0200			987f_ODM
CC-0164			CAL-180
CC-0271			DTM
CC-0180			IDR_7b8
CC-0170			ODM-d80b85
CC-0243			HG66
CC-0103			DRPI-57
CC-0248			CAL-180
CC-0261			ODM-d80b85
CC-0255	06/07/2020	SAMPLE-019	ODM-d80b85

Print Certificate(s)

Open Equipment Record

Open Record Ctrl+O

New Record Ctrl+N

Delete Record Ctrl+D

Select All Rows Ctrl+R

Copy to Clipboard Ctrl+C

Select All in Cell Ctrl+A

Filter Row

Show Fields

GroupBy Box

Row Sizing

Show Top N

Refresh Grid F5

Reset Grid to Default

Print Data Grid Ctrl+P

Export >

Equipment Dialog

Edit Equipment Records

Use the Equipment dialog to edit all data related to the tool / test equipment. Change the name of any field to continue using words that make sense to your organization. If the field name* is in a shaded box, click on it to add new values in the drop-down list of available values (combos).

Equipment Dialog Fields

- **Equipment ID:** The unique identification given to your own tools / test equipment and is often referred to as an Asset Number. This is the only required field for a new record.
- **Serial Number:** The unique identification the manufacturer gives their tools / test equipment. (You may periodically find that some manufacturers do not assign serial numbers or that a serial number label has fallen off.)
- **Model*:** The model number of the tool / test equipment usually assigned by the manufacturer.
- **Size/Range*:** The size or range of the equipment.
- **Description:** Description of the tool model number, which automatically changes to reflect the description of the model field (above). After a Model number is selected, the Description field can be edited.

- **Manufacturer*:** The company that made the equipment/tool. This field is also automatically set if an existing ModelNumber from the Model combo box is chosen.
- **Equip Type*:** Classification system developed by you, the user of Calibration Control. (Examples of Equipment Types include Caliper, 6" Caliper, ThermoTron, Oven, etc.) This field will automatically populate if an existing Model Number is chosen and has a Type field assigned.
- **Alternate ID:** This is a second ID field, the Equipment ID being the first. Use this field to record a company Assetnumber, if it is not the same as the Equipment ID, or a customer's Equipment ID.
- **Site*:** Company field for a geographic Site or Company.
- **Department*:** The department in the organization where the tool is currently located.
- **Location*:** The location within the department where the tool is currently located. (This could be a specific workarea or even an engineer's desk.)
- **Custodian*:** The individual who has custody or possession of the Equipment/Tool.
- **Status*:** The status of the tool / test equipment (e.g., accepted, removed from service, failed calibration, or anyother status needed).
- **Notes:** You can use the Notes field for notes about the Equipment itself or use for whatever purpose needed.

Calibrations Tab of Equipment Record

- **Calibrations Panel Grid:** Click on the [+] button at the bottom of the dialog to create a new Calibration Event.
- **[Calibration] Last:** The date of the last calibration is set automatically when entering a new Calibration Event with aStatus of Pass. The Calibration Last date can be set manually, which automatically updates the Calibration Next field based on the value of the Frequency field.

- **[Calibration] Next:** The date of the next calibration is calculated automatically when entering a new passing calibration event or when manually updating the Calibration Last field. In either situation, the next calibration date is calculated based on the chosen Frequency and Units.
- **Override Calibration Next [check box]:** Manually extend or otherwise override the Calibration Next date by checking the box to the right of the Calibration Next field.
- **Received [for calibration]:** The date when equipment is received for Calibration.
- **[Print Certificate] button:** Click this button to send a Calibration Certificate worksheet straight to your default printer.
- **Frequency / Units:** The frequency between calibrations, like yearly or semi-annual. If choosing a frequency that requires a unit multiplier (e.g., weeks, days, or months), the Units field will enable to adjust the number of units (weeks, days, or months, etc.). The 'Month of' and 'Week of' frequencies set a general due date of an entire month or week.
- **Calibration Standard [check box]:** Check this box if the current equipment is used as a calibration standard for other equipment.
- **Company*:** The organization (Calibration Company) that regularly calibrates this instrument. This is a useful field when wanting to create a report of all the equipment due in a given period for a specific calibration resource.
- **Technician*:** The person who normally calibrates this equipment.
- **Certificate:** Number of the certificate issued by the organization responsible for the last calibration event. Used primarily to aid in traceability back to a national measurement standard, especially if this equipment is a Calibration/Test Standard.
- **Template:** Calibration Template used when creating new calibration events.


The screenshot shows the 'Edit Equipment' window with the 'Calibrations' tab selected. The form includes fields for 'Last' (01/17/2019), 'Next' (04/17/2019), and 'Received' (09/11/2018). It also has dropdowns for 'Frequency' (Months) and 'Units' (3), and text boxes for 'Company' (FICT CO.) and 'Technician' (David Wilcox). A 'Print Certificate' button is visible. Below the form is a table titled '[Calibrations] - 1 Record' with columns: Certificate, Date, As Found, As Left, Cal Company, and Technician. The table contains one record with Certificate 1175, Date 01/17/2019, As Found 'In Tolerance', and As Left 'Pass'. At the bottom are buttons for 'Cancel', 'Cal Label', 'Worksheet', 'Asset Tag', and 'Submit'.

Certificate	Date	As Found	As Left	Cal Company	Technician
1175	01/17/2019	In Tolerance	Pass		David Wilcox

Maintenance


Use the Maintenance tab to add or link maintenance requests to an Equipment record. Keep track of equipment with maintenance due with grid highlights and/or the Equipment Maintenance Due grid (right-click the Equipment grid or press

F12).



Attachments | Documents | Notes | Parts Tab

Use these panel grids to add, link, and edit the Attachments, Documents, Notes, and Part Numbers records related to an Equipment record.



Misc (Miscellaneous) Tab

The Miscellaneous tab contains groups of data that help organize information related to the Equipment record.

- **Miscellaneous Group:**
 - **Owner*:** The Equipment Owner (Person field), if not the organization.
 - **[Equipment] System*:** System that this Equipment belongs to.

- **Barcodes***: Use this field to add additional barcodes to the record that may already be on the equipment. This field is useful when no additional bar-coded labels are needed and Calibration Control will recognize the barcode labels already affixed to the equipment.
 - **Equip Cost**: The purchase price or calibration cost of the Equipment.
 - **Acquired**: Acquisition date of the equipment.
 - **In Service**: In-Service date of the equipment.
 - **Cal Cost**: The standard (normal) cost of calibrating this Equipment.
 - **Time**: The standard amount of time required to calibrate this Equipment.
- ♦ **Check Out Group**: Using the Check In/Out dialog automatically populates the following fields:
 - **By***: The Checked Out By field shows the person who has checked out the Equipment.
 - **Date**: The Checked Out Date field shows the most recent date this Equipment was checked out, and only populates if it is checked out. When not checked out, this field is blank.
 - **Use Count**: Number of times this Equipment was checked out since its last Calibration.
 - ♦ The **Uncertainty Studies** panel grid allows you to add or link and Uncertainty group to the equipment record.

The screenshot shows the 'Edit Equipment' dialog box with the following fields and sections:

- Miscellaneous** section:
 - Owner: James Pinkerton
 - System: [Dropdown]
 - Barcodes: 00024596
 - Equip Cost: 63.00
 - Calibration Cost: 25.00
 - Acquired: [Dropdown]
 - In Service: [Dropdown]
 - Procedural Group: [Dropdown]
 - Personal Property
- Check Out** section:
 - By: Laura Welch
 - Date: 09/02/2019
 - Use Count: 1
- [Uncertainty Studies] - 1 Record** table:

Study ID	Date	Descrip	Nomi	Units
▶ 1	07/30/2			
- Buttons** at the bottom: Cancel, Cal Label, Worksheet, Asset Tag, Submit.

Other Tabs in the Equipment Dialog

There are some extra tabs within the Equipment dialog that help organize data related to Equipment records.

The Custom tab contains extra fields available in text, date, number, and checkbox format to easily rename and use for whatever needed.


The Image tab allows you to select one main image that represents that specific Equipment record.

The Overrides tab, while default labels and reports are defined at a global (application) level in the Options dialog, overrides for the Labels and Reports listed here can be entered for this specific Equipment.

The Meta tab is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.


Calibration Events

In the Equipment dialog, the Calibration summary is displayed from the Calibrations tab. To create or add a new Calibration Event record, click the green [+] icon, or right-click in the Calibrations records gray area to select "New Record" from the context menu.



Create a Simple Calibration Event

The following image shows a simple Calibration Event without any Measurements, as is often the case when recording calibrations handled by an external calibration laboratory.



The Calibration Results (As Found and As Left) can be selected when no Measurement Groups exist. The Technician, Temperature, and Humidity fields are usually not required in the outside lab scenario because this data would be included on the certificate the external organization provided. The calibration certificate provided by the outside calibration company can then be attached to this Calibration Event.


Calibration Event Fields

- **Certificate:** A sequential number is assigned to each Calibration Event. This number is automatically used as the Certificate number, unless another value is entered in this field, and it should be unique. To set the number value to automatically increment for this field, see the Number Masks topic for instructions.
- **Date:** Date the Calibration Event took place. This field defaults to today's date.
- **Company:** The organization that calibrated the equipment, like your own company (Internal) or a calibration laboratory (External).
- **Temp (Temperature):** Enter the ambient temperature of the room where the calibration is performed. If a THUM device (USB hygrometer) is plugged in and configured this field will auto populate.
- **Humidity:** Enter the Relative Humidity of the room where the calibration is performed. This field is also autopopulated with a configured THUM device.
- **Technician:** The individual who performed the calibration. This field is auto populated from the Technician field on the parent Equipment record.
- **As Found:** Record the condition of the equipment when it was received for calibration (e.g., Out of Tolerance or In Tolerance). This field gets auto-populated when the results are tied to the results of the Measurement Groups.
- **As Left:** Record the condition of the equipment at the end of the Calibration Event. This field is auto-populated when using Measurement Groups because the results of the groups are tied to the overall result of the CalibrationEvent.
- **Remarks:** Use this field for any remarks related to the Calibration Event.
- Refer to **Measurement Groups** section.

Attachments | Documents | Notes | Parts Tab

This panel grid contains the Attachments, Documents, All Notes, and Part Numbers grids. Use the Add, Link, Edit, Unlink, and Delete buttons at the bottom of each panel grid to edit the information. Resize the panel grids by clicking and dragging the spacer bars.


- **Attachments:** Attach a certificate PDF provided by an external calibration company, a cal worksheet, or link other existing Attachments to this calibration record.
- **Documents:** Attach your Procedures to this calibration event, if it is not already linked to the Equipment record.
- **Notes:** Link other existing Note records to this calibration or add a new Note.
- **Part Numbers:** Link more information from Part Number records to this calibration.



Standards | Misc. Tab

Calibration Standards can be tied to a calibration event in *three* tiers of traceability for the test/master gages used. You can select Standards at the Calibration record level, the Measurement Group level, and/or at the Nominal level.

This tab allows you to identify the Standards used in the overall Calibration event level.



Miscellaneous information related to the Cal Event is also available in this tab:

- **Use Count:** Shows how many times the equipment being calibrated has been used since its last calibration.
- **Approval:** Checking the box in this tab approves and locks the current Calibration Event record from future editing. If the current user is signed-in at the time the box is checked, that Person's name will be tied to the

approval of the Event and the name and date is signed on the Calibration Certificate stock report.

- **Cal Cost:** The actual cost of this Calibration Event.
- **Cal Time:** The actual time required to perform this Calibration Event.
- **Last Due Date:** This is the most recent Equipment Due for Calibration date, before this current event.
- **Received for Calibration:** Choose the date the equipment was received for calibration. This field auto-populates with the Equipment's Received for Calibration date *if* your setting is defined in Options with the Status Change feature.

Calibration Due for Certificate


At the bottom of the Standards/Miscellaneous tab are settings for the Calibration Certificate.

- **Calibration Frequency:** Choose the frequency of calibration, if this field does not auto-populate from the Equipment record. This displays on the Calibration Certificate stock report.
- **Units:** Choose the appropriate Frequency Units, if this field does not auto-populate from the Equipment record. This displays on the Calibration Certificate stock report.
- **Next Cal:** Select the date of the equipment's next calibration, if this field does not auto-populate from the Equipment record. This displays on the Calibration Certificate stock report.

Measurement Groups: Structure and Fields

- Each Calibration Event can have several Measurement Groups or none at all.
- Each Measurement Group can have any number of Nominal Measurements in it.
- Each Measurement Group can have its own Tolerance, Unit of Measure, and Equipment Standards Used in Calibration.

To illustrate, the diagram below is a single Calibration Event with three Measurement Groups with a different number of nominals in each Measurement Group.



Measurement Group Fields

In the Edit Calibration dialog, under the Calibration tab, use the bottom panel grid to Add [+] Measurement Groups to the Calibration Event. Adding a Measurement Group presents the following terms from the Measurements tab:


- **Group Number:** Measurement Group Number auto-populates; change the number, if needed.
- **Tolerance Type:** Set as None, Percent (%), Manual Entry, Tolerance (+/-). The following is the result of each type:
 - None:** The Upper and Lower Limits are the same as the Nominal.
 - Percent (%):** [Upper Limit = Nominal + (Nominal x Plus Tolerance)] and [Lower Limit = Nominal - (Nominal x Minus Tolerance)]
 - Manual Entry:** Manually enter the Tolerance values.
 - Tolerance (+/-):** [Upper Limit = Nominal + Plus Tolerance] and [Lower Limit = Nominal - Minus Tolerance]
 - Locked:** Once Tolerance and plus/minus values are set, select to Lock the Tolerance Type to avoid accidental editing of upper and lower limits.
- **Plus(+) and Minus(-):** Choose the upper and lower tolerance limits when using 'Percent' or 'Tolerance' as the Tolerance Type. Then set the Tolerance Type to "Locked" to avoid changes.
- **Mask:** The Number Mask is used to set the Integer Spaces (number places to left of decimal) and Significant Digits (number places to right of the decimal). The first mask field controls only the Nominal Mask. The second mask field controls the general mask for all values other than the nominal.
- **Group Name:** Enter a group name for the Measurement, if needed.
- **As Found:** Automatically In Tolerance (green) or Out of Tolerance (pink) dependent on the value of the As Found Result column in the Measurement Grid below (yellow area). If any record has a failing result (red X), the Measurement Group As Found value will be Out of Tolerance.

- **As Left:** Automatically Pass (green) or Fail (pink) dependent on the value of the As Left Result field in the Measurement Grid below (blue area). If any record has a failing result (red X), the Measurement Group As Left value will be Fail.
- **Units:** Select the Unit of Measurement for all number fields (e.g., measurements and tolerances) in the current Measurement Group. (When a different Measurement Unit is required for the same Calibration Event, create an additional Measurement Group in the Cal Event.) Edit all available Measurement Units in the combo box from the main Codes table.
- **Measurement Notes:** Add any notes specific to the current Measurement Group.

	Nominal	Limits		As Found			As Left		
		Upper	Lower	Found	Error	Re	Left	Error	Res
	0.0050	0.0055	0.0045	0.0050	0.0000	✓	0.0050	0.0000	✓
	0.0100	0.0105	0.0095	0.0100	0.0000	✓	0.0100	0.0000	✓
	0.0150	0.0155	0.0145	0.0150	0.0000	✓	0.0150	0.0000	✓
	0.0200	0.0205	0.0195	0.0200	0.0000	✓	0.0200	0.0000	✓
	0.0250	0.0255	0.0245	0.0250	0.0000	✓	0.0250	0.0000	✓

Additional Measurement Group (MG) Features:


- **Show Sequence:** Check the box to enter a custom sequence order for each measurement.
- **Show Description:** Check the box to enter a custom description for each measurement.
- **Auto Fill After:** Check the box to copy As Found values and auto-fill as the As Left field values, too.
- **Show Uncertainty:** Check the box to enter Measurement Uncertainty values to the results. Double-click in the new Other column to enter "Expanded Uncertainty" and/or select an existing Uncertainty Study from the combo box.
- **Show Standards:** (Standards Used for Nominals) Check the box to tie standards to each measurement. Double-click in the new Other column Stds field to Edit Standards and select from available records.
- **Sample Measurements:** Double-click in the As Found or As Left fields for any measurement nominal to enter range of sample measurements and use the average result as the As Found or As Left values.




Standards Used in Measurement Group

Linking Standards in the Measurement Group is the most popular user method of traceability for Calibration Events. (You can link Standards at the Calibration record level, each Measurement Group level, and/or at the specific Nominal levels.)

View the Standards Used tab to see a list of the Calibration Standards used for measurements in the current Measurement Group (MG). To link Standard Equipment in the MG, click on the [Edit Standards] button to bring up a list of available Standards to choose from.



Check the "Use" box(es) of each Standard(s) record, and click OK. Equipment records identified as a Calibration Standard, and have a Calibration Due date greater than the current day, AND with its most recent (if any) Calibration Event as 'Passing', can be selectable from the Standards Used button. Standards that need calibration are shown in red, but cannot be selected.



Additional Tabs in the Edit Calibration Event Dialog


There are three last tabs in the Edit Calibration dialog that help organize and display information related to the current Calibration Event.

- **Custom Fields:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Scatter Plot:** This tab displays a Scatter Plot of the measurements entered for a given Measurement Group. In-Tolerance As Found measurements are displayed Green while Out-of-Tolerance measurements are Red.
- **Meta:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.

Calibration Frequencies

How do Calibration Frequencies Work?

Calibration Frequency is the interval of time between required calibrations of any given piece of test equipment. Set the Calibration Frequency of Equipment in the Calibrations tab of the Equipment dialog (below).



Calibration Fields

There are four fields involved in the Calibration Frequency function. They are:

- **Last [Calibration]:** Date of the last calibration.
- **[Calibration] Frequency:** How often the Equipment is calibrated.
- **Units:** Calculates the frequency automatically based on the unit selected. The Unit field is disabled with Frequencies like Yearly or Monthly because they have an inherent Unit value.
- **Next [Calibration]:** Due date of the next required Calibration.

☰ Calibration Frequencies


The following are all the Calibration Frequencies in Calibration Control where 'n' is the number value in the Units field.

<u>Frequency Name</u>	<u>Frequency Description</u>
Biennial	Every 2 years
Bi-Monthly	Every 2 Months
Bi-Weekly	Every 2 Weeks
Cal Not Required	No Due Date . . . Not Calibrated

Daily	Every Day
Days	Every n Days
Manual	Set the Next Calibration (Due Date) manually
Month End	Every n Months due on the last day of the month
Month of	Every n Months due at any day within the calculated
monthMonth Start	Every n Months due on the first day of the month
Monthly	Every Month
Months	Every n Months
Next Use	No Due Date . . . Next Use
Not Calibrated	No Due Date . . . Not Calibrated
Quarterly	Every 3 Months
Reference Only	No Due Date . . . Reference
OnlySemi-Annual	Every 6 Months
Undefined	No Frequency Defined
Use Count	Calibrate every n uses where n = Units. Use the Check In/Out dialog to automatically increment uses.
Use Days	Calibrate every n days where n = Units. Use the Check In/Out dialog to automatically increment uses.
Week of	Every n Weeks due at any day within the calculated week
Weekly	Every Week
Weeks	Every n Weeks
Yearly	Every Year
Years	Every n Years

“Month of” and “Week of” Frequencies

Unlike the other Frequencies that calculate due dates on specific days, the "Month of" and "Week of" are special Frequencies that calculate the month or week a calibration is due. These Frequencies are also affected by two settings in the Options dialog for "First Day of Week" and "Due Date 'Week/Month of' Offset". By default, the "First Day of Week" is set to whatever is standard on your operating system. The default for the Offset is the Next Month or Week after the calculated due date.



The screenshot shows the "Options" dialog box with the "General" tab selected. The "Advanced" tab is also visible. The "General" tab contains the following settings:

- First Day of Week: System
- Due Date 'Week/Month of' Offset: Current
- Check for Program Updates: Daily
- Measurement System: Metric
- Equipment ID of Default Template: (empty text box)
- Unique Alternate Equipment ID: No
- Language: English (US)
- Label Printer: (empty text box)
- Start Screen: (empty text box)
- Do Not Copy Default Labels, Reports, & Emails on Start
- Use Active Equipment as Default for Asset Transfer
- Dialog Font Size (+/-): 1

Buttons for "Cancel" and "Save" are located at the bottom of the dialog.

"Month of" Example

- **Frequency:** Month of
- **Last Calibration:** 04/27/2019
- **Units:** 12
- **Offset:** Next
- **First Day of Week:** Monday
- **Calibration Due:** 04/27/2020
- **Label Display:** Due Month of Apr-2020

How the math works . . .

1. Calculate 12 Months after Last Calibration, which yields 04/27/2020
2. The Offset is set to Next [Month] so the Due Month jumps to the nextMonth, which yields 05/1/2020

"Week of" Example

- **Frequency:** Week of
- **Last Calibration:** 06/01/2019
- **Units:** 30
- **Offset:** Next
- **First Day of Week:** Monday
- **Calibration Due:** 07/01/2019
- **Label Display:** Due Week of 07/01/19

How the math works . . .

1. Calculate 4 Weeks after the Last Calibration, which results in 06/29/2019.
2. Since the Offset is set to Current [Week] and the First Day of Week is set to Monday, the Monday of the current week results in a due week starting on 07/01/2019.

Referencing Calibration Standards

How to Link Standards to Calibrated Equipment

Linking Standards/Masters to Calibration Events in Calibration Control (our **Calibration Management Software**) is simple. They create a data trail (traceability) between national measurement standards and the Equipment Under Calibration.

ID a Calibration Standard

Before referencing in a Calibration Event, the Standard must exist as an Equipment record with the Calibration Standard checkbox checked. Find this checkbox in the Calibrations tab of the Equipment dialog.

The screenshot shows the 'Edit Equipment' dialog box with the 'Calibrations' tab selected. The 'Calibration Standard' checkbox is checked and highlighted with a red box. Below the form is a table showing one calibration record for certificate 1178, dated 01/30/2019, with 'In Tolerance' and 'Pass' status.

Certificate	Date	As Found	As Left	Cal Company	Technician
1178	01/30/2019	In Tolerance	Pass		David Wilcox

Grid Highlight

By default, Calibration Standards display with a green highlight in the Equipment grid. Refer to the Highlight Colors helptopic for more information.

Equipment					
[Equipment] - 263 Records (8 Past Due)					
Equipment ID	Model	Description	Serial	Manu	
SAMPLE-001	BMS10-83	Type III	112		
SAMPLE-002	IDR-8111	Indicator	2213	East C	
SAMPLE-003	GBS0-5	Gage Block Set			
SAMPLE-004	CRG-101	Cylindrical Ring G	8411		
SAMPLE-005	CRG-101	Cylindrical Ring G	199		
SAMPLE-006		10-32 UNF-3B	1008		
SAMPLE-007	CRG-101	Cylindrical Ring G	195		
SAMPLE-008	CRG-101	Cylindrical Ring G	2002		

Link Calibration Standards

After Equipment records for Calibration Standards are created and identified as Standards, they can be linked to individual Calibration Events. Double-click on a Calibration Event within the Edit Equipment dialog and click on the Standards/Misc tab. From there, click on the [Edit Standards] button to select Standards. Refer to the Calibration help topic for more information on creating a Calibration Event.

Use Count: Check to Approve & Lock Last Due Date:

Cal Cost: Received for Calibration:

Cal Time: **Edit Standards**

Standards used During this Calibration Event							
Equipment ID	Model	Description	Serial	Status	Last Cal	Next Cal	Certificate
SAMPLE-010	PS-d80b5	Pin Set		Active	05/24/2018	05/24/2019	

Calibration Due for Certificate

Calibration Frequency: Units: Next Cal:

Available Standards

Clicking the [Standards] button in the Measurement Group of 'Calibration Details' displays a list of available Standards. Unavailable Standards are those (a) standards currently out of calibration, (b) standards with their last calibration event with an As Left result of "Failed" or "Unknown", or (c) is the current equipment under calibration.

Use	ID	Model	Description	Status	Last Cal	Due Cal
<input type="checkbox"/>	SAMPLE-001	BMS10-83	Type III	Active	11/23/2018	11/23/2020
<input type="checkbox"/>	SAMPLE-009	GPS-68	Gage Pin Set	Active	12/26/2018	12/26/2019
<input checked="" type="checkbox"/>	SAMPLE-010	PS-d80b5	Pin Set	Active	05/24/2018	05/24/2019
<input type="checkbox"/>	SAMPLE-011	PS-d80b5	Pin Set	Active	05/24/2018	05/24/2019
<input type="checkbox"/>	SAMPLE-012	PS-d80b5	Pin Set	Active	05/24/2018	05/24/2019
<input type="checkbox"/>	SAMPLE-014	GPS-68	Gage Pin Set	Active	10/16/2018	10/16/2019
<input type="checkbox"/>	SAMPLE-023	IDCM3	Indicating Micrometer	Active	05/10/2018	05/10/2019
<input type="checkbox"/>	SAMPLE-028	MGBS-878	Master Gage Block Set	Active	06/28/2018	06/28/2019

Results

The process of linking a Standard to a Calibration Event makes a copy of the calibration information for that Standard and stores it in the Calibration Event record. This means that when the calibration information for the Standard changes (after the Calibration Event), a snapshot of the Standard at the time of Calibration is preserved, which helps maintain accurate traceability.

Note: If the Due Cal date has passed the present Calibration date, then the Standard is red and unavailable because it is past due for Calibration.

Calibration Templates


How to Use Calibration Templates

A Calibration Template in Calibration Control (our **Calibration Management Software**) is a collection of Measurement Groups (without final measurement data) and other information such as calibration company, technician, remarks. Therefore, a Template can have multiple groups. A Measurement Group is a collection of information and standard measurements required each time the template is used. For more information on each of the fields within a Calibration Template or Measurement Group, see the Calibration Event help topic.

Default Calibration Templates


There is a setting in program Options: 'Use Last Calibration as Template if No Template Linked' and this box is checked by default. When this option is checked AND if there are no other user-defined Calibration Templates linked to an

Equipment record, Calibration Control will use the most recent Calibration Event record of that Equipment as a Calibration Template when creating new Calibration Events. Therefore, if any history of Calibrations exist in an Equipment record, the next Calibration Event automatically uses the previous Calibration Event's measurements as its template.



Create User-Defined Calibration Templates

Calibration Templates can be created, edited, and managed by clicking on the Calibration Templates icon in the Calibration tab of the ribbon menu.



Alternatively, a template can be created from any existing Calibration Event, by selecting the 'Make Template' button under

the Calibrations tab of the Calibration Event dialog. That way, you can easily create a copy of an existing template by opening the template you'd like to copy and simply click the [Make Template] button. Then, enter a unique template name and Submit. Refresh the Calibration Templates Grid (right-click and click Refresh Grid or press F5). This will show the full template copy as a new template record.

The screenshot shows the 'Edit Calibration CC-0329' dialog box. The 'Make Template' button is highlighted with a red box. Below the form is a table with 3 records for measurement groups.

Group	Units	Tolerance Type	Plus Tolerance	Minus Tolerance	As Found	As Left	Measurements
1	mm ³	Tolerance (+/-)	0.01	0.01	In Tolerance	Pass	3
2	mm ³	Percent (%)	5	5	In Tolerance	Pass	1
3	mm ³	Tolerance (+/-)	0.0001	0.0001	In Tolerance	Pass	4

Calibration Template Structure and Hierarchy


A Calibration Template needs to be linked to an Equipment record to be useful. When linked, new Calibration Events will automatically contain the structure of the linked Calibration Template so that only Measurement data (As Found / Left results) will need to be entered and Calibration Standard(s) will need to be linked for the record. Note: Calibration Worksheets also replicate the measurement structure of user-defined or default templates for consistent data collection.

Link Templates (Hierarchy)

Multiple methods for links exist to provide the choice of how an organization chooses to manage its calibration data. Therefore, choose the method that works best for your team. Referring to the diagram below, a Calibration Template can be linked to an Equipment record in four ways. It is even possible to link four different Templates to an Equipment record although only one template will be used and in the following priority:

1. **Direct (blue line):** Link a Template directly to an Equipment record by selecting the Template in the Calibrations tab of the Equipment record. A Template directly connected to an Equipment record in this way supersedes all other methods. Therefore, this direct link overrides all other indirect links.
2. **Equipment Type (red line):** Link a Template to an Equipment Type record and then link the Type to an Equipment record in the Equipment tab of the Equipment's dialog. This link overrides the two that follow.
3. **Model Number (green line):** Link a Template to a Model Number record and then link the Model Number record to the Equipment record in the Equipment tab of the Equipment's dialog.


4. **Equipment Type – Model Number (yellow line):** Link a Template to an Equipment Type, which is linked to a Model Number, and then is finally linked to an Equipment record. This is the lowest priority.



Attachments


Add Attachments to Equipment and Calibration Records

Add attachments to an Equipment record by double-clicking on a record and navigating to the 'Attachments | Documents | Notes | Parts' tab, where you will find the [Attachments] panel grid. Use this panel grid like you would any other grid: Add, Link, Edit, Unlink, or Delete attachments. To add a new attachment click on the green [+] button and choose whether to add a File, Link to File, add a Folder path, or a URL path.



Add File, Link to File, Folder or URL


[Add File] means a copy of the attachment will be made and saved in the Attachments folder (located in the Files Folder). The [Link to File] selection means the file being attached will stay at its current location and a link to that location will be saved in the database. Adding a link to a specific folder is helpful if more than one file is being attached. The [Add Folder] and [Add URL] buttons allow a user to copy and paste a path to a specific folder (directory) or website.



A dialog box titled "Paste Folder Path" with a close button (X) in the top right corner. It contains a text input field with the text "C:\Users\Ape Software\Pictures\Calibration Equipment". To the right of the input field are two buttons: "OK" and "Cancel".

Edit Attachment

Within an Equipment record or in the Attachments grid, you can double-click on an attachment to edit the record and add information related to it such as title, category, owner (if the image is of a piece of equipment), and a description.




An "Edit Attachment" dialog box with a title bar containing a question mark and a close button (X). The dialog has three tabs: "Attachment" (selected), "Equipment | Calibrations | Companies | People", and "Meta". The "Attachment" tab contains the following fields:

- Title: flashcal.jpg
- Type: File (dropdown menu)
- Attachment: flashcal.jpg (text input field with a pencil icon and ellipsis)
- Category: Image (dropdown menu)
- Link Only:
- Active:
- Owner: Jose Mendoza (dropdown menu)
- Present:
- Date: 09/10/2019 (dropdown menu)
- Description: Image description goes here. (text area)


At the bottom of the dialog are two buttons: "Cancel" and "Submit".

In the Equipment | Calibrations | Companies | People tab, you can view all records linked to the attachment. You can also add, link, and unlink records to the attachment. Similarly, in the Models tab, you can add, link, and unlink Model records to the attachment.



Attachments Grid

The Attachments grid, located in the Data Grids tab of Calibration Control, displays all attachment records information such as file name, size, title, description, etc.



Right click anywhere on the grid to open the context menu and click **Verify All Attachments Exist**. If the Attachments were verified, a green checkmark will appear under the Present column along with the Confirmed Present date. If an attachment was not verified, it likely was moved to a different location, renamed, or is otherwise not accessible to Calibration Control. If *all* attachments were not confirmed present, your Files Folder was likely changed but your Attachments folder was not moved to the new location. When an attachment is verified to be present, you can quickly open the attachment file itself by right clicking a record and select Open Attachment.

Additionally, in the Attachments context menu, you can select different views of the Attachments such as all attachments, attachments linked to Equipment, Calibrations, People, Companies, or attachments not linked to anything.




Image in Equipment Dialog


To add a default image for an equipment record, click on the [Image] tab within a selected Edit Equipment dialog. If no image is linked, you'll see a message "Right-click, Double-click or Drag & Drop to Link an Image." Double clicking will open a grid of existing attachments that are image files. Alternatively, you can drag and drop any image file and Calibration Control will create an Attachment record for that image and link to the Equipment record. Once you have an image linked, you can right-click to open the image, unlink the existing image or link a different an existing image.



Companies


Manage Companies within Calibration Control

Keep track of all the companies you do business with, link them to other records, and manage report overrides using the Company grid and dialog. Find the Companies grid in the Common tab of the menu ribbon of Calibration Control (our **Calibration Management Software**).



Name	Status	Address	City	State
3M	Active	554 3M Circle	Milwaukee	WI
ACME Label Print	Active			
Acme Tools	Approved	996 West Canes	Des Moines	IA
Aerospace, LLC	Active	11342 6th St	Harrisburg	IL
Amazon	Active			
Anderson Windo	Active	3600 Mayo Street	Cincinnati	KY
ARC Manufacturi	Active	721 Wonderful W	Marysville	OR

The grid header shows a count of your number of Company records (or how many filtered records). Right-click in the Companies grid for an alternate grid view to include the Company Types from the context menu. Note the grid header will show a higher count if any records have more than one linked Company Type.



Name	Type
3M	Customer/Client
ACME Label Printers	Supplier
Acme Tools	Customer/Client
Aerospace, LLC	Customer/Client
Amazon	Supplier
Anderson Windows and Doors	Customer/Client
ARC Manufacturing	Prospect
ARC Manufacturing	Manufacturer
B&H Photo	Supplier
Beauty Manufacturing	Manufacturer
Brother	Manufacturer
Calibration 77	Calibration

Company Dialog

Use the Main tab to organize all information related to the company itself.

- **Name:** The name of the Company. This is the only required field for a new Company record.
- **ID Number:** Assign an identification number to a company record.
- **Web:** The company's website URL.
- **Status:** Select the status of the company (Active, Approved, Inactive, etc).
- **Next Onsite:** Select the date that a technician will return on site again.
- **Do Not Auto Email:** Select this checkbox to exclude this company from Auto Email Notifications.
- **Primary and Alternate Addresses:** Enter mailing, billing, or physical address(es).
- **Notes Field:** Enter short notes or remarks related to this Company record.

The screenshot shows a software dialog box titled "Add Company Ape Software (Load Time: 0.0 Seconds)". It features a tabbed interface with the "Main" tab selected. The "Main" tab contains the following fields and controls:

- Name:** Text input field containing "Ape Software".
- ID Number:** Text input field containing "APE".
- Web:** Text input field containing "www.apesoftware.com".
- Status:** Dropdown menu set to "Active".
- Next Onsite:** Date selection dropdown.
- Do Not Auto Email:** Unchecked checkbox.
- Primary Address:** A group of five text input fields for Address, City, State, Postal, and Country.
- Alternate Address:** A group of five text input fields for Address, City, State, Postal, and Country.
- Notes:** A large text area for entering remarks.
- Buttons:** "Cancel" and "Submit" buttons at the bottom.


Email, Phones, People, and Types Dialog

Use these panel grids to keep track of the email addresses, phone numbers, people, and company types associated with the Company record.

- **Email Addresses:** Click on the chain link icon to link an existing email address to a company record or click on the green [+] icon to create a new record. Within the email address dialog, select the Type (i.e., Personal, Work, Other) and choose whether the email is Primary, Opted Out, or Invalid. You are also able to easily link or add Company and People records to the email address record right from the dialog.
- **Phone Numbers:** Add and link phone and fax numbers to a Company record. Within the phone numbers

dialog,select the Type (i.e., Direct, Fax, Office) and choose whether the number is Primary, Opted Out, or Invalid. You are also able to easily link or add Company and People records to the phone number record right from the dialog.


- **People:** Use this tab to link a Person record to a Company. (Use the chain link and green [+] icons, just like the other panel grids.) You can designate someone as the 'Due Cal Contact' at this company for auto-emails. To select a person as the 'Due Cal Contact', open their Person record and within the Main tab check the 'Company Due CalContact' checkbox.
- **Company Types:** Link a Company Type using the 6 checkboxes to specify if it is a Site, Manufacturer, Calibration Service or Maintenance Provider, Client, Supplier, or Other.



Notes and Attachments Tab

Keep track of all notes, certificates, file documents, and/or pictures related to a Company record by using these panel grids.

- **Notes:** Create new or link existing Notes associated with the Company record.
- **Attachments:** Add new or link existing Attachments to a Company record, such as documents or images.



Other Tabs

- **Overrides Tab:** While default labels and reports are defined at a global (application) level in the Options dialog, overrides for the Labels and Reports listed here can be entered for a specific Company.
- **Custom Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these changes occurred.

The screenshot shows a software window titled "Edit Company" with a ribbon menu at the top containing "Main", "Email | Phones | People | Types", "Notes | Attachments", "Overrides", "Custom", and "Meta". The "Overrides" tab is active, displaying two sections:

- Equipment Label Overrides:**
 - Calibration Due: [input field] [X]
 - Asset: [input field] [X]
 - Cal Not Required: [input field] [X]
 - Next Use: [input field] [X]
- Equipment Report Overrides:**
 - Calibration Cert: [input field] [X]
 - Cal Worksheet: [input field] [X]
 - Blank Worksheet: [input field] [X]
 - OOT Form: [input field] [X]

At the bottom of the window are "Cancel" and "Submit" buttons.

Company Types

Edit Company Type Codes

Although any number of Company Type codes can be created, there are six types important to Calibration Control (our **Calibration Management Software**). These six Company Types are: **Site, Manufacturer, Calibration Service (or Maintenance Provider), Client, Supplier, or Other**. These six types will identify *which categories* that a Company record of that Type can show up as combo-box choices in a drop-down selection for other records, (e.g., as a Calibration Company for a Calibration Event, etc.)

Display this grid by clicking on the Company Types icon in the Data Grids tab of the ribbon menu.

Type Name	Site	Manufacturer	Calibration	Supplier
Calibration	<input type="radio"/>	<input type="radio"/>	<input checked="" type="checkbox"/>	<input type="radio"/>
Customer/Client	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Manufacturer	<input type="radio"/>	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>
Other	<input type="radio"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="radio"/>
Prospect	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Site	<input checked="" type="checkbox"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Supplier	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="checkbox"/>
Unknown	<input type="radio"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Company Types Dialog

Create a Company Type name and mark check-boxes to specify if that Type is a Site, Manufacturer, Calibration Service or Maintenance Provider, Client, Supplier, or Other.


Edit Company Type

Type:

Site Client
 Manufacturer Supplier
 Calibration / Maintenance Other

Link to Company Records

Under the Companies tab of the Edit Company Types dialog is a panel grid where you can add, link, edit, or unlink Companies that fall under a specific Company Type.



Other Tabs

- **Custom Fields:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred. This tab is populated automatically for security purposes and is not editable.

Equipment Systems

Groups of Equipment Records

Combine separate equipment records in Calibration Control (our **Calibration Management Software**) into a group (system) for any purpose, such as creating a single test system.

View Systems Grid

Select the Systems feature from the Calibration tab of the menu ribbon.



Equipment System Dialog

From the Systems grid, double-click to open a record or right-click to use the context menu. Use the Equipment System dialog to enter and edit information related to Equipment Systems.

The screenshot shows a software dialog box titled "Edit Equipment System Tool Kit (Load Time: 0.2 Seconds)". The dialog is organized into five tabs: "System", "Equipment", "Notes", "Custom Fields", and "Meta". The "System" tab is currently selected and contains the following fields:


- System Name:** Tool Kit
- System Code:** KOT-1
- Site:** (empty dropdown)
- Department:** Machine Shop
- Location:** Locker
- Status:** Active
- Custodian:** Chris Shaw
- Technician:** Gary Anderson
- Owner:** (empty dropdown)
- Checked Out By:** (empty dropdown)
- Checked Out:** (empty dropdown)
- Certificate:** 559-12

At the bottom of the dialog, there is a large text area for "Notes" and four buttons: "Cancel", "Print Asset Label", "Update Equipment Records", and "Submit".

- **System Name:** Name of the System. This is the only required field for a new Systems record.
- **System Code:** Shorter code name of the System
- **Site:** Company field for geographic Site or Company.
- **Department:** The department in your organization where the System is currently located.
- **Location:** The location within the department where the System can be found. This could be a specific area or work place, even an engineer's desk.
- **Status:** Select the system's current Status from the combo box.
- **Custodian:** Select an individual who has possession or custody of a System.
- **Technician:** The person who normally calibrates this System.
- **Owner:** The System Owner (Person field), if not the organization.
- **Checked Out By:** This field will auto-populate if a Person was required or identified during check out.
- **Checked Out:** This field will auto-populate with the date if the System is currently checked out.
- **Certificate:** Certificate number that was issued by the organization responsible for its last Calibration event. Primarily used as an aid for traceability back to a national measurement standard, especially if this System is a Calibration/Test Standard.
- **Notes:** Add short descriptions to an Equipment System record.

Equipment Tab


To assign Equipment records to an Equipment System, select the Equipment tab in the System dialog. Click the link button to choose which Equipment records belong to this System.



(An alternative method to link Equipment records to an Equipment System is shown in the next section below.)


Notes Tab

Create a new Note record or link existing Notes related to this Equipment System.



Custom Fields Tab

Customizable fields are available to rename and utilize Characters (text fields), Dates, Numbers, and Check boxes (booleans). Visit the Change Field Names help topic to learn how to edit these field names.



Meta Tab

This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.

Mass Update Function

Update Equipment Records: Edit System fields and apply those changes to all related Equipment records tied to that

System with the [Update Equipment Records] button of the System dialog.

The screenshot shows a software dialog box titled "Edit Equipment System Tool Kit (Load Time: 0.1 Seconds)". It features five tabs: "System", "Equipment", "Notes", "Custom Fields", and "Meta". The "System" tab is selected, displaying fields for System Name (Tool Kit), System Code (KOT-1), Site, Department (Machine Shop), Location (Locker), Status (Active), and Notes. The "Equipment" tab shows fields for Custodian (Chare Shaw), Technician (Gerry Anderson), Owner, Checked Out By, Checked Out, and Certificate (659-12). At the bottom, there are buttons for "Cancel", "Print Asset Label", "Update Equipment Records", and "Submit".

From the pop-up dialog, confirm which fields should update for all linked Equipment. All fields are checkmarked by default. Blank fields are empty values that will also update (and clear those fields in) the related Equipment records.

The screenshot shows a software dialog box titled "Update Equipment Records". It contains the text "Check fields to update in related Equipment records." Below this, there are fields for Site, Department (Machine Shop), Location (Locker), Custodian (Chare Shaw), Technician (Gerry Anderson), Owner, Status (Active), and Certificate (659-12). Each field has a checkbox to its right, all of which are checked. A small pop-up dialog titled "Records Updated" is overlaid on the main dialog, showing a blue circle with an "i" and the text "1 Equipment Records Updated", with an "OK" button. At the bottom of the main dialog are "Cancel" and "Update" buttons.

Link System to an Equipment Record

An alternative method to assign an open Equipment record to an existing Equipment System is in the Misc tab of the Equipment dialog. All Equipment Systems will be listed in the combo box to choose a System Name to which the Equipment belongs.

The image shows a software window titled "Edit Equipment" with several tabs: Equipment, Calibrations, Maintenance, Custom, and Attachments. The "Equipment" tab is active, showing a "Miscellaneous" section with the following fields:

- Owner: Walter Volpe
- System: Chamber 001** (highlighted with a red box)
- Barcodes: 0051489
- Equip Cost: 340.00
- Acquired: 01/23/2019
- In Service: 02/13/2019
- Procedure: (empty)
- Personal Property

At the bottom of the window, there are buttons for "Cancel", navigation arrows (|< and <), "Cal Label", and "Works".

Print Asset Label

Configure default System Asset Labels in program Options > Report & Label.

Print System asset labels with barcodes to Use Scanner in the Check Out feature dialog.

Check Out an Equipment System

Enter a valid Equipment System Name in the Check Out dialog to check in and out the group collection of equipment at onetime. The 'Equipment ID' field in the Check Out dialog reflects 'System ID' once a System is validated.

Check Out

Check in or Out?
 Out In Use Scanner

Increment: 0

System ID: Tool Kit

Person: Alberta Mayfield

Reference:

Change Status To:

The Tool Kit System consists of 8 Equipment.

Finished * Required if Checked Submit

The Equipment System record will update to show the recent check out date and/or person, and all related Equipment records linked to that System will update as well.

Equipment Types

Classify Equipment and Link to Calibration Templates

Although the Equipment Type in Calibration Control can be used for any classification, its intent is specifically related to ease of equipment selection and the management of Calibration Templates. For instance, even though generic 6 inch calipers are made by several manufacturers, their use and calibration is usually the same. Therefore, creating an Equipment Type of "6 Inch Caliper" can be a useful grouping.

Additionally, a single Calibration Template (calibration method) for 6 inch calipers can be linked to one Equipment Type record, which in turn is linked to all the Equipment records of the same type. Using Equipment Types and Calibration Templates together like this provides calibration management with the ability to define and apply more consistent calibration controls.

Equipment Types Grid

Display the Equipment Types grid by clicking on the 'Equipment Types' icon in the Data Grids tab of the ribbon menu.

[Equipment Types] - 5 Records					
Type Name	Type Code	Cal Template	Size/Range		P
Caliper	CP		0-6"		
Digital Indicator	DI		0-1"		
Digital Pressure	DPG		0-100 PSI		
Micrometer	MC	TEMP-003	0-6"		
Test Guage	TG				

Equipment Types Dialog

Double-click on a record to open and edit or double-click in the gray space to open a new record dialog box. Right-click inside the grid for the context menu.

- **Type Name:** The name of the Equipment Type or group of equipment. This is a required field.
- **Type Code:** The unique code for the Equipment Type or group of equipment. This is a required field.
- **Size/Range:** Select the Size and Range of the equipment.
- **Calibration Template:** Though not required, Equipment Types can be linked to a corresponding Calibration Template. This may influence the naming convention for the Equipment Type.
- **Procedural Group:** Allows users to define step-by-step procedures that appear on Calibration Worksheets or within Custom Reports.

Add Equipment Type

Type	Equipment	Uncertainty	Overrides	Meta
Type Name:	Micrometer			
Type Code:	MC	Calibration Template:	CT-004	
Size/Range:	0-6"	Procedural Group:	PR-057	
Notes:	Equipment Type notes go here.			

Buttons: Cancel, Submit

Other Tabs

There are four other tabs within this dialog to help manage the Equipment Type records.

- **Equipment Tab:** Use this panel grid to link equipment to an equipment type record.
- **Uncertainty Tab:** Use this panel grid to add or link an Uncertainty Study to an Equipment Type record.
- **Overrides Tab:** While default labels and reports are defined at a global (application) level in the Options dialog, overrides for the Labels and Reports listed here can be entered for a specific Equipment Type.
- **Meta Tab:** Shows a READ-ONLY log of who created the record and when, as well as who updated the record and when. This tab is populated automatically for security purposes and is not editable.

[Equipment] - 4 Records				
Equipment ID	Model	Description	Cal Last	Cal Due
SAMPLE-005	CRG-101	Cylindrical Ring Gage	04/13/2019	04/13/2020
SAMPLE-059	HG66	Height Gage	03/12/2019	03/12/2020
SAMPLE-283	ODM-32321	OD Micrometer	10/12/2018	10/12/2019
SAMPLE-143	SM-4556	Scissor Micrometer	11/03/2018	11/03/2019

Size and Range in Calibration Control

Edit the Size and Range Combo Values of Equipment

Add, edit, or delete the selectable values for the "Size/Range" field for Equipment in Calibration Control. Using these standardized values for size and range provide for more consistent data entry and Equipment searches.

Menu

Find the Size and Range values within the Codes grid, which is located in the Data Grids tab of the ribbon menu.

All Codes

[All Codes] - 25 Records

Code	Active
.0005 +/-	<input checked="" type="checkbox"/>
.0005 +/-	<input checked="" type="checkbox"/>
.0005 +/-	<input checked="" type="checkbox"/>
0-1"	<input checked="" type="checkbox"/>
0-100 PSI	<input checked="" type="checkbox"/>
0-12"	<input checked="" type="checkbox"/>
0-18"	<input checked="" type="checkbox"/>
0-6"	<input checked="" type="checkbox"/>
0-8"	<input checked="" type="checkbox"/>

Group : Sizes & Ranges (9 items)

Group : Status Codes - Company (5 items)

Group : Status Codes - Equipment (12 items)

Dialog

Use the Size/Range dialog to edit all information related to the different sizes and ranges of equipment that will appear under the Size/Range combo box in Calibration Control Equipment Records. This dialog includes custom fields and Meta data.

Add Size/Range

Code Custom Meta

Code: 0-1"

Code Short: 0-1 inches

Group: Sizes & Ranges


Sequence: 4 Is Active

Cancel Submit

Model Numbers

Keep Track of Equipment Model Numbers


Model Numbers are specific to each piece of equipment and the Model Numbers grid keeps them organized within Calibration Control (our **Calibration Management Software**). Find this grid in the Data Grids tab of the ribbon menu. The Model Numbers' dialog can also be accessed within the "Edit Equipment" dialog by clicking on the [Model] button in the Equipment tab.



Model Number Dialog

Double click on a selected record to edit it or double-click within the gray area to create a new record.

Within the dialog, specify the model number, model description, and notes. Add information like Manufacturer, Equipment Type, Measurement Template, Procedure Group, and its Size/Range using the drop down menus. Check or uncheck the [Active] box to show if it is in use or not.



Uncertainty Tab

Add or link an uncertainty study to a model number record.

Study ID	Date	Description	Nominal	Units
1	09/16/2019	Sample Uncertainty Study	0.1	cm

Overrides Tab


Use the Overrides tab to choose a different equipment label or equipment report. Click on a blue hyperlink and select the desired file name then click [Open]. Remove an override by clicking on its corresponding red [X]. Examples of all the labels can be found in the "Labels" section of our Help Topics page.

Equipment Label Overrides	Input Field	Remove (X)
Calibration Due:		X
Asset:		X
Cal Not Required:		X
Next Use:		X

Equipment Report Overrides	Input Field	Remove (X)
Calibration Cert:		X
Cal Worksheet:		X
Blank Worksheet:		X
Custom Report:		X
OOT Form:		X

Custom Tab

The Custom tab allows you to create unique fields that pertain to your company. Edit these using the "Edit Form Labels" feature under the Admin tab of the Program Options dialog. Select the dialog to be edited (i.e., "ModelNumber") from the drop down menu and create your custom content by editing the Control Text. Make sure to click the [Close] button in order to save changes.




The screenshot shows the 'Edit Model' dialog box with the 'Custom' tab selected. The dialog has five tabs: Model, Uncertainty, Overrides, Custom, and Meta. The Custom tab contains the following fields:

- Character 0: [Text Input]
- Character 1: [Text Input]
- Date 0: [Date Picker] (06/06/2019) | Number 0: [Text Input] (0.0000) | Check 0:
- Date 1: [Date Picker] (08/29/2019) | Number 1: [Text Input] (0.0000) | Check 1:

At the bottom, there are 'Cancel' and 'Submit' buttons.

Meta Tab

This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.




The screenshot shows the 'Edit Model' dialog box with the 'Meta' tab selected. The dialog has five tabs: Model, Uncertainty, Overrides, Custom, and Meta. The Meta tab contains the following fields:

- Created: [Text Input] (5/29/2019 3:52:50 PM)
- Created By: [Text Input]
- Updated: [Text Input] (6/6/2019 4:10:02 PM)
- Updated By: [Text Input] (BREANNA/Ape Software)

A red note at the bottom reads: "Note: This page is read-only." At the bottom, there are 'Cancel' and 'Submit' buttons.

Departments

Use the Departments grid and dialog to edit data related to Departments that appear in Calibration Control (our **Calibration Management Software**) Equipment Records. Display the Departments grid by clicking on its icon in the Data Grids tab of the ribbon menu.



Edit Department Dialog

Manage information related to a Department record:

- **Department:** The name of the department. This is a required field.
- **Code:** A short identifier or acronym for the Department. This is a required field.
- **Type:** The department type (i.e., Fabrication, Engineering, etc.).
- **Manager:** The Department Manager; Person to whom Calibration Control can send emails.
- **Contact:** The Department Contact; Person to whom Calibration Control can send emails.
- **Is Active Checkbox:** If checked, this shows that the Department Code is actively being used in the database.
- **Notes:** A notes field for the department record.

The screenshot shows a dialog box titled "Edit Department" with four tabs: "Department", "Notes", "Custom", and "Meta". The "Department" tab is selected and contains the following fields:

- Department: Engineering
- Code: ENG
- Type: Engineering
- Manager: Victor James
- Contact: Penelope Penn
- Is Active:
- Notes: Sample notes.

Buttons for "Cancel" and "Submit" are located at the bottom of the dialog.

Other Tabs

Within the Departments dialog there are three other tabs besides the main Department tab. The first is for *Notes* related to the Department dialog that can be added or linked from other records.

The *Custom* tab contains fields that can be customized to add Characters, Dates, Numbers, and Check boxes. Visit the [Change Field Names](#) help topic to learn how to edit these fields.

The *Meta* tab shows a READ-ONLY log of who created the Department record and when, as well as who updated the record and when. This tab is populated automatically for security purposes and is not editable.

Part Categories

Create Part Category Codes to Assign to Parts

Open the Part Category grid by clicking on its icon in the Data Grids tab of the ribbon menu. Double-click in the gray area to create a new record or right-click and select 'New Record'. Type out a short category code (usually an acronym) and then a short description of the part. Check the Active checkbox if this part category is active, and check the Taxable checkbox if this part is taxable.

The screenshot shows a dialog box titled "Edit Part Category Code" with three tabs: "General", "Custom", and "Meta". The "General" tab is selected and contains the following fields:

- Category Code: LP
- Description: Label Printers
- Active:
- Taxable:

Buttons for "Cancel" and "Submit" are located at the bottom of the dialog.

Other Tabs

- **Custom Fields:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.

Assign a Part Category to Part Number Record

Within the General tab of the Edit Part Number dialog, select the Part Category from the 'Category' drop-down menu.

The screenshot shows the 'Edit Part Number' dialog box with the 'General' tab selected. The 'Category' dropdown menu is highlighted with a red box and set to 'Label Printers'. Other fields include Part Number (PT_9800PCN), Part Name (P-Touch 9800), Description (Desktop Bar Code Network Printer), Standard Cost (389.00), List Price (480.00), and checkboxes for Taxable and Active.

Part Numbers

Manage Part Numbers Related to Equipment

Part Numbers are specific to each piece of equipment and the Part Numbers grid keeps them organized within CalibrationControl (our **Calibration Management Software**). Find this grid in the Data Grids tab of the ribbon menu.

The screenshot shows the software interface with the 'Data Grids' tab selected in the ribbon. The 'Part Numbers' grid is visible, displaying a list of part numbers with columns for Part Number, Part Name, Description, Std Cost, and List.

Part Number	Part Name	Description	Std Cost	List
Warranty 9800		Service & Repair	0	48
PT_9800PCN		Desktop Bar Cod	389	480
PT_P700		PC Connectable L	62	110
PT_P750W		Wireless Label M	106	180
PT_1230PC		PC-Connectable	30	55
Warranty 9700		Server & Repair	0	40

Edit Part Number Dialog

Use the General tab to organize all information related to part numbers.

- **Part Number:** Unique number assigned to the part record.
- **Part Name:** The name of the part record.
- **Description:** A physical description of the part.
- **Category:** The category that a part is assigned to.
- **Standard Cost:** The standard cost of the part.
- **List Price:** The listing price of the part.
- **Discontinued:** If the part has been discontinued, select the appropriate date.
- **Taxable and Active Checkboxes:** Check the corresponding box if a part is taxable or active. Leave the boxes unchecked if it is not.
- **Notes:** Add any notes related to the part number.

The screenshot shows a dialog box titled "Edit Part Number" with a standard Windows window header (minimize, maximize, close) and a help icon. The dialog has six tabs: "General", "Sources", "Equipment | Calibrations", "Work Orders", "Custom", and "Meta". The "General" tab is active. The fields are as follows:

Part Number:	<input type="text" value="PN-003"/>	Standard Cost:	<input type="text" value="62.00"/>
Part Name:	<input type="text" value="PT_P3000"/>	List Price:	<input type="text" value="110.00"/>
Description:	<input type="text" value="PC Connectable Label Maker for"/>	Discontinued:	<input type="text" value=""/>
Category:	<input type="text" value="Label Printers"/>	Taxable:	<input checked="" type="checkbox"/>
		Active:	<input checked="" type="checkbox"/>
Notes:	<input type="text" value="https://www.amazon.com/Brother-PT-P700-Connectable-Label-Maker/dp/B00DSYEB28/"/>		

At the bottom, there are "Cancel" and "Submit" buttons.

Part Number Sources

The Sources tab provides a list of the top 3 suppliers that your company uses to purchase specific parts and their specific part numbers. It also includes a drop-down to choose their manufacturer.

Note: Designate a 'Supplier' by navigating to the Companies grid, opening a record, and clicking on the [Email/Phones/People/Types] tab. Once there, within the Company Types panel grid click on the [Link] icon and select the "Supplier" status.

Edit Part Number [?] [X]

General | **Sources** | Equipment | Calibrations | Work Orders | Custom | Meta

Suppliers

First Choice: Amazon [v] Manufacturer: Brother [v]
 Part Number: 4578948 Part Number: PT-P700

Second Choice: Amazon [v]
 Part Number: 34598625

Third Choice: ACME Label Pr [v]
 Part Number: 23459812

Cancel [Submit]

Equipment & Calibrations

This tab contains panel grids used to link Equipment and Calibrations to a Part Number record. Use the buttons at the bottom of the dialog to link, edit, and unlink records. Click and drag on the middle gray bar to resize the panel grids.

Edit Part Number [?] [X]

General | Sources | **Equipment | Calibrations** | Work Orders | Custom | Meta

[Equipment] - 4 Records		
Equipment ID [v]	Model [v]	Descrip [v]
SAMPLE-004	CRG-101	Cylindric
SAMPLE-052	DSKM-66	Disk Mic
SAMPLE-081	LT	Loctite
SAMPLE-144	CLPR450	Caliper

[Calibrations] - 4 Records	
Equipment ID [v]	Date [v]
SAMPLE-001	09/11/2019
SAMPLE-061	05/16/2019
SAMPLE-019	05/08/2019
SAMPLE-211	04/02/2019

[+] [Link] [Edit] [Unlink] [X]

Cancel [Submit]

Work Orders

Create and manage work orders related to a part number. Click the green [+] icon to add a new Work Order record. Within this new dialog, add information related to a part number's 'Requestor', 'Department', 'Location', which account it's linked to, which category it's in, the priority level, and a space for a description. Link equipment to the work order if necessary.









Equipment	Model	Description
SAMPLE-0	HMIC32	Hole Microm
SAMPLE-0	DBG-72	Dial Bore Ga
SAMPLE-0	CAL-180	Caliper
SAMPLE-1	SM-4556	Scissor Micr

Other Tabs

- **Custom:** Visit the Change Field Names help topic to learn how to customize this tab for your company's needs.
- **Meta:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it.

People

Find the People grid in the Data Grids tab of the ribbon menu. Use the Edit Person dialog in Calibration Control (our **Calibration Management Software**) to edit all information related to the People you interact with for your company and link them to other records using the panel grids.

File	Calibration	Tools	Data Grids	Utilities	Help		
 New Record	 Equipment	 Calibrations	 Documents	 Companies	 People	 Charts	 Calendar
	Equipment		Data Grids			Visual	

All People						
[All People] - 118 Records						
Person	Title	Technician	Active	Department		
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Alberta Mayfield	Manager	<input type="radio"/>	<input checked="" type="checkbox"/>	Production		
Andrew Bradford	Manager	<input type="radio"/>	<input checked="" type="checkbox"/>	Inspection		
Andy Reyes	Supervisor	<input type="radio"/>	<input checked="" type="checkbox"/>	Quality Assuran		
Antonia Jackson	Director	<input type="radio"/>	<input checked="" type="checkbox"/>	Engineering		
Bethany Bitters	Director	<input type="radio"/>	<input checked="" type="checkbox"/>	Shipping		
Betty Niehuhr	Manager	<input type="radio"/>	<input checked="" type="checkbox"/>	Shinninn		

Person Dialog

Add information to a Person record:


- **Salutation/First Name/Last Name:** First and Last Name of the Person, Salutation when applicable. These are the only required fields for a new Person record.
- **Display Name:** Auto-filled by First and Last Name fields.
- **Title:** The Person's job title, if available.
- **Department:** The Department the Person is currently assigned to in their Company.
- **Location:** Use this drop-down when applicable.
- **Status:** The Person's status within the company they work for.
- **ID Number:** The Person's unique company ID, if applicable.
- **Is Active:** Select whether or not the Person has an active account in your records.
- **Primary and Alternate Addresses:** The Person's address(es), if applicable.
- **Notes:** A blank space for notes related to the Person record.

Update Equipment Locations

Click on the [Update Locations] button at the bottom of this dialog to synchronize the Department and Location fields of all Equipment records to include the selections chosen for this Person.

Grids Dialog

- **Email Addresses:** Add or link an Email address to a person record.
- **Phone Numbers:** Add or link phone and Fax numbers to a person record.
- **Companies:** Add or link a person record to an existing company record.
- **Attachments:** Add or link attachments to a person record.




Other Tabs

- **Notes Tab:** Add and link Notes to a Person record.
- **Custom Fields Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.

Documents Dialog

Procedures Are Now Known as Documents


Easily link company Documents (e.g., Calibration Procedures) and related Equipment in Calibration Control (our **Calibration Management Software**). Find the Documents icon under the Data Grids tab of the ribbon menu.



Edit Document Dialog


Double-click on a Document record to open and edit it. The Edit Document dialog contains 4 tabs:

- **Document:** Add information related to the Document including its number, title, active status, type, class, subject, Document URL, and a space for notes. Click on any of the grey shaded boxes to add a new field for that drop-down.
- **Equipment/Calibrations Panel Grids:** Use the panel grids for adding and linking equipment and calibrations to a Document record. These panel grids can be filtered and sorted using the Filter Row at the top of each grid.
- **Custom Tab:** Visit the Change Field Names help topic to learn how to customize this tab for your company's needs.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it.




Link Documents to Equipment Records


- **From the Documents Grid:** Within the Grids tab of the Document dialog click on the 'chain link' icon to link a specific piece of equipment to that Document.



- **From the Equipment Grid:** Within the Grids tab of the Edit Equipment ID dialog click on the 'chain link' icon within the Documents panel grid.




Link an existing Document or Equipment ID by highlighting its row and clicking the [Select] button.



Email Addresses

Organize email addresses related to Companies and People within Calibration Control (our Calibration Management Software). Find the Email Addresses grid in the Data Grids tab of the ribbon menu.




Email Address Dialog

Add a new email address by clicking in the gray area of the grid, or right-click for the context menu and select [New

Record]. Enter the email address and select the 'Type' (i.e. personal, work, or other). Check boxes are available for selections of 'Primary', 'Opted Out', and 'Invalid' which display as green checkmarks in the grid.


Companies Tab

Use the Companies panel grid to link email address records to a specific Company. Use the icons at the bottom of the panel grid to add, link, edit, or break a link (i.e. remove a linked record without deleting it from the database).



People Tab

Use the People panel grid to link email address records to a specific Person. Alternatively, double-clicking in the gray area of the panel grid will open a New Person dialog, and clicking on an existing record will open the Edit Person dialog.




Other Tabs

- **Custom Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it; including the dates and times these occurred.

Phone Numbers

Organize phone numbers related to Companies and People within Calibration Control (our **Calibration Management Software**).




Phone Number Dialog

Add a new phone number by double-clicking in the gray area of the grid, or right click for the context menu and select [New Record]. Enter the phone number with area code in the desired format and then select the number type from the drop-down menu (i.e., Office, Main, Fax). Check boxes are provided to designate the phone number as 'Primary', 'Opted Out', or 'Invalid'.


Companies Tab

Use the Companies panel grid to add and link Companies to a specific phone number record. Double-click in the gray areawithin the panel grid to add a new Account record, or double-click on an existing record to edit.



People Tab

Use the People panel grid to add and link People to a specific phone number record. Double-click in the gray area within the panel grid to add a new Person record, or double-click on an existing record to edit.



Other Tabs

- **Custom Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it; including the dates and times these occurred.


Notes

Create a note within an existing record to organize extra information. Changes such as an asset transfer, asset issue,

asset receive, check out, and status change are also recorded as note records.

Notes Grid


Access the Notes grid by clicking on the notepad icon in the Data Grids tab of the ribbon menu.



	Date	Subject	Note	Type
<input checked="" type="checkbox"/>	09/10/2019	Asset Transfer - S	Changed Site to J	Custody Transfer
	09/10/2019	Equipment history	notes go here.	Observation
	08/07/2019	Check Out	Checked Out	Check Out
	08/07/2019	Check In	Checked In	Check In
	08/07/2019	Check Out	Checked Out	Check Out
	08/07/2019	Check Out	Checked Out	Check Out

Notes Dialog

Use the Notes dialog to organize information related to the Note and its relationship to Equipment, Companies, and People. This dialog also includes a Custom and Meta tab.



Edit Note

Note | Equipment | Companies | People | Custom | Meta


Date: 3/22/2019 Type: Check Out

Subject: Checked Out to Derek Siler


Note: Checked Out to Derek Siler and Status Changed to 'Not for Product Use' for Work Order 202

Cancel Submit


Easily keep track of all the equipment records associated with a specific Note by linking them using the panel grid in the Equipment tab. From within the 'Edit Equipment' dialog (in the Equipment grid), add a new Note to a record and it will automatically be available in the Notes grid.



In the Companies tab, add or link Companies to a Note record. When a Note is added to a Company record it will be available in the Notes grid.



Add and link People to a Note record using the People tab of the Notes dialog. When a Note is added to a Person record it will show up in the Notes grid.




Other Tabs

- **Custom Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it; including the dates and times these occurred.

Application Codes

All of the simple codes (used in drop downs and other locations) are kept in a single Codes grid for easy reference and update. To edit a Code record, click on the Codes icon within the Data Grids tab of the ribbon menu.

By default the Codes grid is in groups for ease of use; click on a [+] icon from the left side to expand each category and double-click on a record to edit it.



Codes List

To view a full list of Code groups, open the Codes grid in the Data Grids tab. You can also visit our Codes help page.

Adding a New Code

Right click in the Codes grid to bring up a new Add Code dialog. Required fields are highlighted **PINK** by default.

- **Code:** The code name that appears in the drop down menu.
- **Code Short:** Acronym or a shortened version of the code name.
- **Group:** The category that the code belongs to, such as Locations.
- **Sequence:** The number order this code will appear in the drop down menu.
- **Is Active** checkbox: Allows you to show or hide a code in the dropdown menu.
- If you select the Locations group, you will have a **Sub Location** dropdown menu to allow you to add a sub location. To do this, create a new code record and choose the Location Types group. Save the record and it will become an option for Sub Locations.
- If you select the Status - Codes Equipment group, selecting the **Is Hidden** checkbox allows you to hide an Equipment record if the status code is marked Is Hidden.

Add Code (Load Time: 0.0 Seconds)

Code Custom Meta

Code:

Code Short:

Group:

Sequence: Is Hidden Is Active

Notes:

Cancel Submit

Custom Tab

Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the [Change Field Names](#) help topic to learn how to edit these fields.

Add Location

Code Custom Meta

Character 0:

Character 1:

Character 2:

Character 3:

Character 4:

Date 0: Number 0: Check 0:

Date 1: Number 1: Check 1:

Date 2: Number 2: Check 2:


Date 3: Number 3: Check 3:

Date 4: Number 4: Check 4:

Cancel Submit

Meta Tab

This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it, including the dates and times these occurred.



The screenshot shows a dialog box titled "Edit Location" with three tabs: "Code", "Custom", and "Meta". The "Meta" tab is selected and displays the following information:


Created:	11/28/2017 6:08:47 PM
Created By:	
Updated:	12/7/2018 12:31:34 PM
Updated By:	APELAPTOP

Below the table, there is a red note: "Note: This page is read-only." At the bottom of the dialog, there are "Cancel" and "Submit" buttons.

Program Options

Modify Application Settings

You can modify Calibration Control (our **Calibration Management Software**) program options by going to the File menu and selecting Options or access the Options dialog in the Utilities tab of the ribbon menu.




General Tab

- **First Day of Week:** By default, Calibration Control uses the System defined first day of the week. Otherwise, you can define the first day of the week using the dropdown. This setting determines the due date for calibrations using

the 'Week of Frequency.

- **Due Date 'Week/Month of' Offset:** When calculating the calibration due date for the 'Week of' and 'Month of' frequencies, the Calibration Due Date is calculated by adding the Week of/Month of units to the Last Calibration and then choosing the first day of the Current OR Next week/month.
- **Check for Program Updates:** Calibration Control can notify you of available software updates on a Daily, Weekly, or Monthly cycle when the application is first opened. By default, you will receive daily updates.
- **Measurement System:** Use the drop-down to switch the default measurement system between U.S. and Metric.
- **Equipment ID of Default Template:** Enter the Equipment ID of an existing Equipment record to be used as a template when creating new Equipment records.
- **Unique Alternate Equipment ID:** Determine whether Alternate Equipment ID's should be unique or not, or be unique within a given Site. Note that Equipment ID's are required to be unique.
- **Language:** Select the language displayed. Changing the language requires a program restart. Alternatively, you may select a language for each specific User.
- **Label Printer:** When multiple label printers are installed on a single computer, identify the default label printer in this field.
- **Start Screen:** Select your default start screen when opening Calibration Control. By default and when no selection is made, the Equipment grid will be the start screen. To have no start screen, select "None" from the dropdown.
- **Do Not Copy Default Labels, Reports & Emails on Start:** If any default labels or reports are missing, CalibrationControl will add the missing files to your files folders. Check this box to prevent the action.
- **Use Active Equipment as Default for Asset Transfer:** If this box is checked, any selected Equipment record will automatically be populated in the Asset Transfer dialog.
- **Dialog Font Size:** Define the size of all dialogs with 0 being the default. You can also use negative values to shrink the dialog size.



Advanced Tab

The Advanced tab is split up into separate dialogs for each category. Click each link below to learn more about each option.

Calibration

Check Out

Email

Folder & Path

Highlight Colors

Reference Records - Add Missing

Reference Records - Refresh All

Report & Label

Security


Switch to Sample Database*

Label Edit
Masks - Date
Masks - Number

Spell Checker
Visibility

* When using the sample database, this button will show "Refresh Dates in Sample DB" which will update sample calibration due dates to more recent dates.

- **Allow Multiple Equipment Browse Grids:** When this option is selected, you have the ability to open additional equipment tabs at once.
- **Remember Grid Layouts:** When this option is selected, your layouts for all grids are saved and stored as a layout file. When you sign back in to Calibration Control, your grid layouts will remain as you left them, otherwise, they are reset to default.
- **Show and Hide Fields in Grids:** When this option is selected, you can right click any grid and select the Show Fields dialog which allows you to show or hide fields in the specified grid, otherwise, the Show Fields item will be greyed out and inactive.



Options - Calibration

Within the Option dialog's Advanced tab, modify Calibration Control's Calibration settings to fit your company's standards.

- **Status on Calibration Pass:** Select a status that will be applied to an Equipment record when a new Calibration event has Passed.
- **Status on Calibration Fail:** Select a status that will be applied to an Equipment record when a new Calibration event has Failed.
- **Received for Calibration Status:** When using the Status Change dialog, if the selected status is set, three changes are made to the related Equipment dialog: **(a)** the Status is changed to the selected status, **(b)** the Received field is set to TODAY, and **(c)** a note record is created describing the event.
- **Reminder Lead Days:** The number of days that determine which equipment is "soon due" and will appear in the DueCal report and when browsing equipment that is due for calibration.
- **Max Measure Groups to Submit at Once:** The maximum number of measurement groups you can submit at once. By default, the value is 0 which indicates there is no max.
- **Second Pause between MG Submit:** The number of seconds to pause between each measurement group being submitted. By default, the value is 0 which indicates there is no pause between submitting measurement groups.
- **Manually Set Due Dates in Equipment Dialog:** Normally, the Calibration Due (Next Cal) is automatically calculated based on the Last Calibration and the Calibration Frequency. Therefore, checking this option for manual set will disable the automatic calculation.
- **Ask to Print OOT Worksheet When Received OOT:** With this option checked, the User who saves a Calibration Event with the Received Out-Of-Tolerance (OOT) condition will be asked if the OOT Worksheet should be printed. (Works with the next option below.)
- **Null Cal Received Date on Successful Cal:** (Usually used in conjunction with the above Received for Calibration Status feature) With this option checked, when a new Calibration is saved with a 'Passed' Status then the Received [for Calibration] field in the Equipment record is automatically nulled.
- **Use Last Calibration as Template if No Template Linked:** With this option checked, creating a new Calibration Event uses the most recent Calibration Event as the Measurement Template if no Template is otherwise linked for that Equipment.
- **Prevent Entry of Future Date for Calibration Dates:** Prevent the manual entry of a calibration due date in the future.
- **Work Days:** Select which days of the week are considered work days for calibrations.

Calibration Options

Status on Calibration Pass: Active

Status on Calibration Fail: Inactive

Reminder Lead Days: 30

Received for Calibration Status: In Calibration

Manually Set Due Dates In Equipment Dialog

Ask to Print OOT Worksheet When Received OOT

Null Cal Received Date on Successful Cal

Use Last Calibration as Template if No Template Linked

Prevent Entry of Future Date for Calibration Events

Block Printing of Due Cal Labels for Failed Calibrations

Work Days

S M T W T F S

Cancel Save

Options - Check Out

Configure settings for the Check Out dialog in Calibration Control and CC Terminal.

- **Prevent Check Out if already Checked Out:** Users will be required to Check In an equipment before attempting to Check Out.
- **Default Use Count Check In/Out Quantities:** Define the Default Use Count quantity for Check In and Check Out using the corresponding fields.
- **Auto Assign Authenticated User in Check Out Dialog:** When using User Authentication (Sign-in mode), the user currently signed in will automatically populate the Person field.
- **Do Not Allow Check Out of Hidden Equipment:** Prevent Equipment with a hidden Status from being checked out.
- **Only Current User Can Check Out Equipment:** When using User Authentication (Sign-in mode), the user currently signed in will be the *only* user allowed to Check Out equipment.
- **Require Fields:**
 - Person Field During Check Out
 - Work Order Field During Check Out

◦ Status Field During Check In AND Out

Check Out Options

Prevent Check Out if already Checked Out

Default Use Count Check Out Qty:

Default Use Count Check In Qty:

Auto Assign Authenticated User in Check Out Dialog

Do Not Allow Check Out of Hidden Equipment

Only Current User Can Check Out Equipment

Require Person Field During Check Out

Require Work Order Field During Check Out

Require Status Field During Check In AND Out

Cancel Save

Options - Email

Use the fields in the Email tab to setup a connection to your company or personal email server. Calibration Control makes use of these settings in the Auto Notify feature that sends scheduled Email notifications or messages to People identified in the Calibration Control.

- **SMTP Server:** An SMTP server address usually looks something like SMTP.YOURDOMAIN.COM .
- **Email/Password:** This is the email address (or user name) and password of the email account Ape will be sending mailfrom.
- **Seconds Between Emails:** Enter the number of seconds to wait between sending each email. By default, this value is zero.
- **Port:** By default, the port number is 0. Otherwise, enter the port number required by your email server.
- **Disable CSS in Emails:** The default Email templates use embedded cascading style sheets (CSS) to format text with colors, backgrounds, decorations, etc., especially in the table of Equipment. If the email services you send emails to blocks some CSS, like Gmail, check this box to get a list of some minimal coloring back into the Equipment grid.
- **Connection Settings:**
 - Connect Asynchronously
 - Secure Socket Options

- **Use Legacy Email (pre 9.2):** Check this box if you require the "Use Secure Connection" checkbox previously used in versions prior to 9.2. Alternatively, you can try using a Secure Socket Option of **SSL on Connect**.
- **Output Folder:** Optionally, select a folder path for a copy of all emails sent out. The Enabled checkbox must be checked for emails to be copied to the selected path.
- **Email File Mask** The default file mask can be slightly edited to your liking. {RND32} will contain the unique identifier of each email.
- **Send Test Email to:** Solely for the purpose of testing your email settings, enter an email address to send a test message to. Alternatively, you can choose to Send to Email Output Folder.


Options - Folder & Path

The **Files Folder** contains the three important sub folders of Attachments, Labels, and Reports. While the Attachments folder is the location where Calibration Control stores the actual folders, the Labels and Reports folders store templates (how the labels and reports are printed).

- **Location of Files:** It is important to move the Files Folder to a location where it will be backed up regularly and can be accessed by other Calibration Control users. Click the Edit button to browse to your desired location. Note: If you cannot find your desired location through this dialog, you can manually enter the path in your database. Take a look at the Folders help topic for more detailed information.
- **Location of Auto Generated Website:** This is the location where a simple website listing all equipment grouped by

Department. It can be created manually by clicking on the 'Publish Website' icon in the Utilities tab of the ribbon menu or scheduled via the Auto Notify utility.

- **THUM Database:** The THUM (Temperature and Humidity USB Monitor) database is the location where the THUM device stores its data. If the device is used and the database exists, then Ape needs to know where this database is



so temperature and humidity data can be added automatically to Calibration Event records.


Change Field Names

Change the Text of Field Labels


Rename any field of Calibration Control dialogs or forms for better usability.

Edit Form Labels

From Options (located in the Utilities or File tabs) navigate to the 'Advanced' tab and select the [Label Edit] button.




From the Edit Forms Labels dialog (shown above), use the 'Select Dialog' drop-down to select which dialog screen to modify. In the example below, selecting Company will display all fields within the Company dialog.



Edit Fields and Labels

After finding the label to modify, manually edit the **Control Text** (what the field says), **Updated By** (who modified it), and **Record Updated** (when it was modified).




Change the 'Help Tip' field to edit the text shown while the cursor hovers over the control. (Note: Not all controls have the ability to display 'Help Tip' text, like column headers in data grids.)

Reset Changes

If a mistake is made, set all the values back to default by clicking the [Reset THIS Dialog to Default] for the single dialog displayed, or by clicking [Reset ALL Dialogs to Default] for all dialogs. Additionally, double-click on any single record and right-click to select the [Reset THIS Record] option.

Modifications of Form Labels are highlighted green to show the custom changes.



Custom Fields

The Custom Tabs in dialogs are very useful. Rename these available fields with the same steps. For example, here is the default Custom Fields tab in the 'PersonEdit' dialog:

The screenshot shows the 'Edit Person Record' dialog box with the 'Custom Fields' tab selected. The form is organized into several sections:

- Character Fields:** Character 0 through Character 5 are represented by text input boxes. Character 6 and Character 7 are represented by larger, multi-line text input boxes.
- Date Fields:** Date 0 through Date 6 are represented by dropdown menus.
- Number Fields:** Number 0 through Number 6 are represented by text input boxes.
- Check Fields:** Check 0 through Check 6 are represented by checkboxes.

At the bottom of the dialog, there are 'Cancel' and 'Save' buttons.

And after renaming the Custom fields, it could display as shown below:


The screenshot shows the 'Edit Person Record' dialog box with the 'Custom Fields' tab selected after renaming. The form is organized into several sections:

- Text Fields:** Available, Conflicts, Requester, Emergency, Relation, and Phone are represented by text input boxes.
- Notes and Other:** Notes and Other are represented by larger, multi-line text input boxes.
- FTA Fields:** FTA 0 through FTA 6 are represented by dropdown menus.
- Number Fields:** Number 0 through Number 6 are represented by text input boxes.
- Check Fields:** Contract, Freelance, Local, Transport, Dependable, On-Call, and DND are represented by checkboxes.

At the bottom of the dialog, there are 'Cancel' and 'Save' buttons.


Options - Reference Records

Calibration Control's reference tables include: Auto Emails, Charts, Code Groups, Codes (only required codes), Company Types, Coverage Factors, Grid SQL (non-custom), Permissions & Roles. You can add missing reference records or refresh them (replace with default) with a simple click through the Options dialog's Advanced tab.




Add Missing

Any missing reference records will be added. If one or more records is accidentally deleted from tables like Language (used for field names), Roles (defines Admin, Supervisor, etc.), or Permissions (which Roles are required for each Permission), then pressing this button will re-add all missing records. You'll see a message that states how many records were added. Otherwise, you'll see a message that all reference records are up to date.



Refresh All

Refreshing all reference records will delete and readd the contents of all reference tables (listed in the screenshot below). Note the warning that this may overwrite some customizations in those tables.



Reports & Labels


Calibration Control uses the default Reports and Labels, however you have the ability to set any report as your default certificate, worksheet, etc. You can copy and customize any report or label and use it as an override with a specific equipment record or set it as a default report/label to use with all records.

Default Reports

These settings determine which reports are used by default. While most reports can be previewed using the Reports PrintMenu, *the Cal Certificate includes specific data so that it can only be previewed through an Equipment record.*

- **Calibration Cert:** Template used when printing a Calibration Certificate.
- **Due Cal Report:** Report template used when printing the Calibration Due Report from the ribbon menu.
- **Cal Worksheet:** Calibration Worksheet that prints when the Equipment record has a Measurement Template or a previous Calibration Event.
- **Blank Worksheet:** Calibration Worksheet that prints when the Equipment record has neither a Measurement Template nor a previous Calibration Event.
- **Uncertainty:** Default report for the Measurement Uncertainty Budget.
- **OOT Investigation:** Designates the Out-Of-Tolerance (OOT) Investigation worksheet used when Equipment received for calibration is found OOT.
- **Custom Report:** Custom report for a single Equipment record, which can print from the context menu in the Equipment grid.


Note: All reports can be edited using Ape software's built-in Report Designer. Custom changes made to default templates should be renamed, (e.g. "CalibrationCertificate - COMPANY").



Default Labels Tab

- **Default Equipment Labels:** These settings determine which label templates are used by default when using the 'Print Label' button from the Equipment dialog. (Default label settings can be overridden for a specific Equipment, Type, or Model record.)
 - **Calibration Due:** Prints when the calibration frequency generates a due date and uses labels with the 'Due Cal' string in their file name.
 - **Asset:** File templates are identified as Asset labels by the 'Asset' string in the file name.
 - **Cal Not Required:** Prints in place of the Calibration Due label when the Calibration Frequency is 'Cal Not Required' and is identified as CNR label by the 'CNR' string in the file name.
 - **Next Use:** Print in place of the Calibration Due label when the Calibration Frequency is 'Next Use' and the labels are identified by a 'Next Use' string in their file name.
 - **System:** Prints with information relating to an Equipment System.

Note: The number at the beginning of the file name indicates the width of the label tape in millimeters. The number at the end of the file name indicates the sequence. Any custom labels should be given a *unique* sequence. All labels can be edited using the Brother P-touch Editor software that comes free with your Brother P-touch label printer. For a preview of all default label templates included in Ape Software, visit our Label Fields topic.




Options - Security

Modify Calibration Control's advanced Security settings in the Options dialog. Easily activate User Authentication (sign-in mode) with just a click of a button, disable CC Terminal features, and more.

- **Activate/Deactivate Sign-In Mode:** Enable or disable the User Authentication Mode (Sign-In Mode). At least one User must have Administrator privileges to Enable User Authentication and, by default, only an Administrator can disable Authentication.
- **Password Security:** Edit User password strength, rules, and expiration times.
- **Inactivity Exit:** Adjust how many minutes pass before Calibration Control prompts a user to sign-in again.
- **Sign In with Windows User Name:** For ease of use, use a Windows account to sign in to Calibration Control.
- **Only Current Users Can Check Out Equipment:** Prevents anyone but the currently signed in user from checking out Equipment.
- **Prevent Users From Saving Passwords on Client PCs:** When the Authentication mode (sign-in) is enabled and a User signs in, they are given the option of saving their User Name and Password for future automatic sign-ins. Checking THIS option prevents the ability to save their User Name and Password.
- **Require Re-Entry of Password When Signing In:** Requires users to enter their password twice in order to successfully sign in.
- **No Authentication for Auto-Notifications:** Allows email notifications to be sent without prompting the user sign-in authentication dialog.

Disable Terminal Features for All Users

- **Disable Terminal Equipment Grid:** Prevents any user from seeing the main list of Equipment from the Terminal dialog.
- **Disable Check In/Out Dialog:** Prevents the use of the Check In/Out dialog when Terminal mode is activated.
- **Disable Status Change Dialog:** Prevents use of the Status Change dialog when Terminal mode is activated.
- **Disable Asset Transfer Dialog:** Either disable the Asset Transfer dialog completely or disable one or more of the fields that can be edited.




Feature Visibility Options

Simplify Calibration Control's ribbon menu by hiding features your company is not using.


Make Features Invisible

Feature Visibility is found under the Advanced tab in the program Options. Find the Options dialog under the Files tab or under the Utilities tab of the ribbon menu.




Click on Feature Visibility to view the available features you can hide. Select a feature (or group of features) to become invisible on the ribbon menu.

Note: Restart Calibration Control to apply these changes.




In the above image, notice how the feature group 'Advanced Tools' is marked to be invisible.

The ribbon with the 'Advanced Tools' feature group *visible* looks like:




The ribbon with the 'Advanced Tools' feature group *hidden* looks like:



When it is made hidden or invisible, a feature or a feature group is simply removed from view. To make a hidden feature visible again, return to the Feature Visibility dialog in Program Options and remove the feature's checkmark, then restart the program.

Minimize Ribbon


For a better view of the open windows, right-click on the ribbon and select the 'Minimize the Ribbon' option. Bring it back the same way, by right-clicking where the ribbon would normally appear.



Options - Spell Checker

You can enable spell check mode in Calibration Control through the Spell Checker dialog (Advanced Options). There are several spell check modes to choose from plus further configuration options in the checkboxes shown below.

- **Mode:** Choose from the following spell checker options:
 - **None:** By default, Calibration Control does not use the Spell Check mode.
 - **As You Type:** Underlines incorrectly spelled words in red as you type. This option does not offer suggested correct spelling.
 - **Dialog On Validating:** Upon submission of any dialog, a spell check window will appear if any spelling errors are found, which includes correct spelling suggestions.
 - **Dialog On Validating And As You Type:** This option shows a red line under incorrectly spelled words and a spell check window upon submitting any dialog with spelling errors.
- Enable the User dictionary when spell check mode is active.
- Enable the option to ignore words in UPPERCASE, Mixed CaSe and words that contain numbers.
- Enable spell check on Compound Words and Hyphenated words
- **Maximum Word Length to Check:** Enter the maximum word length to check. By default the value is 80.
- **Language:** Select the language for the spell checker. Available languages: English, Dutch, French, German, Italian, Portuguese or Spanish.




Measurements in Metric or Imperial

Display Metric or Imperial as Default MeasurementUnits


Use the Options dialog to change the default measurement system to either Metric or Imperial (US) units in CalibrationControl (our Calibration Management Software).

Go to File > Options and in the General tab of the Program Options dialog, select either 'Metric' or 'US' as the defaultmeasurement system from the drop down menu.



Asset Transfer Dialog

Transfer Assets quickly to new Sites, Departments, Locations, Custodians, or change the Status in Calibration Control (our **Calibration Management Software**). Administrators can control which fields are shown by adjusting permissions in the Options dialog under the Security tab. Transactions made using the Asset Transfer dialog are automatically logged in that Equipment dialog's Notes tab.



Open the dialog by clicking on the Asset Transfer icon in the Tools tab of the ribbon menu. This dialog can also be opened directly from the Equipment grid by right-clicking on the Equipment record to be changed and selecting 'Asset Transfer'.

Use the Asset Transfer Dialog

First type an Equipment ID, Master ID, or Serial Number in the 'Equipment ID' box or select the [Use Scanner] button first when using a barcode scanner. Then select the desired fields for a new Site, Department, Location, Custodian, or

Status and click [Validate], then the [Submit] button. Fields left blank will not be changed.

Asset Transfer

Use Scanner

Equipment ID: SAMPLE-047

BM-64 (Blade Micrometer)

Transfer To

Site: Jensen and Watts

Department: Machine Shop

Location: Locker

Custodian: Gerry Anderson

Status: Active

Historical Note

Dated: 09/10/2019 04:33 PM

Cancel Validate Submit

Asset Issue & Receive Dialogs

Asset Issue Dialog

The Asset Issue dialog quickly issues equipment using specified values and saves a history note to that specific equipment's Notes tab. Open the Asset Issue dialog by clicking on the red arrow icon in the Tools tab of the ribbon menu. For convenience, the Site, Department, Location, Custodian, and Status selections will be remembered by the software until they are changed.



3-Step Process to Issue Equipment

- **STEP 1:** Select the Site, Department, Location, Custodian, or Status that the equipment is going to be issued to. Unused selections may be left blank.
- **STEP 2:** Check the box next to the desired Print and Update selections.

- **STEP 3:** Type the equipment ID in the provided box or check the 'Scanner' box to scan an asset tag, and then click[Submit].

The new information will be displayed in the info box once the Asset Issue has been submitted.

Asset Issue

STEP 1: Select Transfer To Values

Site: Jensen and Watts

Department: Engineering

Location: Cabinet

Custodian: Andrew Bradford

Status: Limited Use

STEP 2: Worksheet Print Preferences

To Printer Print Preview

Update First Update After

STEP 3: Scan Asset Tags

Equipment ID: Scanner

Equipment ID SAMPLE-303. Changed Site to Jensen and Watts from Blank, Department to Engineering from Manufacturing, Location to Cabinet from Test Area, Custodian to Andrew Bradford from James Pinkerton, Status to Limited Use from Active

Close Submit

Asset Receive Dialog

The Asset Receive dialog has all the same features as the Asset Issue dialog. Open the Asset Receive dialog by clicking on the green arrow icon in the Tools tab of the ribbon menu.



3-Step Process to Receive Equipment

- **STEP 1:** Select the Site, Department, Location, Custodian, or Status that the equipment is going to be received to. Unused selections may be left blank.
- **STEP 2:** Check the box next to the desired Print and Update selections.
- **STEP 3:** Type the equipment ID in the provided box or check the 'Scanner' box to scan an asset tag, and then click[Submit].

The new information will be displayed in the info box once the Asset Receive has been submitted.


Status Change Dialog

Quickly Change Status of Multiple Equipment

Change the status of multiple pieces of equipment quickly and accurately in Calibration Control. For instance, when receiving several pieces of equipment for calibration (or any other status), change the status of each tool with a single barcode scan after setting the "Change to What Status?" field.

Status Change in Ribbon

Open the Status dialog from the Status Change option.



Using the Status Change Dialog

Select a new status code for a piece of equipment in the "Change to What Status?" field combo-box. If entering the Equipment ID manually, uncheck the "Use Barcode Scanner for Data Entry" checkbox, and manually type in the Equipment ID and click [Submit]. Repeat with as many pieces of equipment that require the new status. When all necessary status changes are completed, click [Finished] to exit the Status Change dialog.

If using a barcode scanner, ensure the "Use Barcode Scanner for Data Entry" box is checked and scan the barcode on the asset label previously printed from Calibration Control. There is no need to click the [Submit] button when using a barcode scanner (due to carriage return setting in scanners). Repeat the equipment scan, or data entry, until all necessary status changes are complete. Click [Finished] to close out of the Status Change dialog.

Received for Cal with Status Change Dialog

By default, using the Status Change activity dialog will **(a)** update the Status code in the Equipment record and **(b)** save a Note in the Equipment record and Change Log describing that this status change took place.

Additionally, if a specific Status code is defined in the Options dialog (Advanced tab/Calibrations), then the Received Date auto-populates in the Equipment record when this Status code is used in the Status Change dialog. Check the box to null the Received (for Cal) Date field in the Equipment record whenever a new Calibration Event is saved with a 'Passed' Status.

Alternatively, if the Equipment Status code is changed *from* its designated Received status using the Status Change dialog, then the Received (for Cal) Date automatically clears in the Equipment record.

Equipment Check In/Out

Use the Equipment Check Out feature to quickly issue out items while automatically creating an audit trail, tracking UseCount or Use Days, linking to a Person, and connecting to a Work Order.

Opening the Check Out Dialog

Select the Check Out feature icon from the Tools tab of the ribbon menu. This option is best when you are using a barcodescanner.



Alternatively, you can right-click on an Equipment record from the Equipment grid and select [Check Out/In] from the context menu. Checking Out equipment directly from the Equipment grid will automatically generate the Equipment ID number in the Check Out dialog.

Check Out

To check out Equipment, select the 'Out' checkbox, enter a valid Equipment ID, and click the [Submit] button.

Enable the checkboxes next to optional fields to track the Person, Work Order, or Status Change related to the current Equipment Check Out.

A screenshot of a 'Check Out' dialog box. At the top, it asks 'Check In or Out?' with 'Out' selected (checked) and 'In' unselected. There is also a 'Use Scanner' checkbox which is unselected. An 'Increment' field is set to '0'. Below this, the 'Equipment ID' field contains 'SAMPLE-172'. The 'Person' dropdown menu is set to 'Trevor Williams' and has a checked checkbox next to it. The 'Reference Number' and 'Change Status To:' fields are empty, each with an unselected checkbox to its right. A large 'Response Area' is located at the bottom of the dialog. At the very bottom, there are three buttons: 'Finished', '* Required if Checked', and 'Submit'. The 'Submit' button is highlighted with a red rectangular box.

Confirmation of Check Out will appear in the Response Area of the dialog:

Check Out

Check In or Out?
 Out In

Use Scanner

Increment
0

Equipment ID:

Person:

Reference Number:

Change Status To:

Equipment ID SAMPLE-172
Checked Out to Trevor Williams

Finished * Required if Checked Submit

Check In

To return the Equipment, select the 'In' checkbox, enter a valid Equipment ID, and click [Submit]. Confirmation of Check In will appear in the Response Area of the dialog:

Check Out

Check In or Out?
 Out In

Use Scanner

Increment: 1

Equipment ID: _____

Person: _____

Reference Number: _____

Change Status To: _____

Equipment ID SAMPLE-172
 Checked In
 (Location: Test Area)

Finished * Required if Checked Submit

Check Out by Equipment System Name

Equipment System names can be used in place of an Equipment ID, making it appear like an Equipment Master ID, and it does NOT have the same name of an Equipment ID. This will check in/out all equipment in that system.

Check Out with Barcode Scanner

Enable the checkbox for "Use Scanner" to enable barcode label scanning. Using a scanner for field input will disable manual entry.

Check Out

Check In or Out?
 Out In

Use Scanner

Increment: 1

Equipment ID: _____

Using Scanner for Barcoded Equipment ID

The default barcode used in Calibration Control is the Equipment Master ID (a unique sequential number for each Equipment).

However, the Equipment Master ID can be overridden in barcode labels for an Equipment ID or Equipment System Name

using the following format:

- Upper Case A-Z (Capital Letters)
- Integers 0-9 (Numbers)
- Special Characters: Plus (+), Minus/Dash (-), Period (.), Forward Slash (/), Space ().
- **Note:** Barcodes must begin and end in alpha-numeric form, surrounding special characters.

Barcode Recognition Order

Sometimes scanning barcodes in Check Out can reveal data conflicts, so it is helpful to know the barcode recognition order in Calibration Control. When a barcode is scanned, Calibration Control performs database tests to determine what information (string) the barcode contains. The string is searched for in the following order:

1. Is it an Equipment Master ID? (Default in Calibration Control)
2. Is it an Equipment ID?
3. Does the string appear in the Barcodes field of an Equipment record?
4. Is it an Equipment System Name?
5. Is it the ID of a Person? (if a Person entry is required)
6. If none of the above, the string will be placed in the Work Order field.


Equipment Use Count or Use Days

Some equipment or tools require calibration after a certain number of uses or days used, as opposed to a time interval (e.g., by units of months or years). We refer to tools with this kind of Calibration Frequency as 'Equipment Use Count' or 'Equipment Use Days'.

Tracking Uses: The Increment field in the Check Out dialog only affects Equipment Use Count and Equipment Use Days, and only if its Next Calibration date is not overridden.

Increment Value: By default, the Use Count Increment during Check Out equals zero (0), and Check In equals one (1).

Manual Override: Increasing the Increment that automatically appears in the Check Out dialog will override the default values to apply the new value to that Equipment Use Count.



Default Increments in Check Out

The Response Area of the Check Out dialog will recognize Equipment Use Count / Days.

The screenshot shows a 'Check Out' dialog box with the following elements:

- Check In or Out?:** Radio buttons for 'Out' (checked) and 'In'. A 'Use Scanner' checkbox is also present.
- Increment:** A spinner box set to '1'.
- Equipment ID:** An empty text input field.
- Person:** A dropdown menu with an empty selection and a checkbox.
- Reference Number:** An empty text input field with a checkbox.
- Change Status To:** A dropdown menu with an empty selection and a checkbox.
- Central Display:** A large box containing the text:
 - Equipment ID SAMPLE-123
 - Checked Out
 - 19 Uses Left Before Calibration Is Required. (highlighted with a red border)
- Buttons:** 'Finished', '* Required if Checked', and 'Submit' (highlighted with a blue border).

Equipment Use Count (+1 per Use)


The default Increment during Check Out is zero (0), increasing the Use Count by nothing. Equipment Use Count updates during Check In to increase by one (1), the default value. In other words, return the tool or check it out again (if applicable) to add one count of usage in the Equipment record.

Equipment Use Days (+1 per Day)


Check In increases the Use Count (Days) by one (1) Increment for each day the Equipment was checked out. The nature of this calibration frequency is time sensitive so tracking days of usage also occurs during Check Out and may update the Use Count (Days):

- If the tool is checked out multiple times (before Check In), then Calibration Control will apply its automated algorithm to count the number of days this tool has been out and used.
- If the previous Check Out occurred on another date, then the Use Count (Days) adds one (1) Increment for each day, just like Check In would.
- If the previous Check Out occurred on the same date, then there is nothing to change.
- Note that if Check In occurs on the same day as Check Out, by default, it counts as one day of use.

Options: Check Out options are available in the Advanced tab of the Program Options dialog.



Verify Check Out Status








Current Status: In the Equipment dialog, select the 'Misc' tab to verify the current Check Out status with the Person (if applicable) and date. The Use Count is also visible, which resets to 0 after calibration.

Override Count: You are able to override any of these fields to include current and accurate information.

Check Out History: All Equipment Check Out activity is automatically logged in the Notes panel grid of the Equipment dialog located under the Grids tab. Double-click any note to view or edit.

[All Notes] - 3 Records

Date	Subject	Type
06/17/2019 02:47 PM	Check Out	Check Out
06/17/2019 02:47 PM	Check In	Check In
06/17/2019 02:46 PM	Check Out	Check Out




























Check Out Log

In Version 9 and later, all check in/out history is available in the Check Out log in the Data Grids tab. By default, the top 100 records are shown. Use the dropdown to change the number of records that appear. In addition, you can use the searchbar plus any of the filter rows to show only the data you want to see.

Ape Calibration Control (version 9.2.3)

File Calibration Tools **Data Grids** Utilities Help

Where and Who Work Flow Miscellaneous Logs Metadata

Equipment **Check Out Log** Equipment Systems

Select Top 100

[Check Out Log] - 19 Records

Equipment ID	Check Out	Increment	Check Out To	Reference	Status Changed	Scanner Used	Created
SAMPLE-205	<input checked="" type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:51 AM
SAMPLE-158	<input checked="" type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:51 AM
SAMPLE-194	<input checked="" type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:51 AM
SAMPLE-002	<input checked="" type="checkbox"/>	1				<input type="checkbox"/>	02/11/2021 10:51 AM
SYS: 777777	<input checked="" type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:51 AM
SAMPLE-194	<input type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:46 AM
SAMPLE-002	<input checked="" type="checkbox"/>	1				<input type="checkbox"/>	02/11/2021 10:46 AM
SAMPLE-158	<input type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:46 AM
SAMPLE-205	<input type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:46 AM
SYS: test2	<input checked="" type="checkbox"/>	0				<input type="checkbox"/>	02/11/2021 10:46 AM
SAMPLE-072	<input type="checkbox"/>	0	Ron Harrington		In Calibration	<input type="checkbox"/>	02/11/2021 09:57 AM
SAMPLE-024	<input checked="" type="checkbox"/>	1			Active	<input type="checkbox"/>	02/11/2021 09:56 AM
SAMPLE-024	<input type="checkbox"/>	0	John Reilly		In Calibration	<input type="checkbox"/>	02/11/2021 09:56 AM
SAMPLE-018	<input type="checkbox"/>	1				<input type="checkbox"/>	02/11/2021 09:55 AM
SAMPLE-018	<input type="checkbox"/>	1				<input type="checkbox"/>	02/11/2021 09:54 AM
SAMPLE-009	<input type="checkbox"/>	1				<input type="checkbox"/>	02/11/2021 09:54 AM
SAMPLE-005	<input type="checkbox"/>	1			Active	<input type="checkbox"/>	02/11/2021 09:54 AM

User Authentication [admin (Alberta Mayfield) Signed In] Using ACE Database Drivers Equipment Systems Load: 0.1 Seconds Sample DB

Recording and Researching Jobs

Record Equipment Used on Jobs for Traceability

The following dialog ('Record Job' in Tools tab) is available for quick recording of equipment used on Jobs or Work Orders and can be used with barcode scanners for increased data accuracy. Any string can be entered in the Job field while the Equipment field will recognize the label barcodes printed by Ape Software or the Equipment ID.

EquipmentID	Description	Job
ID-101	Printer	Job-101

Browse Jobs

The "Job Browse" grid is now the Check Out Log grid, which contains past Job entries. You can filter jobs by EquipmentID, Job, Dates, and several other fields. The results of a search / filter can then be printed or exported to Excel or PDF. To filter for Jobs, simply type "Job" in the search bar or filter for Job in the Check Out To column.

Equipment ID	Check Out	In	Check Out T	Reference	Status	Chan	Scanner U	Machine	Created	Created By
ID-101	<input checked="" type="checkbox"/>		0	Job	Test Job		<input type="checkbox"/>	BIANCA	03/18/2021	Alberts Mayfield

Auto Notify Utility

Send Emails and Desktop Notifications Automatically

Configure one or several Notification Events for scheduled reminders through email, desktop, or website publishing in Calibration Control. Send email notifications to specific groups of people such as Custodians, Technicians, Owners, Department Managers, etc. Receive a desktop notification of Equipment due for calibration.

Note: Auto notify only needs to be set up on ONE computer in order to avoid duplicate emails.

Email Settings

Before sending automatic Emails, you will need to add your email settings in the Options dialog . The database will need additional related data discussed at the bottom of this page.


User Authentication

With user authentication turned on, when an Auto Email is scheduled, Calibration Control will require you to enter your credentials and will not send the scheduled emails until you sign in. To prevent this from happening, go to File > Options > Advanced tab > Security. Then, ensure the checkbox for No Authentication for Auto Notifications is checked.

Auto Notify Grid

The Auto Notify utility is located in the Tools tab in the ribbon menu. Currently, 17 Auto Notify events are defined by default.

- **DueCal_[Person]:** Auto email notifications for the Person identified after the "DueCal_" string.
- **DueCal_Desktop:** Auto desktop notifications on the current computer.
- **DueMaint_[Person]:** Auto email notifications for the Person identified after the "DueMaint_" string.
- **Website_General:** Auto publication of a simple website listing all equipment in Department groups. Set the publish location in the Folders tab of the Options dialog.



Event	Email Subject	Frequency	Start	Hour	Minute	Active	Host
DueCal_CheckedOutBy	Equipment Due	Weekly	06/24/2020	8	30	<input type="radio"/>	
DueCal_Custodians	Equipment Due	Weekly	06/24/2020	8	30	<input type="radio"/>	
DueCal_DepartmentCo	Equipment Due	Weekly	06/24/2020	8	30	<input type="radio"/>	
DueCal_DepartmentMa	Equipment Due	Monthly	06/24/2020	8	30	<input type="radio"/>	
DueCal_Desktop		Weekly	06/24/2020	8	0	<input type="radio"/>	
DueCal_Owners	Equipment Due	Monthly	06/24/2020	8	30	<input type="radio"/>	
DueCal_Site	Equipment Due	Weekly	06/24/2020	8	30	<input type="radio"/>	
DueCal_Technicians	Equipment Due	Weekly	06/24/2020	8	30	<input type="radio"/>	
DueMaint_CheckedOut	Equipment Maint	Weekly	06/24/2020	8	55	<input type="radio"/>	
DueMaint_Custodians	Equipment Maint	Monthly	06/24/2020	8	55	<input type="radio"/>	
DueMaint_DepartmentC	Equipment Maint	Weekly	06/24/2020	11	30	<input type="radio"/>	
DueMaint_DepartmentM	Equipment Maint	Weekly	06/24/2020	11	30	<input type="radio"/>	
DueMaint_Owners	Equipment Maint	Weekly	06/24/2020	11	30	<input type="radio"/>	
DueMaint_Site	Equipment Maint	Weekly	06/24/2020	11	30	<input type="radio"/>	

Auto Notification Dialog

- **Email Subject:** The text that appears in the Subject line of the Auto Email Notifications.
- **Signature File:** Defines which HTML file is used as the email signature (Optional).
- **BCC:** Whenever an Auto Email is sent, a copy of the email can be sent to an additional email or list of emails. If more than one email is added to the BCC field, separate those emails with a comma, semicolon, or a hard return.
- **Schedule:** There are three frequency options available: Daily, Weekly and Monthly. For a Monthly frequency, select which Day(s) of the month from the Calendar. For a Weekly frequency, select which Day(s) of the Week.

Additionally in this section, Select the Is Active checkbox to automatically create a Task in your Task Scheduler which you can now easily navigate to using the Open Task Scheduler button.

- **Date & Time:** Select the starting date for this event and the specific hour and minute at which you would like the emails to be sent.

The screenshot shows a software dialog box titled "Edit Auto Notification DueCal_Custodians (Load Time: 0.1 Seconds)". It has several tabs: "Notification", "SQL", "Body Template", "Notes", "Custom", and "Meta". The "Notification" tab is selected. The dialog contains the following fields and controls:

- Event Name:** DueCal_Custodians
- Email Subject:** Equipment Due Cal
- Signature File:** (empty field with a red 'X' icon)
- BCC:** (empty field)
- Notes:** (empty text area)
- Schedule:**
 - Frequency:** Weekly (dropdown menu)
 - Is Active:**
 - Open Task Scheduler:** button
- Day(s) of Month:** A calendar grid with days 1 through 31. A "Clear" button is below it.
- Day(s) of Week:** A grid of checkboxes for days of the week: Mon, Tue, Wed, Thu, Fri, Sat, Sun. A "Clear" button is below it.
- Date & Time:**
 - Starting On:** 06/24/2020 (dropdown menu)
 - Hour:** 8 (spin box)
 - Minute:** 0 (spin box)
- Buttons:** Cancel, Send All Message(s), Send All Test Message(s) to BCC(s), and Submit.

Send Message(s)/Send Test Message(s)

- **Send Message(s):** Pressing sends all the email notifications to the defined Person group (e.g., Custodian, Technician, Manager), displays the desktop notification, or creates the default website.
- **Send Test Message(s) to BCC:** Sends the email notifications to the email(s) listed in the BCC field only as a test run.


SQL tab

SQL Statements

Previously located in your database, all required SQL statements can be viewed/edited here.


Table Previews

View a preview of the current email data. Double-click each user from the Recipients list to view their Equipment Due that will be included in their emails.



Body Template tab

Previously located in your database, all HTML Email Templates can be viewed/edited here.



Note: Prevent specific users from editing SQL/HTML by defining these skills/permissions in the User dialog (Advanced Skills section).

Email Activity Log

Under the Data Grids tab, in the Logs section, you'll find the Email Activity log. Here you can review details of all emails that were sent out.

Troubleshooting

The Auto Email Notifications require related data in the following tables:

- **People Fields in Equipment Record:** Identify the Custodian, Technician, Owner, or Checked out By fields so Calibration Control knows which equipment is related to who.
- **Department Field in Equipment Record:** Identify the Department in the Equipment dialog and then, in the Department record, identify the Manager and Contact people.
- **Site in Equipment Record:** Identify the Site (e.g., Customer or Company Site) in the Equipment record and then, in the People records, connect them to Companies.
- **Emails in People Records:** Ensure the People identified in other locations (e.g., Equipment or Department records) have valid email addresses and ensure they do NOT have "Do Not Email" checked in their person record.

Equipment Emails

In the tools tab, go to the **Equipment Emails** grid. Right-click and select which Auto Notify event you are troubleshooting from the bottom of the menu.

Easily filter all equipment that won't be included in that email notification by utilizing the **Passed** checkbox column. Additionally, all pink highlighted text indicates why an email is not being sent out.

Equipment ID	Description	Custodian	Status	Passed	Email Address	Email Date
SA001-001	Tool Kit	SA001-001	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-002	Indicator	SA001-002	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-003	Open Record	SA001-003	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-004	New Record	SA001-004	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-005	Close Record	SA001-005	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-006	Select All Rows	SA001-006	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-007	Copy to Clipboard	SA001-007	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-008	Select All in Cell	SA001-008	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-009	Filter Rows	SA001-009	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-010	Show Fields	SA001-010	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-011	Group by Size	SA001-011	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-012	Row Sorting	SA001-012	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-013	Refresh Grid	SA001-013	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-014	Reset Grid to Default	SA001-014	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-015	Print Data Grid	SA001-015	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-016	Export	SA001-016	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-017	Equipment with Due Dates	SA001-017	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-018	Equipment to Email for Custodians	SA001-018	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-019	Equipment to Email for Checked Out	SA001-019	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-020	Equipment to Email for Owners	SA001-020	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-021	Equipment to Email for Calibration Technicians	SA001-021	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-022	Equipment to Email for Maintenance Technicians	SA001-022	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-023	Equipment to Email for Department Managers	SA001-023	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-024	Equipment to Email for Department Contacts	SA001-024	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-025	Equipment to Email for Site Contacts	SA001-025	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-026	Practice Setup Test	SA001-026	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-027	Drop Indicator	SA001-027	Active	<input checked="" type="checkbox"/>	SA001@Example.com	
SA001-028	Drop Indicator	SA001-028	Active	<input checked="" type="checkbox"/>	SA001@Example.com	


Also, take a look at the SMTP Test Procedure help topic if you are having difficulties with sending Auto-Notify emails.

Terminal Mode

Minor Tasks in Ape Terminal Mode

Calibration Control in Terminal Mode (CC Terminal) is useful for those who require limited access. The simplified Terminal allows your Production employees or Shop personnel to make minor changes via feature shortcuts. Terminal users can *only* use CC Terminal. They will receive an error message upon starting Calibration Control stating they can only use the Terminal Mode.

Both the main Ape program and Ape Terminal install together. The Terminal is the same application executable (apecal.exe), except it is started with the switch -t1 (i.e., apecal.exe -t1).



Limited Use

Currently featured in the Ape Terminal dialog:

- Check Out/In
- Status Change
- Asset Transfer

- Equipment Grid
- Systems Grid
- Calibration History
- Print Menu

Further Limitations

Disable any Terminal features by going to Options > Advanced tab > Security. In the Disable Terminal Features for AllUsers section, check the features you'd like to disable.

The screenshot shows the 'Security Options' dialog box. At the top, there are two buttons: 'Deactive Sign-In Mode' and 'Password Security'. To the right is the 'Inactivity Exit (when Signed In)' section, which includes two input fields: 'Inactivity Minutes' and 'Message Seconds', both containing the value '0'. Below these are four checkboxes: 'Sign In with Windows User Name', 'Prevent Users From Saving Passwords on Client PCs', 'Require Re-Entry of Password When Signing', and 'No Authentication for Auto Notifications' (which is checked). The 'Disable Terminal Features for All Users' section contains a grid of checkboxes: 'Check In/Out', 'Status Change', 'Equipment Grid', 'Systems Grid', and 'Calibration History' are unchecked, while 'Asset Transfer', 'Site', 'Department', 'Location', 'Custodian', and 'Status' are checked. At the bottom, there are 'Cancel' and 'Save' buttons, with the 'Save' button highlighted by a blue dashed border.

Projects

Keep track of projects within your organization using the Projects grid and dialog. Find the Projects icon in the Calibration tab of the ribbon menu.



Projects Dialog


The 'Add/Edit Project' dialog keeps all information related to a project organized.

Project Dialog Fields

- **Project ID and Title:** Choose a project ID and give the project a title.
- **Description:** Add a description about the project.
- **Type & Priority:** Select the project type and its priority from the drop-down fields.
- **Requestor/Sponsor/Manager:** Choose which person is requesting the project, which person is sponsoring the project, and which person is the project's manager.
- **Status & Site:** Select the status of the project's completion and its location.
- **Start Est. & End Est.:** Select the estimated start and end dates for the project.
- **Start Act. & End Act.:** Select the actual start and end dates for the project.
- **Hours Est. & Hours Act.:** Choose the estimated and actual hours spent on completing the project.
- **Cost Est. & Cost Act.:** Choose the estimated and actual cost spent on completing the project.


Equipment/People/Parts/Work Orders

Use these panel grids to add, link, and edit records for a project's Equipment, People, Part Numbers, and Work Orders.



Notes/Attachments/Documents

Use these panel grids to add, link, and edit records for a project's Notes, Attachments, and Documents.



Other Tabs

- **Custom Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it; including the dates and times these occurred.

Work Orders

Submit | Service | Panel Grids | Other Tabs

Create and manage work orders within Calibration Control (our **Calibration Management Software**) and link them to specific equipment records.

Find the Work Orders icon in the Calibration tab of the ribbon menu.



Work Orders Dialog

The 'Edit Work Order' dialog contains 5 tabs useful in organizing each record.

Submit Tab

Create a work order number and add information like who requested it, which department and location it's at, which account it's related to, what category it falls under, and a priority level. Add a description of the work order and then link related equipment using the panel grid on the right side of the dialog. Click on any of the shaded boxes to quickly add new fields for its drop-down box.

The 'Edit Work Order' dialog box is shown with the 'Submit' tab selected. The dialog has a title bar with a question mark and a close button. The main area is divided into several sections. On the left, there are several input fields with labels and dropdown menus: 'Work Order:' (text box with 'WO-0006'), 'Requestor:' (dropdown with 'Idell Norton'), 'Department:' (dropdown with 'Machine Shop'), 'Location:' (dropdown with 'Shelf'), 'Account:' (dropdown with 'Production'), 'Category:' (dropdown with 'Repair'), and 'Priority:' (dropdown with 'Moderate'). Below these is a 'Short Description:' text box with 'Work order description....' and a larger 'Full Description:' text box. On the right, there is a table titled '[Equipment] - 4 Records' with columns 'Equipment I /', 'Model', and 'Descriptio'. The table contains four rows of data. At the bottom of the dialog, there are 'Cancel' and 'Submit' buttons.

Equipment I /	Model	Descriptio
SAMPLE-06	HMIC32	Hole Micro
SAMPLE-07	DBG-72	Dial Bore G
SAMPLE-09	CAL-180	Caliper
SAMPLE-14	SM-4556	Scissor Mic

Service Tab

Add work order information related to its status, due date, received and completed date, as well as the technician in charge, the work order's estimated labor time, labor cost, and material cost. This tab also includes a space to write out the actions taken to complete the work order.

The screenshot shows a software window titled "Edit Work Order" with a standard Windows-style title bar (minimize, maximize, close buttons). The window has a tabbed interface with the following tabs: "Submit", "Service" (which is the active tab), "Documents | Notes | Parts | Attachments", "Custom Fields", and "Meta".

Under the "Service" tab, the form contains the following fields:

- Status:** A dropdown menu with "In Process" selected.
- Technician:** A dropdown menu with "David Wilcox" selected.
- Due:** A date field with "3/7/2019" entered.
- Labor Time:** A numeric field with "9.2500" entered.
- Received:** A date field with "2/22/2019" entered.
- Labor Cost:** A numeric field with "221.0000" entered.
- Completed:** An empty date field.
- Material Cost:** A numeric field with "124.0000" entered.


Below these fields is a checkbox labeled "Check to Approve & Lock" which is currently unchecked.

At the bottom of the form is a text area labeled "Action Taken:" with the placeholder text "Write a description of the actions taken to complete the work order." The text area is currently empty.

At the bottom of the window are two buttons: "Cancel" on the left and "Submit" on the right.

Documents/Notes/Parts/Attachments Tab

Within these 4 panel grids you can add and link documents, notes, part numbers, and attachments to a work order. Filter and sort each panel grid using the gray shaded Filter Row.



Other Tabs

- **Custom Fields Tab:** Custom fields are available for adding Characters, Dates, Numbers, and Check boxes. Visit the Change Field Names help topic to learn how to edit these fields.
- **Meta Tab:** This is a READ-ONLY information tab used as a reference to view which user created the record and which user was the last to edit it; including the dates and times these occurred.

Custom Barcodes

Use Unique Barcodes Already on Equipment

If Equipment already has a unique barcode (e.g., asset or serial number), that number can be entered into its equipment record and used to locate and open its record. Make sure that this barcode is only used once per equipment record because there is no way for Calibration Control (our **Calibration Management Software**) to differentiate between two equipment records with the same barcode number.

Custom Barcode Field

Open an Equipment record and select the Misc tab. Enter the unique custom barcode in the 'Barcodes:' field in the Miscellaneous category (left side). The most accurate way of entering a custom barcode is to place the computer's cursor into the custom Barcode field and scan the barcode into the field.

Miscellaneous

Owner: John Reilly

System:

Barcodes: 0001574586

Equip Cost: 68.99

Acquired: 09/26/2018

In Service: 06/17/2019

Procedure:

Personal Property

Barcode Scan Dialog

Close the Equipment dialog, returning to the Equipment grid. Press the [F3] button to bring up the following dialog. When this dialog (below) is visible and a barcode is scanned, Calibration Control will know to look in the custom barcode field in addition to the default field.

Barcode Scan

Scan Barcode

... printed by CC, that represents a serial number, or is in barcode field of equipment record.

Cancel

Date Masks in Calibration Control

A mask is a way to format a field to make the data in that field look a specific way. Calibration Control gives you the ability to format some of the more important date, string, and number fields. Find 'Masks - Date' in the Advanced tab of the Program Options dialog.

Date Masks

Date/Time Masks & Formats for GRIDS

Date Mask: [LOC]mm/BB/yyyy

Date Format: MM/BB/yyyy

Date/Time Mask: [LOC]mm/BB/yyyy hh:mm S

Date/Time Format: MM/BB/yyyy hh:mm S

Date/Time Masks & Formats for CONTROLS

Date Mask: [LOC]mm/BB/yyyy

Date Format: MM/BB/yyyy

Date/Time Mask: [LOC]mm/BB/yyyy hh:mm S

Date/Time Format: MM/BB/yyyy hh:mm S

Label/Date/Email Date Masks

Standard Mask: MM-dd-yyyy

Month of Mask: MM/BB

Week of Date Mask: yyyy MM dd

MAKES - When cursor is NOT in control [vs.] FORMAT - When cursor IS in control (EXT masks)
DATE/TIME - Controls that show date and time [vs.] DATE - Controls that show only date

Cancel Reset All Masks & Formats Submit

The date formatting methods used in CC are software industry standards used by MS. For example, a date mask likeMM/dd/yyyy, would format August 23, 2022 into 08/23/2022.

Here are some sample date mask results for the date of August 23, 2022 . . .

 Date Masks	
Mask	Result
MMMM	August
MMM	Aug
MM	08
M	August 23
dddd	Monday
ddd	Mon
dd	23
d	8/23/2022
yyyy	2022
yy	22
y	August, 2022
[M]	M
Mask	Result
[Anything in Brackets]	Anything in Brackets

Each section in this dialog controls the settings for that specific date mask. They give users the ability to define a standard mask for all new records in that category. The date masks are sorted into 3 categories - masks for grids, masks for controls, and label date masks. Each mask uses a drop-down menu to select the desired format (i.e., MM/dd/yyyy). An example of each mask is shown in orange and changes in real time. To reset all masks to default settings, click on the [Reset All Masks & Formats] button at the bottom of the dialog. When done, click [Submit] to save changes.

GRIDS

Date/Time Masks & Formats for CONTROLS

Label/Auto Email Date

Date Mask: {LOC}mm/dd/yyyy

03/18/2021

Date Format: MM/dd/yyyy

03/18/2021

DateTime Mask: {LOC}mm/dd/yyyy h

03/18/2021 12:36 PM

DateTime Format: MM/dd/yyyy hh:mm t

03/18/2021 12:36 PM

Standard Mask:

Month of Mask:


Week of Date Mask:

MASK = When cursor IS in control | vs. | FORMAT = When cursor is NOT in control (EDIT mode)
DATETIME = Controls that show date and time | vs. | DATE = Controls that show only date

Reset All Masks & Formats

Number Masks in Calibration Control

Calibration Control (our **Calibration Management Software**) gives you the ability to format some of the more important date and number fields. Find the 'Masks - Number' dialog in the Advanced tab of the Program Options dialog.



Number Masks

Each section (Equipment ID Auto Increment, Calibration Certificate Number Auto Increment, etc.) controls the settings for that specific number mask. This gives users the ability to define a standard mask for all new records in that category. The first box controls the appearance of the mask, with a sequence of letters followed by a dash and zeros. This field should only contain zeros and no other numbers or this will cause an issue with the current increment. The 'Current Increment' is the next number CC will use in the mask when the next new record is created. An example of each mask is shown in orange and changes in real time. Activate or deactivate each mask feature by checking or unchecking the 'Enable' checkbox.

The Calibration Certificate Number Auto Increment textbox at the bottom of the dialog is an example of what the PDF format will look like when saved. Select (highlights blue) and then right-click at least one Equipment record. Then, click Print PDF.

Reserved Characters

In the Maintenance Activity Number Increment section, notice the **[M]** that is used in the Mask field. The character **M** is a reserved character for date masks which represents the Month. If you'd like to use a reserved character(s), simply placing brackets around the character allows you to use it in your number mask. Not sure what else is a reserved character? Take a look at the Date Masks table in the Date Masks help topic to learn more.

The screenshot shows a dialog box titled "Number Mask Options" with a grid of settings for different number masks. Each setting includes a label, a text input for the mask, a numeric input for the current increment, a preview of the next value, and an "Enable" checkbox.

Section	Mask	Current Increment	Next Value	Enabled
Equipment ID Auto Increment	EID-0000	1	EID-0001	Yes
Calibration Certificate Number Auto Increment	CC-0000	1	CC-0001	Yes
Project Number Increment	P-0000	1	P-0001	Yes
Work Order Number Increment	WO-0000	1	WO-0001	Yes
Maintenance Activity Number Increment	[M]A-0000	1	MA-0001	Yes
Uncertainty Study Number Increment	US-0000	1	US-0001	Yes

At the bottom, there is a "Calibration Certificate Number Auto Increment" section with a text input containing a complex mask: `Certificate_{@Certificate}_for_{EquipmentID}_{@printdat` and a "Mask Sample" button. Below this are "Cancel" and "Submit" buttons.

Leading Zeros Mask

Leading zeros can be important in number formats if all numbers require the same minimum number of characters. Instead of a counting number sequence of 1, 2, 3, etc. an LZ Mask of '000' results in 001, 002, 003, etc.

In the case where a counting number has more characters than the mask, the counting number wins. For instance, with a mask of 000 and a number of 5348, the result will be 5348 even though the mask has only 3 characters.

Calibration Measurement Group Mask

The Mask field in the Measurement Group of the Calibration dialog sets the leading and trailing zeros of all the numbers in the given group. Unlike the Leading Zeros Mask, the number of leading and trailing places of these fields are strictly limited by the mask. For instance, a mask of 00.000 will only allow two places to the left and 3 places (rounded) to the right.

Therefore, if a number of 100 or higher is required, the mask must be modified to include three zeros to the left 000.000.

Group Number: _____ 1 Group Name:

Tolerance: Tolerance (+/-) As Found:

Plus (+): _____ 0.0010 As Left:


Minus (-): _____ 0.0010 Units:

Masks: 0000.0000 00.0000

	Nominal	Upper	Lower	Found
▶	0.0500	0.0510	0.0490	0.0
	0.2000	0.2010	0.1990	0.2
	0.5000	0.5010	0.4990	0.5
	0.7500	0.7510	0.7490	0.7
	0.9500	0.9510	0.9490	0.9

Troubleshooting Number Masks

If you are receiving an error message when creating a new record, your number masks may not be set up properly. Go to the Options dialog > Advanced tab > Masks - Number. Ensure the mask field contains no other numbers besides zeros (e.g. EID-0000). Then, enter the next increment number in the Current increment field (e.g. 25).




Procedural Steps and Groups

Step-by-step Procedures for Worksheets

The Procedural Steps feature allows users to define step-by-step procedures that appear on Calibration Worksheets or

within custom reports in Calibration Control (our **Calibration Management Software**). Add multiple procedure steps to form a Procedural Group. If you don't see Procedural Groups or Steps in your Data Grids tab, follow along with the FeatureVisibility help topic to learn how to show those features.



Concept

Each *Step* in a Procedure is a Procedural Step while each *Procedural Group* (i.e., step-by-step procedure) is a collection of Procedural Steps in a specific order. To eliminate the need for duplication, each Procedural Step can be linked to multiple Procedural Groups.

Procedural Groups can be linked to Equipment Type, Model Number, and Equipment records. When linked, the Steps for a given Group will appear on the default Calibration Worksheet for Procedural Steps.

Procedural Groups

Create a new Procedural Group record to represent a new Procedure.

A screenshot of the 'Add Procedural Group' dialog box. The dialog has a title bar with a question mark and a close button. It has four tabs: 'Procedural Group', 'Procedural Steps', 'Custom', and 'Meta'. The 'Procedural Group' tab is active. It contains the following fields: 'Group Name' with the value 'Calibration', 'Code' with the value 'CAL', 'Group Type' with a dropdown menu set to 'Calibration', and an 'Active' checkbox which is checked. There is a 'Notes' field with the placeholder text 'Notes go here.' and a scroll bar. At the bottom, there are 'Cancel' and 'Submit' buttons.

- **Group Name:** Name that represents a collection of Procedural Steps
- **Code:** Procedure Group Code, if needed
- **Group Type:** Used to classify the group
- **Active:** Currently in use (active record)
- **Notes:** Keep track of additional details about the group

Procedural Steps

Create a Procedural Step record from its main grid or from within the Procedural Steps tab of the Procedural Group dialog.


- **Description:** Short description of the actual step (action)
- **Code:** Procedural Step Code, if needed
- **Step Type:** Assigning a Step Type is a good way to categorize Steps
- **Response Type:** While currently not in use, this field will designate what kind of response (in the software) will be required from the user
- **Active:** Procedural Step is still in use when checked
- **Notes:** Add any additional notes required for a Step

Add Procedural Steps to a Procedural Group

The Procedural Steps tab of the Group dialog contains the steps linked to the Group. Use the buttons at the bottom of the grid to:


- **[Add]** a new Step
- **[Link]** an existing Step
- **[Edit]** the selected Step
- **[Unlink]** the selected Step from the current Group without deleting it
- **[Delete]** the selected Step from ALL Groups

To add sequence numbers to the procedural steps simply type the desired number into the Sequence column of the 'Procedure Steps' panel grid.




Linking Procedural Groups to Equipment

Procedural Groups can be linked to the Equipment Type, Model Number, or directly to the Equipment record. Therefore, any Equipment record linked to an Equipment Type or Model Number will also be linked to that Type or Model Number's Procedural Group. In the instance that a different Procedural Group is linked to all three records (i.e., Type, Model, & Equipment), only the highest priority link will be displayed on the Calibration Worksheet. The record priority is **(1)** Equipment record, **(2)** Model Number, and **(3)** Equipment Type.



Equipment Record

Model Number




The 'Edit Model' dialog box contains the following fields:

Model	Uncertainty	Overrides	Custom	Meta
Model Number:	505-675	Measure Template:	TEMP-003	
Description:	0-6" Dial Caliper	Procedure Group:	Calibration	
Manufacturer:	ARC Manufacturing	Size/Range:	0-6"	
Equipment Type:	Caliper		<input checked="" type="checkbox"/> Active	
Notes:	Notes go here.			

Buttons: Cancel, Submit

Equipment Type



The 'Edit Equipment Types' dialog box contains the following fields:

Type	Overrides	Meta
Type Name:	Micrometer	
Type Code:	MC	Calibration Template: TEMP-003
Size/Range:	0-6"	Procedural Group: Calibration
Notes:	Notes go here.	

Buttons: Cancel, Submit

Calibration Worksheet

Calibration Worksheet

Equipment Information

Equipment ID:	Sample 101	Site:	California Plant
Description:	Pressure Tank	Department:	Production
Model No.:	505-716	Location:	Final Inspection
Manufacturer:	NEWPORT	Custodian:	Mick Moose
Size / Range:	0-100 PSI	Technician:	Dawn Ducking
Serial No.:	SN 101	Status:	Undefined

Calibration Summary

Calibrated:	01/18/2016	Frequency:	Temp:	As Found:
Next Cal:	01/18/2017	Humidity:	Result:	

Remarks: _____

Measurement Group 1

<u>Tolerance</u>	<u>As Found</u>	<u>Result</u>	<u>Procedures</u>
± 2% in			SP-1001, SP-1002

Desc.	Nominal	Limits		As Found			As Left (Cal Status)		
		Upper	Lower	Actual	Error	Result	Actual	Error	Result
	10.0000	10.2000	09.8000						
	50.0000	51.0000	49.0000						
	90.0000	91.8000	88.2000						

Notes: _____


Procedural Steps

Sequence	Description
1	Inventory All Parts Present
2	Inspect for Shipping Damage
3	Replace Batteries
4	Perform Functional Test Per Procedure RI-1013
5	Repair As Necessary
6	Ensure Chassis Connections Are Firm
7	Wipe Down Chassis
8	Complete Activity Report
9	Supervisor Approval

Approved By: _____ Date: _____

Set the Cal Worksheet with Procedural Steps as the Global Default Report

Within the Program Options dialog, go to the Reports tab and click on the 'Cal Worksheet' hyperlink. Select the[CalibrationWorksheetProcedureSteps.rpx] file and click [Open].




Override Instructions at Equipment Level

To override the procedural steps at the equipment level, navigate to the Overrides tab of the Edit Equipment dialog for a piece of equipment. Click on the 'Cal Worksheet' hyperlink and choose the [CalibrationWorksheetBlankProcedureSteps.rpx] file. Click [Open], then [Submit].

Required Fields

Available in Version 9.1 and higher, you can define which dialog fields must require data input.

Required fields in a dialog are highlighted light pink (by default) when the required field is empty or blank, and highlighted light green (by default) once the required field is populated. An error message will display when a required field is missing data when you attempt to save that record. Follow these instructions to set up your Required Fields.



Make Procedural Groups feature Visible

By default, the Procedural Groups grid is hidden from the toolbar. You'll need to make it visible by going to:

File > Options dialog > Advanced tab > Visibility > Data Grids tab

Un-check the checkbox for **Procedural Groups** in order to make the feature visible. Restart Calibration Control to see the Feature Visibility changes reflected in the application.

Create a New Procedural Group

Open the Procedural Groups grid (located under the Data Grids tab in the Work Flow section) and add a New Record.

1. **Group Name:** Enter a unique Group name (i.e. Required fields for Equipment)
2. **Code:** The Code field is where you will enter the actual name of the dialog. Below are some dialog names for your reference.
 - Equip
 - EquipmentSystem
 - EquipmentType
 - Calibration
 - MeasurementGroup
 - Attachment
 - Document
 - Note
 - Department
 - Person
3. **Group Type:** Required Fields
4. **Active:** Ensure the Active checkbox is marked.
5. **Notes:** This is where you'll add your required field(s) information. The syntax for a required field is:

**Field Name|Custom field name|Tooltip message
(Custom error message)**

- **Field Name:** This is the actual field name that is in your database. You can open up your database to find the specific field name. Below is a list of some Equipment dialog field names for your reference:
 - SerialNumber
 - ModelNumber
 - CompanySite
 - DepartmentGuid

- ModelDesc
 - SizeRange
 - ManufacturerGuid
 - EquipmentTypeld
 - Alternateld
 - LocationGuid
 - Custodian
 - StatusGuid
 - CertificateNumber
 - Technician
- **Custom Field name:** This custom field name will be used in your error message instead of the actual fieldname.
 - **Tooltip message:** Optional. This message will appear when hovering over a required field.
Note: You can add more than one required field, by adding required fields on the next line.
 - **Custom error message:** Optional. You can add a custom message that will appear above the default error message when a user attempts to submit a form with a missing required field. You can do this by adding an additional line with your message inside parentheses.

Custom Highlight Colors

The Required fields feature uses Highlight colors to indicate whether a required field has been filled out. By default, a required field is highlighted light green when it has a value or light pink when it is empty. Optionally, you can change these colors in **Program Options > Advanced tab > Highlight Colors > Good Response & Bad Response**.

Neutral:	White	Past Due (Cal):	LightPink
Calibration Standard:	LightGreen	Past Due (Standard):	LightPink
Pass:	LightGreen	Past Due (Maint):	LightPink
Fail:	Pink	Soon Due (Cal):	Yellow
Limited Calibration:	Yellow	Soon Due (Standard):	Yellow
Out of Tolerance:	Red	Soon Due (Maint):	Yellow
In Calibration/Maint:	Gray	Good Response:	LightGreen
		Bad Response:	LightPink

Buttons: Cancel, Reset, Save

Get Ape Support

If you need assistance setting up the Required Fields feature, contact us and we can do this together.


Publish a Custom Website

Publish Website of Equipment Grouped by Departments

Instantly publish a static Website of all equipment in Calibration Control to your company network so everyone has access to online lists of equipment sorted by Department. Schedule the website to be refreshed regularly using the Auto Notify Utility.

One-Click Website Publishing (Create)

Select the 'Publish Website' menu item in the Utilities tab.



Website will be automatically published from the available Equipment Records. A dialog box will pop up with an option to go to the Home Page, which has links to all Department pages. Here is an example of the type of information available with the Publish Website feature.

Manufacturing


[back](#) | (Last updated: Thursday, June 20, 2019 2:02 PM)

ID	ModelNumber	Description	Serial	Status	Frequency	Last Cal	Due Cal
Test Area							
SAMPLE-006		10-32 UNF-3B	1008	Active	Yearly	7/3/2018	7/3/2019
SAMPLE-025	RGS263e	Radius Gage Set	N27627	Active	Yearly	5/3/2019	5/3/2020
SAMPLE-040	DPM	Depth Point Micrometer	2060	Active	Yearly	4/17/2019	4/17/2020
SAMPLE-043	010-MM1	Micrometer	34247-1093	Active	Yearly	4/16/2019	4/16/2020
SAMPLE-049	DTM	Deep Throat Micrometer	2021	Active	Yearly	4/15/2019	4/15/2020
SAMPLE-056	BM-64	Blade Micrometer	109	Active	Yearly	4/9/2019	4/9/2020
SAMPLE-068	DRPI-57	Drop Indicator	2098	Active	Yearly	4/17/2019	4/17/2020
SAMPLE-097	IDR-87abd	Indicator	2247	Active	Yearly	3/4/2019	3/4/2020
SAMPLE-114	MM-67	Micrometer	345	Active	Yearly	4/11/2019	4/11/2020
SAMPLE-118	IDR_7bf8	Indicator	2194	Active	Yearly	4/11/2019	4/11/2020
SAMPLE-123	CMMA	CMM ARTIFACT	187	Active	Yearly	9/19/2018	9/19/2019
SAMPLE-125	SM-079e	Scissor Micrometer	7126751494	Active	Yearly	10/3/2018	10/3/2019


Auto Publication with Auto Notify Utility

Publish (refresh) the website automatically at a set schedule using the Auto Notify utility. Begin by selecting 'Auto Notify' from the Utilities group of the Calibration tab.

From the Auto Notify grid, open the WebSite_General record (e.g., double click to open).




Within the Notification tab, **(a)** change the Frequency (i.e., 'How Often'), **(b)** the start date and time of day (i.e., 'Date & Time'), **(c)** ensure the Is Active checkbox is checked, and **(d)** submit by selecting the [OK] button.



Change Location Where Site is Published

The default location where the site is published is the Settings Folder. Change the Website publication directory to a network folder so other network users have access. Go to File > Options and select the 'Advanced' tab. Click 'Folder & Path' and select the desired network location in the Location of Auto Generated Website section.



Measurement Uncertainty Studies

Perform Uncertainty Analysis of a Measurement

The Measurement Uncertainty Studies grid (formerly called Uncertainty Budgets) in Calibration Control (our **Calibration Management Software**) is a simplified budget where Ape Software assumes no responsibility or liability for any unintentional errors. It is the User's responsibility to ensure calculations and data entry are accurate and to confirm the suitability for any particular purpose.

This budget is useful in situations where **(a)** each Component of Uncertainty is independent of the others, **(b)** the major Components are estimated with a high degree of accuracy, and **(c)** measurement results are measured directly or calculated with a linear formula using only multiplication or division (e.g., $W=V \times I$).

Note that each Measurement Uncertainty Budget is conducted for a single Nominal measurement for a unique piece of

Equipment.

Uncertainty Study Dialog

- **Study ID:** Add a name or number as the ID for an uncertainty study.
- **Date:** Date the Uncertainty Budget calculations were performed.
- **Nominal:** Target value for the overall Uncertainty Budget analysis.
- **Measure Units:** Measurement unit (e.g., inches, degrees, etc.) of the Nominal value for the Uncertainty Budget analysis.
- **Significant Digits:** Number places to the right of the decimal.
- **Description:** Description or title of the Uncertainty Budget analysis.
- **Notes:** Any notes relevant to the Uncertainty Budget analysis.
- **Combined [Standard] Uncertainty (uc):** Summation in quadrature (i.e., root sum of the squares) of all the Sources of Uncertainty in the grid.

$$uc = \sum_{k=1}^n u_k^2$$

- **Coverage Factor (e.g., k=2):** Used in the calculation of the Expanded Uncertainty and has a default value of k=2. The Coverage Factor is a confidence level of the interval accounts for the variation in the estimate of the residual standard deviation, and is based on the assumption that the random errors have a normal distribution.
- **Expanded Uncertainty (U):** The Expanded Uncertainty (U) assures a high level of confidence by re-scaling the Combined Uncertainty (uc) expressed with a different confidence level, e.g., 95 percent. The re-scaling is accomplished by multiplying by the coverage factor (i.e., $U=k \times uc$). In many cases, the coverage factor may be 2, which is typically used to approximate a 95% level of confidence.
- **Expanded Uncertainty Rounded (Uc):** The expanded uncertainty (Uc) rounded to two significant digits.

Study ID: Test 003 Description: Sample Uncertainty Study

Date: 10/10/2019 Notes: Notes go here.

Nominal: 0.10

Measure Units: cm

Significant Digits: 2 Active

Source	Va	Units	Type	Distribution	Divisor	Freedom	Coef	Std Un
Sample Uncertainty Component	0.1	cm	A (Statistical)	Norm @ 99%	2.58	0	0	0.04

Combined Uncertainty: 0.04 Coverage Factor (K*): 2 Expanded Uncertainty: 0.08

Expanded Uncertainty (Rnd): 1.0

Buttons: Calculate Uncertainty, Cancel, Submit

Uncertainty Detail Dialog

Use the Detail dialog to describe each source of uncertainty (i.e., Component of Uncertainty) to be included in the overall Measurement Uncertainty Budget.

- **Source Name:** Descriptive name for the Source of the Component of Uncertainty, which will be combined with other Components to create a single Uncertainty Budget.
- **Value (a):** Numeric value related to the Unit of Measure.
- **Unit of Measure:** Units measured in such as inches, Fahrenheit, mph, etc.
- **Type A/B:** Method used to estimate uncertainty either Type A, if estimated statistically, or Type B, if inferred through non statistical methods such as historical measurements, manufacturer specifications, calibration certificates, published data, mathematical formulas, or general experience.
- **[Probability] Distribution:** The Probability Distribution selected dictates the Divisor (next field). For example, a Normal 1s distribution has a Divisor of 1 where a Normal 2s distribution has a divisor of 2. The following values for each Probability Distribution are defaulted into the Divisor field, which can be overridden.
 - No Divisor
 - Normal 1s = 1
 - Normal 2s = 2
 - Normal 3s = 3
 - Normal @ 99% = 2.58
 - Rectangle = 3 = 1.7321
 - Triangle = 6 = 2.4495
 - U-Shaped = 2 = 1.4142
 - Rectangle = 12 = 3.4641
- **Divisor (d):** The numeric value related to the Probability Distribution chosen. This value is automatically set when choosing the Probability Distribution (previous field) but can be edited independently after choosing the Distribution.
- **Degrees of Freedom (v):** For Type A uncertainties, use the d.f. for the corresponding standard deviations. For Type B sources refer to calibration certificates or published reports. In some instances, where the standard deviation must be estimated from scientific judgment or partial data, an infinite d.f. is normally assumed.
- **Sensitivity Coefficient (c):** Used to describe how the Component of Uncertainty contributes to the Combined Uncertainty when a simple functional relationship does not exist between the input quantities and the measurement result. For example, the unit of measure for the study may be in length units where an Uncertainty Component maybe expressed in temperature or angular units. Use this field to change the weighting of the Component from the default of 1 (one).
- **Standard Uncertainty (u):** Calculated by dividing the Value (v) by the Divisor (d) and then multiplying by the Sensitivity Coefficient (c):
$$u = \frac{v}{d} \times c$$
- **Assumptions Explanation:** Use this field to Explain any Assumptions made for the current Uncertainty Component.

Learning about Measurement Uncertainty

Here are a few online resources to get started with learning about Measurement Uncertainty.

- How to Calculate Uncertainty - *Dr. Jody Muelaner*
- Uncertainty Budgets and Sensitivity Coefficients - *NIST, Information Technology Laboratory*
- Calculating an Uncertainty Budget for a Measurement - *wikiHow*
- Example of Uncertainty Budget - *NIST*
- Policy on Estimating Measurement Uncertainty for Construction Materials & Geotechnical Testing (P103d) - *The American Association for Laboratory Accreditation*
- A Beginner's Guide to Uncertainty of Measurement - *Stephanie Bell, National Physical Laboratory*
- An Introduction to Expressing Uncertainty in Measurement - *Mr. Ouellette, National Research Council Canada*

Here are a few spreadsheet examples of Measurement Uncertainty Budgets:

- Uncertainty Budget Template - *NIST, Laboratory Metrology Program*
- Simplified Uncertainty Budget Template (xls) - *National Research Council Canada*
- Measurement Uncertainty Budget Template - *American Society of Crime Laboratory Directors / Laboratory Accreditation Board*

References

- Degrees of Freedom - *Wikipedia*



THUM Temperature and Humidity Device

Automatically Record Temperature and Humidity


Calibration Control is compatible with the Temperature and Humidity USB Monitor (THUM) sold by Practical Design Group. With a THUM device plugged in and set up on a computer, Ape can automatically retrieve temperature and humidity data for calibration events.



Setup

Follow these instructions to integrate the THUM with Calibration Control:

1. **Setup THUM:** Follow the setup instructions that arrive with the THUM device.
2. **Move THUM Database:** The default location of the THUM database (thum.mdb) is in the Program Files folder. This will only work if the signed in user has read/write access to the Program Files folder. Therefore, always move the database out of the Program Files folder to a location where everyone has read/write access, like the CC Settings folder.
 - a. Move the THUM database by opening the THUM application and selecting 'Show Options' from the Options drop-down menu.
 - b. Update the 'Database location' field to reflect the new location, like the CC Settings folder.
3. **Tell Ape Where THUM Is:** Open the Folders tab of the Options dialog in Calibration Control and click the Edit link in the 'Location of THUM . . .' field. Navigate to the same location as above.
4. **Test Setup:** Test the setup by creating a new Calibration Event record (Calibrations tab of Equipment Edit dialog) and confirming that the Temp and Humidity fields are automatically populated.



Troubleshooting

- **Insufficient Permissions:** Attempt to start the THUM Monitor software and it reports that the current user has insufficient permissions to make changes. Resolve this by running the THUM monitor as an Administrator. Do this by right-clicking on the THUM icon and selecting "Run as Administrator".
- **THUM Service Will Not Start:** If the THUM service cannot be started from the THUM Monitor, one or more of the following steps should fix the problem. For a complete reset, follow each step:
 1. **Close THUM Monitor:** Resets the software.
 2. **Delete THUM Database:** Delete the "thum.mdb" file that the THUM Monitor is pointed to. Ensure this is not the database in the Program Files folders, which should not be used for live data. Deleting the database resets the database file to ensure no corruption and that (thum.mdb) is in the proper version format for the THUM service. Therefore, when opening the "thum.mdb" file with MS Access, do not allow Access to convert the "thum.mdb" file to a newer version of MS Access.
 3. **Reset THUM Device:** Ensure the electronic THUM device itself is reset by unplugging it from the computer for at least 5 seconds and plugging it back in.
 4. **Restart THUM Monitor:** Remember to start the THUM Monitor with Administrator Privileges, if necessary.
- **Still Need Help?** Contact us if the above steps do not work and we can do a screenshare together.

Best Label Printers for Ape Software

Which label printers will work with Ape's Software?

Calibration Control is compatible with the **PC-connectable Brother P-touch label printers.**

Brother P-touch Label Printers

The printer we recommend most is the PT-P900W (replacement to the PT-9700PT), which is a high volume desktop unit intended for industrial environments. The PT-P900W will print up to 32mm (1.26 inch) width tape. If the PT-P900W is too expensive, maybe the portable PT-P700 would be a better starting choice. The PT-P700 uses label tape up to 24mm (1 inch) in width. Otherwise, either printer will work with Ape Software applications.

You can find the PT-P900W and the PT-P700 at Amazon.com and many other online retailers.

TZe Label Tape

All of the P-touch label printers use the same label tape (model number starts with TZe). This tape is resistant to temperature, UV, scratching, and spills.

Here are some of the common sizes along with their Amazon links:

- 6 mm (1/4 inch) Black on White (TZe211)
- 9 mm (3/8 inch) Black on White (TZe221)
- 12 mm (1/2 inch) Black on White (2 pack - TZe2312PK)
- 18 mm (3/4 inch) Black on White (TZe241)
- 24 mm (1 inch) Black on White (TZe251)
- 36 mm (1 1/2 inch) Black on White (TZe261)

Extra Adhesive Label Tape

Brother also sells label tape with extra adhesive for high-use environments or extra smooth surfaces. Before choosing the extra adhesive, we suggest trying the normal adhesive formula first (above labels) because the additional adhesive can gum up label printers faster than the normal level of adhesive.

- 6 mm (1/4 inch) Black on White (TZeS211)
- 9 mm (3/8 inch) Black on White (TZeS221)
- 12 mm (1/2 inch) Black on White (TZeS231)
- 18 mm (3/4 inch) Black on White (TZeS241)

Here are a couple of links to B&H Photo for Cleaning

- Tape: 36 mm (1 1/2 inch) Good for PT-P900W
- 18 mm (3/4 inch) Good for PT-P700

Free Label Editing Software and Drivers

Free Brother Printer Drivers and Label Editing Software

Ape Software uses Brother P-touch label printers for the highest quality laminated thermal transfer label printing.

Use the following quick links to find the drivers and software from the Brother support site.

Printer Drivers

- Download the PT-9800PCN Printer Drivers
- Download the PT-9700PC Printer Drivers
- Download PT-2430PC Printer Drivers
- Download PT-P700 Printer Drivers
- Download PT-P900W Printer Drivers

P-touch Editor - Label Editing Software

Download the Brother P-touch Editor for editing the labels that come with Calibration Control or creating your own.

Labels Quick Start Instructions

Sample Labels

Check out the sample labels for Calibration, Asset, and CNR (Calibration Not Required). They are a good resource for determining which labels to use.

Important Points

Keep in mind when printing labels with Calibration Control:

1. **Brother Label Printers:** Calibration Control will only print labels using one of the Brother P-touch (PT) label printers. Click on the link to view a list of compatible models.
2. **Printer Drivers:** The printer drivers MUST BE INSTALLED in order to use a label printer with Calibration Control, these drivers are accessible with the CD included with the label printer (or downloaded from Brother). Follow the instructions *before* plugging in and turning on the printer.
3. **Label Files:** Calibration Control uses label layout files (*.lhx) which can be accessed by clicking on the File tab at the top left corner of Calibration Control and selecting [Open Files Folder]. The label files *must* reside in the Labels folder for Calibration Control to find them. Watch the video on Calibration Control folders to understand the folder layout.
4. **Label File Number Prefixes:** The 150+ label files that come with Calibration Control begin with the millimeter prefixes of 12 (1/2 inch), 18 (3/4 inch), or 24 (1 inch) to identify the width of label tape required.
5. **Set Default Labels:** Use the Options dialog to set the default labels for **Calibration Due**, **Asset**, and **Cal**

NotRequired. Be sure to select labels with the same width as the tape in your P-touch printer.

6. **Print a Label:** Open an equipment record and click the [Asset Tag] button in the bottom-right corner of the dialog.

Access All the Label Fields

Understanding Labels Fields in Calibration Control

Follow these guidelines to get the most out of printing labels with Calibration Control (our Calibration ManagementSoftware).

Need To Know:


1. **Label Basics:** Watch this label quick start video for an introduction to [the basics of label printing](#) in CalibrationControl.
2. **Files Folder Basics:** Watch this video for an overview of the Calibration Control folders.
3. **Sample Labels:** Check out these screenshots of sample labels included in Calibration Control for **Calibration**, **Asset**, and **CNR (Calibration Not Required)**. They are a good resource for choosing which labels you want to use or modify.
4. **Label Naming Requirements:**
 - Printing the equipment labels or **Asset** tags must have '**Asset**' in the label file name, (e.g. 24-Asset-05).
 - Printing **Calibration** labels must have '**DueCal**' in the label file name, (e.g. 18-DueCal-02).
 - Printing **Maintenance** labels must have '**DueMaint**' in the label file name, (e.g. 24-DueMaint-01).
 - Equipment **System** labels must have '**System**' in the label file name, (e.g. 36-System-01).
 - Due Cal Labels for **Calibration Not Required** Frequency must have '**CNR**' in the label file name, (e.g. 12-CNR-01).
 - Due Cal Labels for **Next Use** Frequency must also use '**CNR**' in the label file name for current version 9, (e.g. 18-CNR-09). In older version 8, it must have '**NextUse**' in the file name, (e.g. 24-NextUse-03).
 - **Note:** The format of the label file name begins with the millimeter (mm) size of your tape width.
5. **P-touch Editor:** The P-touch editor software that comes with P-touch printers (or download from Brother) is required to edit the 150+ labels that come with Calibration Control or to create custom labels.
6. **Editing Label Templates:** You can edit Calibration Control label templates in the free **P-touch Editor software**. The label fields below are the designated **Object Names** in the Textbox or Barcode Object Properties.
7. **Using Label Fields:** The label fields below correspond to the specific data fields used in Calibration Control. Copy the exact label field format when modifying any given template.
8. **Label Fields:** All of the available label fields are shown below.

Take a look at our Label Fields help topic to see a full list of label field names.

Sample Due Cal Labels


Below are the sample Calibration Due (Due Cal) Label templates included in Calibration Control. These labels can be edited and additional labels can be created as needed.

The number at the beginning of each label file refers to the width of the label tape (in mm) or height of the label as it's presented below. The number at the end of the label is simply an id.






18-DueCal-10



18-DueCal-11



18-DueCal-12




18-DueCal-13




18-DueCal-14



18-DueCal-15



18-DueCal-16



18-DueCal-17



18-DueCal-18



18-DueCal-19



24-DueCal-01



24-DueCal-02




24-DueCal-03




24-DueCal-04




24-DueCal-05




24-DueCal-06




24-DueCal-07




24-DueCal-08




24-DueCal-09



24-DueCal-10




24-DueCal-11



24-DueCal-12



24-DueCal-13




24-DueCal-14



24-DueCal-15



24-DueCal-16



24-DueCal-17



24-DueCal-18


Print Future Calibration Labels

Print Calibration Labels to Apply at a Future Date

Using Calibration Control, this feature is useful in a scenario where a Technician must calibrate Equipment at a remote site but does not have the ability to print labels while at that site. With this feature, a Technician can print a range of labels with a future Calibration date and apply those labels on site assuming the Unit Under Test (UUT) passes. Printing labels with this feature does not update the Calibration dates of the Equipment.

Printing Future Labels

Print Future Calibration Labels by selecting a group of records, right-clicking to get the pop-up menu, and selecting Print->Future Calibration Labels. Use the Filter Row in the Equipment grid to filter only the rows required.



Set the future calibration date and the appropriate technician for the labels.


Future Calibration Date

Calibration Date: 07/10/2019

Technician: Tina Williams

Cancel OK

Depending on the type of label selected, different information will be displayed in different styles. This label shows the Technician name, the Future Calibration Date, and the date it will be Due for Calibration. The customizable fields are shown in red for reference.



Sample Asset Labels

Below are the sample Asset Label templates included in Calibration Control. These labels can be edited and additional labels can be created as needed. The number at the beginning of each label file refers to the width of the label tape (in mm) or height of the label as it's presented below. The number at the end of the label is simply an id



12-Asset-01



12-Asset-02




12-Asset-03




12-Asset-04




12-Asset-05



12-Asset-06



12-Asset-07



12-Asset-08



12-Asset-09



12-Asset-10



12-Asset-11




12-Asset-12



12-Asset-13



18-Asset-01



18-Asset-02



18-Asset-03



18-Asset-04



18-Asset-05



18-Asset-06



18-Asset-07



18-Asset-08



18-Asset-09



18-Asset-10



18-Asset-11



18-Asset-12



18-Asset-13



18-Asset-14



Calibration Not Required Labels

Below are the sample Calibration Not Required (CNR) Label templates included in Calibration Control. These labels can be edited and additional labels can be created as needed.

The number at the beginning of each label file refers to the width of the label tape (in mm) or height of the label as it's presented below. The number at the end of the label is simply an id.



12-CNR-01




12-CNR-02



12-CNR-03




12-CNR-04



12-CNR-05




12-CNR-06



12-CNR-07



12-CNR-08




12-CNR-09



12-CNR-10




12-CNR-11




12-CNR-12



12-CNR-13




12-CNR-14




12-CNR-15




18-CNR-01




18-CNR-02




18-CNR-03



18-CNR-04




18-CNR-05



18-CNR-06



18-CNR-07




18-CNR-08



18-CNR-09



18-CNR-10



18-CNR-11



18-CNR-12



Chain Printing Labels

How to Chain Print Labels to save Label Tape

There are two ways to print labels using the Brother label printers and each utilizes a small portion of blank label tape used as a "grab-tab" to help peel away the label from its backing. The first method by default prints them individually with a small piece of "grab-tab" left on the end of each scored piece. The second, using this method, prints and scores them all in a row with one piece of "grab-tab" at the beginning of the strip.


Individually Printed Labels

By default, labels are printed individually with a small "grab-tab" at the front of each:



Chain Printed Labels

From the File tab at the far left of the ribbon menu, choose the 'Options' button and navigate to the 'Labels' tab. At the bottom of the dialog there is a 'Chain Printing (not for all models)' checkbox. Check this box to enable chain printing. Uncheck the box in order to go back to printing them individually.



The chain printed labels will look like this:



Labels & Reports Overrides


Override Global Defaults for Reports and Labels

By default, the Global (application-wide) settings in Calibration Control for file templates that are used for Labels and Reports are defined in Options > Advanced tab > 'Report & Label'.

In situations where the Global template is not acceptable, users can override the template in the Equipment, Site (Client), Model Number, or Equipment Type. For example, if a 12mm label is set at the Global level and a different size (larger or smaller) is required for a specific Model Number or Equipment, a different label template can be selected at the Model or Equipment level.

Override Settings

In the Overrides tab of the Equipment dialog, you can override Labels and Reports. The other dialogs that have Label and Report Override ability have similar fields. The template values can be cleared by clicking on the red [X] buttons.



Label Fields

- **Calibration Due:** Prints when the calibration frequency generates a due date and uses labels with the "Due Cal" string in their file name.
- **Asset:** File templates are identified as Asset labels by the "Asset" string in the file name.
- **Cal Not Required:** Prints in place of the Calibration Due label when the Calibration Frequency is "Cal NotRequired" and is identified as a CNR label by the CNR string in the file name.
- **Next Use:** Print in place of the Calibration Due label when the Calibration Frequency is "Next Use" and the labels are identified by a "NextUse" string in their file name.

Report Fields

- **Calibration Cert:** Template used when printing a Calibration Certificate.
- **Due Cal Report:** Report template used when printing the Calibration Due Report.
- **Cal Worksheet:** Calibration Worksheet that prints when the Equipment record has a Calibration Template or a previous Calibration Event.
- **Cal History:** Calibration History of an Equipment record.
- **OOT Investigation:** Designates the Out-Of-Tolerance (OOT) Investigation worksheet used when Equipment received for calibration is found OOT.


Overrides apply to Labels and Reports via four paths of priority.

1. **Direct:** Override directly through an Equipment Record on the 'Overrides' Tab. An override through an EquipmentRecord supersedes all other methods.
2. **Company:** Override labels or reports through a Company record connected to an Equipment Record. This override holds priority over those below.
3. **Model Number:** Override labels or reports through a Model Number record connected to an Equipment Record. This override holds priority over those below.
4. **Equipment Type:** Override labels or reports through an Equipment Type record connected to an Equipment Record. This override holds priority over those below.
5. **Equipment Type — Model Number:** Override labels or reports by linking an Equipment Type to a Model Number, then linking the Model Number to the Equipment Record.
6. **Options Dialog:** Override labels or reports through the Options Dialog.

Print Reports

Printing the Standard Reports is Easy

Click on the [Print] icon in the Calibration tab of the ribbon menu of Calibration Control (our **Calibration Management Software**), then select which report to print. If the report requires additional information (i.e., parameters, due date, or other filters) to print, enter the required parameter and then click the [Submit] button to




generate the report.

Ape software also offers a Report Designer to modify any of the existing reports or to create custom reports.

Calibration Due Report

Quick Access to the Due Cal Report

Print the Calibration Due Report by selecting the Due Cal Report option from the Calibration tab or from the context menu within the Equipment grid under Print options.



The following Due By dialog will appear.

A screenshot of a dialog box titled 'Report Parameters'. It has a close button (X) in the top right corner. The main field is labeled 'Calibration Due By:' and contains the date '06/17/2019'. At the bottom left is a 'Cancel' button, and at the bottom right is a 'Submit' button.

If a date range is required, click on the [PRINT] button in the Calibration tab of the ribbon menu to find the Calibration Due Date Range report. The Date Range report (file name CalibrationDueDateRange.rpx) can also be set as the default report. The range dates are set by default to the first date of the current month and the last date of the current month. See the Program Options help topic instructions on how to change the default report.

A screenshot of a dialog box titled 'Report Parameters'. It has a close button (X) in the top right corner. There are two date fields: 'Cal Due From:' set to '06/01/2019' and 'Cal Due To:' set to '06/30/2019'. At the bottom left is a 'Cancel' button, and at the bottom right is a 'Submit' button.

Note: The default versions of these Due Cal reports do not include equipment with non-calibration Frequency codes, like 'Cal Not Required'. Also, equipment with Status Codes marked as Hidden will not be displayed.

Due Cal: Print Preview

Here is a sample Calibration Due Report ready to print.

Calibration Due

Equipment ID	Model	Description	Serial	Status	Last Cal	Next Cal	Department
SAMPLE-006		10-32 UHF-3B	1003	Active	10/20/2015	10/20/2016	Manufacturing
SAMPLE-010	PS-03006	Pin Set		Active	9/2/2015	9/2/2016	Quality Assurance
SAMPLE-011	PS-03006	Pin Set		Active	9/2/2015	9/2/2016	Quality Assurance
SAMPLE-012	PS-03006	Pin Set		Active	9/2/2015	9/2/2016	Quality Assurance
SAMPLE-015	FRC	Fuel Resistant Coating	2145	Active	10/20/2015	9/2/2016	Shipping
SAMPLE-017	CMH-421	CMH	185	Active	10/22/2015	10/22/2016	Inspection
SAMPLE-018	20MG-KT		0015	Active	11/17/2015	11/17/2016	Inspection
SAMPLE-023	DRP-67	Drop Indicator	2119	Active	10/18/2015	10/18/2016	Shipping
SAMPLE-023	IDMS	Indicating Micrometer	2165	Active	9/19/2015	9/19/2016	Production
SAMPLE-028	MG85-878	Master Gauge Block Set		Active	10/7/2015	10/7/2016	Quality Assurance
SAMPLE-031	14 - 28 UNJ-3B		2930	Active	10/6/2015	10/6/2016	Machine Shop
SAMPLE-063	DBS-72	Dial Sore Gauge	2059	Active	12/10/2015	12/10/2016	Stockroom
SAMPLE-065	PTM8111	Point Micrometer	1027827	Active	11/16/2015	11/16/2016	Machine Shop
SAMPLE-071	DBS-55	Dial Sore Gauge	2023	Active	11/12/2015	11/12/2016	Production
SAMPLE-085	BI-67	Block Ink	103	Active	11/16/2015	11/16/2016	Stockroom
SAMPLE-098	RC-4005	Reconditioner	M041734	Active	11/9/2015	11/9/2016	Production
SAMPLE-099		4 Kg Weight	101	Active	11/20/2015	11/20/2016	Stockroom
SAMPLE-103	H-CLPR-22	Caliper	140	Active	11/13/2015	11/13/2016	Stockroom
SAMPLE-104	DM-11	Depth Micrometer	2049	Active	12/29/2015	12/29/2016	Production
SAMPLE-111	H-CLPR-22	Caliper	140	Active	11/19/2015	11/19/2016	Production
SAMPLE-116	IDR663	Indicator	2190	Active	11/20/2015	11/20/2016	Quality Assurance
SAMPLE-119	DR_7018	Indicator	307	Active	11/20/2015	11/20/2016	Production
SAMPLE-123	CM84	CHMARTIFACT	187	Active	12/23/2015	12/23/2016	Manufacturing
SAMPLE-161	32-MH	Micrometer	301	Active	11/11/2015	11/11/2016	Stockroom
SAMPLE-162	32-MH	Micrometer	301	Active	11/11/2015	11/11/2016	Production
SAMPLE-163	32-MH	Micrometer	303	Active	11/11/2015	11/11/2016	Manufacturing
SAMPLE-164	32-MH	Micrometer	324	Active	11/11/2015	11/11/2016	Quality Assurance
SAMPLE-165	32-MH	Micrometer	305	Active	11/11/2015	11/11/2016	Shipping

12/27/2016 10:31 AM

[Calibration Due]

Page 1 of 3

Calibration Worksheets


Record Calibration Data While Away from a Computer

Calibration Worksheets are a convenient way to record Calibration Event data without being near a computer. The Worksheets in Calibration Control have the proper measurements pre-determined to make calibration data easier to record. The Worksheets are easily configured for the needs of the organization.

When printing a Calibration Worksheet:

- If the Equipment has a Calibration Template attached to it, the Worksheet will replicate the Template structure on the Worksheet so all the technician needs to do is fill in the blanks.
- If no Calibration Template is attached to the Equipment, Calibration Control will use the measurement structure from the most recent calibration of the Equipment.
- If neither a Template nor past calibration exists, Calibration Control will print a blank Worksheet with several lines for recording measurement data.

Right-click for the Context menu to print.



The worksheet will look like this in the print preview. Notice that some of the fields are already filled-in with information from the record:

Calibration Worksheet

Equipment Information

Equipment ID: SAMPLE-016 Site:
 Description: Caliper Department: Inspection
 Model No: Location: Shelf
 Manufacturer: Custodian: David Wilcox
 Size / Range: 0-12" Technician:
 Serial No: 0796 Status: Active

Calibration Summary

Calibrated: 03/19/2019 Frequency: Yearly Temp: _____ As Found: _____
 Next Cal: 03/19/2020 Humidity: _____ Result: _____
 Remarks: _____

Measurement Group 1


Tol Type Units As Found Result Procedures
 Tolerance (+/-)

Desc	Nominal	Limits		As Found			As Left (Cal Status)		
		Upper	Lower	Actual	Error	Result	Actual	Error	Result
00.0000		00.0010	-00.0010						
02.0000		02.0010	01.9990						
05.0000		05.0010	04.9990						
08.0000		08.0010	07.9990						
11.0000		11.0010	10.9990						

Notes: _____


Approved By: _____ Date: _____

The blank worksheet looks like this when neither a Template nor past calibration exists. All of the fields are blank to accommodate necessary information.




Define Report Data Source

If you click on the [Edit Data Source...] link in the bottom-right corner of the Report Designer a blank Report Data Sourcedialog will appear. The Connection String field defines the connection to the database being used. The Query field defines the fields and records needed from a specific table and in a defined order.




Connection String - Choosing the Provider

The easiest way to change the connection string is to select the Connection String icon from the ribbon menu. This will show the same connection string the Ape application is using to connect to the database. Click the [Copy to Clipboard & Close] button and paste the result into the Connection String field of the Report Data Source dialog.



Otherwise, build the connection string by clicking the [Build] button to the right of the Connection String text box and the Data Link Properties will be shown. Click on the Provider table if it is not already displayed. Assuming you want to connect to an Access database (i.e., apecal.mdb), select the 'Microsoft Office 12.0 . . .' provider and click the [Next] button. If connecting to an instance of SQL Server, select the SQL Server OLE DB Provider.

If this does not work, ask your database administrator (DBA) for assistance in choosing the correct provider.




Connection String - Setting the Connection

Continuing the assumption that you need to connect to the apecal.mdb file, paste the path to the apecal.mdb file in the Data Source field. With help on finding your database, read the Locating Database File help topic. The entire path will look something like this:


C:\Users\Public\Documents\Ape Software\Calibration Control\apecal.mdb

After entering the path to the database, click the [Test Connection] button and a 'Test connection succeeded' response should show. If the connection is unsuccessful, repeat the above steps until it is.



Writing a SQL Select Statement

When the SQL SELECT statement has been created, it will look something like the following picture. Note that the Calibration Due field is set to a specific filter of less than 6/18/2019. This date is called a parameter (report variable), which can be changed with every printing. See the help topic on Report Parameters to learn how to add Parameters to SQL statements.



Adding Fields, Labels, and Report Info

In the image at the top of this page, I performed the following actions to create the Calibration Due report:

1. Dragged the bound fields (e.g., Equipment ID, Model & Description) from the right side of the page (Fields - Boundtree) to their current location in the Detail band and resized them to fit their contents.
2. Selected the Date fields and edited their properties (bottom-right corner of screen) so that the OutputFormat = 'M/d/yyyy'.
3. Dragged Label objects from the left side of the screen to locations above each field in the page header, resized them to fit their corresponding fields, bolded, and underlined them. The page title (Calibration Due Report) was created the same way.
4. Report Info objects were dragged from the left side of the screen to the left and right side of the footers. I clicked on each object and changed their properties (bottom-right) to . . .
 - Set the Format String property of the date/time field to '{RunDateTime:M/d/yyyy}'
 - Set the Format String property of the page number field to 'Page {PageNumber} of {PageCount}'
 - Clicked the right-align button (top of page) for the page number field.
5. Other minor tasks involving bolding, underlining, aligning, and positioning were made to make the report look the way I wanted it to look.

Menu Settings

Assuming the custom report is stored in the Ape application Reports folder, the default file name of the new report in the report menu will be the file name. Override the file name displayed by entering a preferred name in the User Data field of the report properties. Also, to place the new report under the first or second node of the report tree, place a 1 or 2 before the other text in the User Data field. For example, entering '1 My New Report' in the User Data field will place the name 'My New Report' under the first node.

Preview and save the Report Layout

Preview the report by clicking the 'Preview' tab at the bottom of the Report Designer window. Save the report layout by selecting the [Save] button at the top of the Report Designer Window.

Choose a save file location and name your report 'my due cal report'. The file location defaults to the Report folder under the Files Folder (select Open Files Folder from the File drop down menu) for the Ape application to find it.

Report Parameters

Pass Data to SQL Statements at Run Time

Use Report Parameters in Calibration Control to pass information to the report's SQL string at run time. A Parameter dialog will then prompt the user for the required input when generating reports.

The general syntax for a report parameter is:

<%ParameterName|PromptString|DefaultValue|Type%>

Where:

- **Parameter Name:** Must be a unique Parameter name not used for any other parameter or field in the report.
- **Prompt String:** The text displayed in the Parameter dialog asking for input (e.g., Calibration Due By:)
- **Default Value:** The default value of the Parameter.
- **Type:** The Parameter type code of the parameter (see below).

Parameter Type Codes

- **AS:** String
- **AD:** Date
- **AB:** Boolean
- **AI:** Integer
- **ADB:** Double
- **AC:** Combobox (i.e., drop-down lists)

The following is an SQL statement that uses a 'hard coded' date, which means the SQL statement must be edited eachtime the default date needs to change.

```
SELECT *
FROM qryEquipmentMaster
WHERE CalibrationDue < #2018-10-15#
```

Conversely, the following sample shows the same SQL statement using a Parameter, which allows users to enter a datevalue at run time.

```
SELECT *
FROM qryEquipmentMaster
WHERE CalibrationDue < #<%CalibrationDue|Due Date:|AD%>#
```

Date Default Values

When using the AD (Date) Parameter Type code, several options are available for the Default Value:

- **Specific Date:** Exact date with the syntax of #YYYY-MM-DD# for MS Access or 'YYYY-MM-DD' for MS SQL Server.
- **Number of Days:** Positive or negative integer (counting number) indicating the number days added (or subtracted)from the current date.

- **BOM:** Beginning of *current* month.
- **EOM:** End of *current* month.
- **BONM:** Beginning of *next* month.
- **EONM:** End of *next* month.
- **BOPM:** Beginning of *previous* month.
- **EOPM:** End of *previous* month.

The following example uses the EOM Default Value code for the End of the Current Month:

```
SELECT *
FROM qryEquipmentMaster
WHERE CalibrationDue < #<%CalibrationDue|Due Date:|EOM|AD%>#
```

Combobox Parameter Type

Using the AC (Combobox) data type requires the use of the DefaultValue Parameter field to configure the Combobox. Within DefaultValue exactly six sub-fields, each separated by a slash ("/"), are required. This also means that exactly 5 slashes must be present in the DefaultValue field.

The six sub-fields are:

1. **Table or SQL *:** Table or SQL for the records used in the combobox.
2. **Value Member *:** Field name for the value, usually an ID field.
3. **Display Member *:** Field name for the text displayed in the combobox.
4. **Filter:** SQL filter phrase that limits the records displayed (e.g., "Active = -1") when using a Table name instead of a SQL statement
5. **Sort:** SQL sort phrase that places the list of choices in a specific order (e.g., "Name"). When left blank, the DisplayMember is the default sort field. This field is only valid when using a Table name instead of a SQL statement.
6. **Default Value:** One of the values within the list created by the Value Members field (e.g., 55ec4215-7f9b-4e9f-b583-56be9871b895 without apostrophes or quotes).

Although not all fields are REQUIRED (*), their places must be defined with slashes like in the following **DefaultValue** example. Note that there are still five slashes ("/") even when only three sub-fields are used.

```
tblPeople/PersonId/DisplayName///
```


Report Settings

Find and Change Settings of Custom Reports

In Calibration Control within the 'Report Designer' grid (found in the Tools tab of the ribbon menu), find and modify the settings of custom reports. Adjust paper size, orientation, gutters, styles, fonts, and more.


Open Settings

Find most of the settings for the custom reports by double-clicking the Settings node in the top-right corner of the ReportDesigner.




Page Setup

Adjust the Margins and Gutter settings.




Printer Settings

Set paper size, orientation, and other page settings.




Styles

Define new or edit existing font styles.



Global


Modify the Report Designer settings for layout, units, maximum number of preview pages, and word-wrap in the ScriptEditor page.



Update Default Reports


Quickly Update Default Report and Label Templates


In Versions 9.0.4.2 and higher, go to the Utilities tab and click **Refresh Templates**. All of your reports and labels will be removed to the Backup folder and replaced with the most up to date default templates. All of your custom reports/labels will need to be MANUALLY updated. If you need help updating your custom reports, contact us for a screen share and we can show you how to update your reports. If you only added some text and a logo to your custom report, for example, it would be simple to recreate the custom report using the default reports instead of manually updating it.



Manually Update Default Reports Using a 2-Step Process

To find the Reports folder, click on File and select "Open Files Folder", then open the Reports folder to find all of the default reports. When a report is customized, make sure to create a custom name for the file, adding "*_custom*" or "*_yourcompanyname*" works best. During version updates for CC, default reports may be updated and the custom reports will be overridden if the original file name is not changed.



<input type="checkbox"/> Name	Date modified	Type
 Backup	12/6/2018 1:48 PM	File folder
<input type="checkbox"/> AllUserRolePermissions.rpx	11/14/2017 3:15 PM	RPX File
<input type="checkbox"/> AsFound.rpx	6/23/2015 8:25 PM	RPX File
<input type="checkbox"/> AsFoundSql.rpx	6/23/2015 8:25 PM	RPX File
<input type="checkbox"/> CalibrationCertificate.rpx	2/12/2016 9:58 PM	RPX File
<input type="checkbox"/> CalibrationCertificate_Sub.rpx	11/12/2015 4:24 PM	RPX File
<input type="checkbox"/> CalibrationCertificateSimple.rpx	10/26/2015 6:06 PM	RPX File
<input type="checkbox"/> CalibrationCertificateSql.rpx	2/12/2016 10:02 PM	RPX File
<input checked="" type="checkbox"/> CalibrationDue_custom.rpx	11/13/2017 1:13 PM	RPX File
<input type="checkbox"/> CalibrationDueCustodianFilter.rpx	4/19/2018 3:12 PM	RPX File
<input type="checkbox"/> CalibrationDueDateRange.rpx	10/30/2015 3:44 PM	RPX File
<input type="checkbox"/> CalibrationDueDateRangeDeptFilter.rpx	3/6/2018 9:17 AM	RPX File
<input type="checkbox"/> CalibrationDueDateRangeSiteFilter.rpx	11/15/2017 5:59 PM	RPX File
<input type="checkbox"/> CalibrationDueDepartmentGroup.rpx	11/10/2017 5:35 PM	RPX File

Steps to Update Default Reports

1. Delete the Default reports in the Reports folder.
2. Restart Calibration Control.

If the customized reports do not have a new custom name then they will be overridden with the new version update. To accommodate for this, Ape creates a **backup** folder within the Reports Folder and moves all of the old reports into that folder. This means that if the custom reports you had in the main Files Folder are gone, first check the backup folder and rename them before replacing.

SQL SELECT Statement

Useful with Calibration Management Software Reports

SQL stands for Structured Query Language. SQL is the basic language of most common databases, including MS Access and MS SQL Server, the two databases that Calibration Control can use. For the purpose of report writing, the part of the SQL language we need to focus on is the SELECT statement, which retrieves information from Calibration Control.

Although SELECT statements are used by reports, web pages, on screen displays, and even when moving data between applications, we will need SQL primarily for reports.

A Primer on Tables

Before we get started, let's make sure we understand the source of our data when we use a SELECT statement. Within a database, information is stored in tables that look something like a spreadsheet with columns and rows. Unlike a spreadsheet, a table uses records (horizontal rows) and fields (vertical columns).

Imagine each record in the database as a photocopied standard form used for keeping track of test equipment and imaginethat these forms are kept in a file. The file may have a hundred forms, each with information describing a specific piece of test equipment. Just like a table, these forms can be sorted in different orders and records can be scanned by a single fieldon each form; it just takes a bit longer with paper compared to a database table.

The Equipment Master View

Within Calibration Control, you will probably derive most of your reports from the qryEquipment view (also called a query) because it has most of the required fields for your test equipment records. The fields in the view will also be easier to readthan their corresponding codes in the root table. Think of a view as a way to pre-package part of the SQL complexity that makes data easier to work with.

As an example, the contents of the tblEquipmentMaster table (where the data is actually stored) contain fields like ModelNumberID with meaningless numbers in the fields (see image below). On the other hand, look at the qryEquipment and you will see the actual model numbers and descriptions that look familiar. This is because the view has SQL code in itthat looks up and displays the meaning of the ModelID code so you don't need to.

The image shows two screenshots of a database interface. The top screenshot displays the 'tblEquipmentMaster' table with columns: EquipmentID, ModelDesc, SerialNumber, FrequencyUni, and Ca. The bottom screenshot displays the 'qryEquipment' view with columns: MasterGuid, EquipmentID, Alternateld, SerialNumber, and PersonalPro.

EquipmentID	ModelDesc	SerialNumber	FrequencyUni	Ca
SAMPLE-184	OD Micrometer	2344	0	
SAMPLE-301	Acme 1" Micrometer	1015	0	
SAMPLE-220	Caliper	135	0	
SAMPLE-081	Loctite	321	2	
SAMPLE-180	Indicator	2209	0	
SAMPLE-131	Drop Indicator	2096	0	
SAMPLE-275	Cylindrical Ring Gage	2003	0	
SAMPLE-101	Depth Micrometer	2055	0	
SAMPLE-007	Indicator	2247	0	

MasterGuid	EquipmentID	Alternateld	SerialNumber	PersonalPro
4fa23f81-b95a-	SAMPLE-057		2152	(
88bf5641-bf43-	SAMPLE-058		2150	(
07e827a0-f9de-	SAMPLE-017		186	(
a26c469b-3e4c-	SAMPLE-059		2148	(
74b113c4-59a9-	SAMPLE-060		2153	(
3b446cf0-f83f-	SAMPLE-061		2156	(
f979d4c4-a740-	SAMPLE-062		2092	(
65ede61e-4eac-	SAMPLE-063		2069	(

The SQL SELECT Statement

There are four main parts of the SELECT statement that we need to cover:

- ♦ SELECT – (IDs the fields from a table that will be included)
- ♦ FROM – (IDs the table where the fields come from)
- ♦ WHERE – (Defines the filter that includes only the records you want to include)
- ♦ ORDER BY – (Defines the sort order of the records)

While the SELECT and FROM clauses are always required, only the WHERE and ORDER BY clauses need to be included when filtering and sorting. Otherwise, the results will include all records and in no particular order other than the physical order of the underlying table.

Here is an example of a common SQL statement used to return all the fields from the qryEquipment view:

```
SELECT *
FROM qryEquipment;
```

Pretty easy, right? Notice the asterisk (*)? The asterisk is a wildcard that includes all fields. With the relatively easy SELECT statements that we need in our Calibration Management Software, using the wildcard is usually the best because you don't need to worry about forgetting a field when designing the report.

Now let's assume that we want to filter and sort our results in the following example:

```
SELECT *
FROM qryEquipment
WHERE DepartmentCode = 'QA'
ORDER BY Location;
```

You can even add multiple filters and multiple sorts with the following **MS Access**:

```
SELECT *
FROM qryEquipment
WHERE (DepartmentCode = 'QA') AND (CalibrationDue <> #5/1/2019#)
ORDER BY Location DESC, SerialNumber;
```

Note that the date value has number signs (#) around it rather than the single quotes of the text values? You need to use the # sign when using MS Access and the single quote (') when using SQL Server. Although both databases use SQL Server, there are still slight differences.

Here's the same code for **MS SQL Server**:

```
SELECT *
FROM qryEquipment
WHERE (DepartmentCode = 'QA') AND (CalibrationDue <> '5/1/2019')
ORDER BY Location DESC, SerialNumber;
```

The final example deals with filtering for numeric values and is compliant with both **MS Access** and **MS SQL Server**. Note that the value (1) uses neither the single quote (') or the number sign (#)?

```


SELECT *
FROM qryEquipment
WHERE FrequencyUnits = 1;

```

You can learn more about the SQL SELECT statement through an Internet search or go directly to Microsoft for SQL Server SELECT Statements or MS Access SELECT Statements .

Dashboard Charts

Under the Calibration tab in the Visual group, you'll find a new Dashboard including Column charts, Pie charts and more. You can change the look and feel of the dashboard by customizing each chart.



Refresh the Dashboard and the randomized chart colors by right-clicking a grey area and clicking Refresh. Additionally, you have Options to modify the alignment and spacing of the charts.

Edit Dashboard Charts


To edit each chart, double click the chart or right-click and select Edit.

Chart tab

If you want to remove a chart from the Dashboard, deselect the Active checkbox. You can also edit all Dashboard Charts in the Tools tab of Calibration Control.


Sort Order determines the order of the charts in the Dashboard.

Height and Width: The default height of 0.5 refers to 50% and the Width of 0.2 is 20%. Any number above 0.9 is in pixels. Click the **Refresh Chart** button to see the reflected changes while editing.



SQL | Grid tab


Whether you are using a MS Access or SQL Server database, the SQL for the appropriate database type will be automatically selected. The SQL script produces the chart data and you can view the Results in grid format by clicking **Refresh Grid**.



Legend tab


First, make the Legend visible by checking the Visible checkbox and Refreshing the Chart. Then, you can use the

available settings to modify the legend.




Pie | Doughnut tab

If you are editing a Pie or Doughnut chart, use these available settings to modify the chart.



Column | Bar | Line tab

Column or Bar spacing settings are located on the left. Settings for Line or Spline charts are located on the right.



Edit Titles

By default, all charts have a Top title already visible. You can edit each available title in this dialog but only the titles marked as Visible will appear on the chart.

Equipment by Cal Frequency

Property	Top	Bottom	Left	Right
Text	Equipment by Cal Frequency			
Clip Text	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Best Fit	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Font Size	11.00	8.00	8.00	8.00
Font Style	B / U S	B / U S	B / U S	B / U S
Font Color	Black	Black	Black	Black
Extent	25	25	25	0
Horizontal Align	Center	Center	Near	Near
Vertical Align	Center	Center	Near	Near
Margins	5 5 5 5	5 5 5 5	5 5 5 5	5 5 5 5
Orientation	Horizontal	Vertical Left Facing	Vertical Left Facing	Vertical Left Facing
Orientation Angle	0	0	0	0
Visible	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Edit XY Titles

You can edit XY chart titles here. The Visible titles depend on the type of chart you are editing. For example, Stack Columncharts use X and Y titles.

The screenshot shows a dialog box titled "Edit Chart Labels Upcoming Calibrations by Month (Load Time: 0.1 Seconds)". It contains four columns of settings for different chart elements: "x", "x Series", "y", and "y Series". Each column has a "Visible" checkbox at the bottom. The "x" and "y" columns have "Visible" checkboxes that are unchecked, while the "x Series" and "y Series" columns have "Visible" checkboxes that are checked. The settings include Extent, Range Type, Range Min, Range Max, Font Size, Font Style (Bold, Italic, Underline, Strikethrough), Font Color, Format, Horizontal Align, Vertical Align, Orientation, Orientation Angle, and Padding. The "Format" field for "x" is "<ITEM_LABEL>", for "x Series" it is empty, for "y" it is "<DATA_VALUE:0>", and for "y Series" it is empty. The "Horizontal Align" for "x" and "y" is "Center", for "x Series" is "Center", and for "y Series" is "Near". The "Vertical Align" for "x" and "y" is "Center", for "x Series" is "Center", and for "y Series" is "Near". The "Orientation" for "x" and "y" is "Horizontal", for "x Series" is "Custom", and for "y Series" is "Vertical Left Facing". The "Orientation Angle" for "x" and "y" is "0", for "x Series" is "320", and for "y Series" is "0". The "Padding" for "x" and "y" is "2", for "x Series" is "2", and for "y Series" is "5". There are "Cancel" and "Submit" buttons at the bottom.

Displaying Charts

Use Pie, Column, and Bar Charts to Visualize Your Data

Display the Charts grid by clicking on the Charts icon in the Calibration tab of the ribbon menu.






Chart Data Tables

Use the Data Table combo box to select what information is shown in the Chart.




Customize a Chart by editing the default Title and by changing the settings for Chart Type, Color Model, and Legend position. Use the Data Table combo box to choose which department information is used.




Saving and Printing Charts

Right-click the Chart screen to Save (several formats) or Print the Chart. 'Save as PNG – Transparent Background' is recommended for use in reports and presentations with backgrounds other than white.




3D Charts

Change the viewing angle on 3D Chart Types by holding down [Alt] and then clicking and dragging the mouse to the bestangle. Zoom in and out of the Chart by rolling the mouse roller back and forth while keeping the [Alt] button pressed.



'Min Others Pie Slice' (for pie charts only)

Use 'Min Others Pie Slice' to group very small percentages into an 'Other' category. This helps to declutter the pie chart when there are too many slices.



Due Cal Calendars

Visualize Due Cal Equipment in Calendars

Display a Calendar of Equipment Calibration Due dates by clicking the Calendar icon in the Calibration tab of the ribbon menu.



Toggle between Month, Week, or Day views using the tabs at the top of the Calendar screen.

Click on specific dates in the Calendar to view which tools are due to be calibrated on that date. The 'Equipment Due Cal' label automatically displays how many tools are due for calibration on that date. Click in the default label to customize.

Switch between viewing month, week, or day by using the tabs shown in the top left corner of the dialog below.

Month	Week	Day
Sunday		Monday
Nov 25	26	27
8 Equipment Due Cal		
2	3	4

Quickly view Equipment Due for Calibration dates in different months by clicking through the mini calendar view in the upper left-hand corner of the screen. The current day is highlighted in a red box while the day selected is shown in orange. The **bold** dates are ones with at least 1 equipment that will be due for calibration.

Adjust the number of weeks visible in the Calendar screen by clicking up or down in the 'Visible Weeks' box (Month view only).

◀ December 2018 ▶

S	M	T	W	T	F	S
25	26	27	28	29	30	1
2	3	4	5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31					


▶ January 2019 ◀

S	M	T	W	T	F	S
	1	2	3	4	5	
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30	31	1	2
3	4	5	6	7	8	9


Visible Weeks: ▲▼

Printing Calendars

Print a 5-Week, 1-Week, or 1-Day Calendar by viewing dates you want to be printed in the Calendar screen and clicking [Print].



Customize printing Paper Size and Source, Orientation, and Margins by clicking on [File], then [Page Setup...].




Security Methods - How They Work

Database Security

When using a MS SQL Server database in Calibration Control, the built-in SQL Server and Windows authentication will keep the database protected. Otherwise, when using MS Access the database file (apecal.mdb) can use a database password.

Activating User-Authentication (Sign-In Mode)

Any user can activate User-Authentication mode after ensuring *at least one active User* has Admin privileges. Only a User with Admin privileges can deactivate User-Authentication.



What an Administrator Can Do

Administrators have full access to all parts of the application that require any level of security. An Admin can create Users and change the privileges of any other User, including other Admins. Administrators cannot remove themselves from being an Admin while User-Authentication mode is activated; one Admin must remove the Admin privileges of another. This is a safety feature to ensure there is *at least one Admin* while the application is in User-Authentication mode.

Administrators and Passwords

When an Administrator creates or changes the password for any User other than themselves, that User will be required to change their password the next time they sign in. Use the Password Security Dialog page for help in implementing and adjusting Password Security methods.

Pre-Defined Roles

There are six pre-defined security roles, which are: Administrator, Super User, Supervisor, Technician, Layouts and No Role (read-only user). Although these Roles cannot be changed, Admins and Super Users can change which Role has access to which Permission. In addition to a role, a user can be a Terminal User. This means they can **only** use Calibration Control's Terminal Mode.

Hierarchy of Roles

A User with no Roles assigned can see, print, and export any information within the database. Each of the other Roles has the privileges of a User with no Role in addition to:

- **Administrator:** Can do anything that requires security.


- **Super User:** Can do anything an Administrator can do *except* add or edit Users.
 - **Supervisor:** Same authority as Technician and Production, can also edit Technician Name and Status of calibration records, and can remove relationships between Equipment and Jobs.
 - **Technician:** Can create and edit equipment records. Cannot edit the Technician Name (added automatically) or the Status of calibration records.
 - **Layout:** Has the authority to make most application layout changes relating to the look and feel of the database.

Permissions

There are pre-defined Permissions for specific areas of the application that require a given level of security. Although each Permission has a default minimum Role, Administrators can change the Role of any of the Permissions. Non-admins can only view the Permissions grid but cannot make changes.

Change Log

All field and label changes are tracked in the 'Change Log' grid found in the Data Grids tab of the ribbon menu. It records time and date, user (if User-Authentication mode is activated), the screen where the changes were made, machine (computer) name, and the detail of the change. The change detail includes field names and before/after data.



Users

Creating a User

Select the **Users** icon under the **Utilities** tab in the ribbon menu and create a new User record by double-clicking in the gray area or right-click for the context menu and choose [New Record].

User Tab

The only required fields are 'User Name' and 'Password'. When creating a password, select from the various checkboxes below the password field. By default, the password will be marked as Temporary and the new user will be marked as Active.

If a Person record already exists for this User, select that record in the Person dropdown.

You can choose a language for the user in the Language dropdown. The language is unique to each user.

Roles

In the Roles section, select a role for this User. To view the Permissions each role has by default, go to the Permissions grid located in the Utilities tab. Users with **No Role** (read-only rights) are free of charge.

In addition to the user's role, selecting **Terminal User** will allow the user to connect to Calibration Control's Terminal Mode but the user will NOT be allowed to use the full version of Calibration Control.

The screenshot shows a dialog box titled "Edit User admin (Load Time: 0.0 Seconds)". It has four tabs: "User", "Advanced", "Notes", and "Meta". The "Advanced" tab is active. The "User" section contains: "User Name:" with the text "admin"; "Password:" with a masked field of six dots; "Show Password (new passwords only)" with an unchecked checkbox; "Ignore Password Rules" with an unchecked checkbox; "Temp Password (change during next login)" with an unchecked checkbox; and "Is Active" with a checked checkbox. The "Person:" dropdown menu shows "David Wilcox". The "Language:" dropdown menu shows "English (US)". The "Comments:" field is a large empty text area. The "Roles" section contains: "Administrator (can do anything)" with a selected radio button; "Super User (edit anything but other users)" with an unselected radio button; "Supervisor (can approve calibrations)" with an unselected radio button; "Technician (can create calibration records)" with an unselected radio button; "Layouts (can change layouts)" with an unselected radio button; and "No Role (read-only user)" with an unselected radio button. There is also an unchecked checkbox for "Terminal User (can only use terminal)". At the bottom, there are "Cancel" and "Submit" buttons.

Advanced tab

Lockouts

Here you can Remove Lockouts, view the number of failed log in attempts and the last locked out date.

Logins & Password

Calibration Control keeps track of the user's last log in date and the number of log ins since the user record was created.

Domain

A Windows User name can be used instead of a username and password.

Advanced Skills

In this section, you can give users the ability to edit SQL statements or HTML in the software (i.e. Auto Emails & DashboardCharts). By default, these checkboxes are not selected.

The screenshot shows the 'Edit User admin' dialog box with the 'Advanced' tab selected. The dialog is divided into several sections:

- User:** Contains a dropdown menu for 'Advanced'.
- Lockouts:** Includes 'Failed Login Attempts' (input field with '0'), 'When Locked Out' (dropdown), a 'Locked Out' checkbox, and a 'Remove Lockout' button.
- Logins:** Includes 'Last Logged In' (dropdown with '05/10/2020 01:40 PM'), 'Login Count' (input field with '10'), and 'Password Last Changed' (dropdown).
- Domain:** Includes 'Windows User' and 'User Domain' (input fields) and a 'Use Current Machine and Account' button.
- Advanced Skills:** Includes two checked checkboxes: 'Can Edit SQL Statements' and 'Can Edit HTML'.
- Buttons:** 'Cancel' and 'Submit' buttons are located at the bottom.

Notes Tab

Add or Link Note records to each User record.

Sign-In Mode (Enabling Authentication)

Set up and Turn on User Authentication


Utilize User Authentication to keep track of changes made in Calibration Control. Activating Sign-In mode is simple. All you need is at least one user with an Administrator role.

Creating an Admin User

Create a User record and select Administrator as the role.

Turning On User Authentication

After creating at least one user with an Administrator role, open the Options dialog > Advanced tab > Security and select **Activate Sign-In Mode**.



Click Yes and Calibration Control will restart. Then, the Sign-In dialog will open.


Note: If you select Remember Sign-In, you can undo this in Calibration Control by going to File > Sign Out and restarting the program.



Password Security

Using the Password Security Dialog

Go to the 'Options' dialog and go to the Advanced tab, then Security and select **Password Security**.



Use the Password Security dialog to set expiration times, password strength, and naming rules.

A default value of zero (0) means that a setting is not active.

Maximum Settings

Days Before Expiration:

Unused Days Before Expiration:

Invalid Password Attempts:

Minimum Characters

Upper Case:

Lower Case:

Special (e.g., ~@#\$\$%^):

Length:

Check Settings

Different From Temp

Cannot Contain User Names

Cancel Submit

Maximum Settings

Set the maximum limits for password usage:

- **Days Before Expiration:** Number of days before a User's Password expires and must be changed.
- **Unused Days Before Expiration:** Number of days a User can go without signing in before their Password expires and must be changed.
- **Invalid Password Attempts:** Number of times an incorrect Password can be entered for a given User before that User is locked out. Entering a correct Password for the same User sets the number of incorrect password attempts for that user back to zero (0).

Minimum Characters

Minimum Character password rules are enforced when a User creates a *new* password after signing in with a *temporary* password. The password assigned by an Administrator in a User dialog is temporary and is not subject to these rules.

- **Upper Case:** Minimum number of Upper Case characters required in a Password.
- **Lower Case:** Minimum number of Lower Case characters required in a Password.
- **Special (e.g., ~@#\$\$%^):** Minimum number of Special Characters required in a Password.
- **Length:** Minimum character length required for a Password.


Check Settings

Additional safety settings for passwords.

- **Different From Temp:** When checked, Passwords cannot be the same as a Temporary Password set by an Administrator.
- **Cannot Contain User Names:** When checked, Passwords cannot contain the First, Last, or User name of registered User.

Reset Login by Admin

When a User is locked out due to invalid password attempts, an Administrator can reset the User Login information from the Advanced tab of the 'Edit User' dialog.



The screenshot shows the 'Edit User' dialog box with the 'Advanced' tab selected. The 'Lockouts' section is highlighted with a red box. It contains the following fields and controls:

- Failed Login Attempts: 0
- When Locked Out: 10/27/2015 12:00 AM (dropdown)
- Locked Out
- Remove Lockout button

Other sections visible include:

- Logins: Last Logged In: 07/28/2020 08:51 AM (dropdown), Login Count: 8
- Password: Password Last Changed: 07/08/2020 12:00 AM (dropdown)
- Domain: Windows User, User Domain (text boxes)
- Advanced Skills: Can Edit SQL Statements, Can Edit HTML
- Buttons: Cancel, Submit

Database Password Encryption


Encrypting an MS Access Database with a Password

Part of the data security within Calibration Control includes the password encryption of MS Access databases. This meets the standards of the FDA's 21 CFR Part 11 requirements.


You can take the following steps to add a password to your MS Access database:

In order to encrypt your database with a password, you must first open the database in **Exclusive Mode**.

- Start Microsoft Access.
- Click Open, then Browse.
- Browse to your database file location and select the database by clicking it once.
- Click the downwards arrow next to Open, then click **Open Exclusive**.



- In MS Access, click **File** then **Info**.
- Click **Encrypt with Password**. Your version of MS Access may say "Set Database Password"
- Type in your password, type it again to verify and click OK.
Note: Use a strong password and store the password in a secure place. If you forget your password, Microsoft won't be able to retrieve it.
- Start Calibration Control and you'll be prompted to enter your database password.



If you do not have MS Access installed on your computer, or you could not complete the above steps successfully, contact us and we can help.


Permissions

Change the Roles Required for Any Permission


Permissions are now in the **Utilities tab** (Security group). Any user can view the Permissions grid but only Administrators (by default) can edit Permissions.

****Pre Version 9.2**

First, go to File > Options. For Version 9.1.2, go to the Advanced tab > Security > Edit Permissions. For Versions 9.1.1 and older, go to the Security tab > Edit Permissions.



Permissions can be edited by double clicking each record. You can change each permission's required Minimum Role by selecting from the dropdown. In the reference below, anyone with a Technician role or higher has Permission to Add Attachments.



Reset Options


There are two reset options, one for *each* Permission and one for *all* Permissions. You can reset EACH permission within each Permission dialog and you can reset ALL permissions by right-clicking on the Permissions grid and selecting Reset AllPermissions.

Change Log

Audit Changes Made in Calibration Control

The Change Log grid shows a complete list of all the changes made within the program, including time, date, user, and technical information. Find this feature in the Data Grids tab of the ribbon menu.

Determine how many records are shown using the combo box on the left-hand side and search for a specific record using the search bar on the right-hand side. Filter options are the same as other grids, and can be found in the row under the column headers; enable or disable the Filter row feature by right-clicking in the grid to bring up the context menu and toggling the [Filter Row] option. Each column has a default search option of "Contains", and hover the cursor over the row to reveal a drop-down box and the clear filter button.



Columns of the Change Log Grid

- **Change Time:** Date and time that the change took place.
- **User:** When sign-in mode is enabled, this will show the name of the user who made the change.
- **Dialog:** The dialog box in which the change occurred.
- **Change Detail:** An auto-generated description of the change that was made in the database.
- **Person:** The full name of the User profile who made the change.
- **Computer:** Which computer the change was made on.

FDA 21 CFR Part 11 Compliance

FDA Software Compliance

Ape Software's calibration management software (Calibration Control) complies with FDA 21 CFR Part 11 Rule on electronic records and signatures by containing the following functionality:

Authentication at User Level

Calibration Control has the ability to individually authenticate each user for the purpose of audit trails and electronicsignatures.

File Integrity

The default MS Access back-end database is secured through password encryption, which includes tamper detection. If an MS SQL Server is used as the back-end database, instead of MS Access, the necessary security features are already part of SQL Server.

Change Management (Audit Trails)

All user activities (e.g., record changes & calibrations) are automatically recorded in the Change Log to provide an audittrail.

Electronic Signatures

Re-entry of password required upon electronically signing (approving) Calibration Event (optional). Also, inherent in the Change Log, all user activity relating to equipment maintenance and calibration is recorded (electronically signed).

Password Strength

Ability to define the strength required for passwords including numbers of upper, lower, number, and special characters. Ability to define expiration intervals and other features.


Inactivity Settings

Ability to define the number of inactivity minutes in a usage session before the application times out and users are requiredto sign in again.

Sessions

See Who is Using Calibration Control

See which users or machines (i.e., computers) are currently using or have used the software. Time and date stamps show when the software was used (start and end times) and also shows if the session was interrupted due to network problems.Find the Sessions grid in the Utilities tab of the ribbon menu.




Sessions Grid & Dialog

The Sessions grid shows a READ-ONLY log of sessions and displays the following fields:


- **Machine:** Shows the name of the machine (computer) the session was conducted from.
- **User:** Displays the name of the User (e.g., Admin) of the Session. A blank 'User Name' field indicates that theSession was conducted while User Authentication was disabled (i.e., Machine Authentication).

- **Person:** Displays the name of the person from the Person table record linked to the User record at the time of theSession. If this field is blank while the User Name field is not, it means that the User record was not linked to a Person record.
- **Session Start:** When the user signed in to the software.
- **Session End:** When the user signed out of the software. If this field is blank the session is still active or was interrupted (e.g., network interrupted). When the user for an interrupted session signs in again, the previous sessionwith a blank 'Session End' field will automatically close (receive a date and time) and the 'Interrupted' checkbox will be ticked to indicate an ungraceful exit.
- **Interrupted:** This box is checked when the session was ended ungracefully (e.g., software not shut down normallyor network interrupted).
- **Session Bumped checkbox:** Check this box to terminate a session for the selected user. This is helpful when theuser count exceeds the software license agreement and is preventing registered users from using the software.



View All or Current Sessions


Toggle between viewing All or Current sessions by right-clicking in the grid for the context menu and selecting 'All Sessions' or 'Current Sessions'. Sort the sessions grid using the filter row to easily find the record you are looking for.



Machines


Control Computer Access in the Machines Grid

The **Machines grid** shows a list of all the machines (computers) that connect to the Calibration Control database and the number of times they have connected. Find the Machines grid icon under the Utilities tab of the ribbon menu.




The **Machines dialog** contains a name field which is auto-populated by the software when a machine (i.e., computer) user is signed in, a space for notes to be added, and an 'Is Blocked' checkbox. The only fields that can change in the Edit dialog are the 'Is Blocked' checkbox and the Notes field. Only check the 'Is Blocked' checkbox if you need to remove a machine from the license count by blocking that machine from using the software.

The **Meta Tab** displays the meta-data pertaining to the selected record and IS A READ-ONLY INFORMATION tab showing the date and time the record was created, the user that created it (if they are signed in), as well as the most recent date and time the record was updated, and the machine that made the change (if no user is signed in).



In addition to the IsBlocked checkbox, you can quickly block ALL machines besides the one you are current using. Right-click the grid and at the bottom of the context menu, select **Block All but Me**.



Licensing and Authentication

Ape Software applications can be used in "Machine" or "User Authentication" modes. The default mode is Machine (no sign-in) although User Authentication (sign-in) can be enabled for additional security controls and more accurate changerecords (e.g., who changed what and when).

Users can only have one active session at a time.

Machine Authentication (default)

Machine Authentication means that User Authentication (sign-in mode) is **not** enabled. Therefore, the license count is managed by the names of the individual computers (i.e., machines) that connect to the database. This is the default method of authentication.

User (Machine Perspective)

When using Machine Authentication all users who connect to the database have full access, so everyone (each machine) is considered a user for the purposes of license use.

Per Seat Licensing

With a Per Seat license and User Authentication **not** enabled, the number of machines that connect to the database cannot exceed the license count.

For example, with a 2-user license no more than two unique machines can connect to the database.

Concurrent (at the same time) Licensing

With a Concurrent license, and User Authentication **not** enabled, the number of machines that connect to the database *at the same time* cannot exceed the license count.

For example, with a 2-user concurrent license, unlimited machines can connect to the database, but no more than two of them *at one time*.

User Authentication

User Authentication means that at least one User was created, at least one of those Users was given Administrator-level access, and User Authentication **is** enabled. In this scenario, license count is managed by the number of Users created within the software who have some level of write access within the software (e.g., Technician, Supervisor, or Administrator).


User (Authenticated)

When User Authentication is enabled only users that have any level of write access (i.e., have a role checked in their Userrecord) count as a user for the purposes of license use.

Per Seat Licensing

With a Per Seat license, and User Authentication **is** enabled, the number of authenticated users that connect to the database cannot exceed the license count.

For example, with a 2-user license there should be no more than two authenticated users (listed in the "Users" grid) that have a role defined.



Concurrent (at the same time) Licensing

With a Concurrent license and when User Authentication is enabled, the number of authenticated users who connect to the database *at the same time* cannot exceed the license count.

For example, with a 2-user concurrent license any number of authenticated users can exist in the User Grid but no more than two of those users can connect *at the same time*.

Simple Troubleshooting

Fix Most Setup or Configuration Problems

Most configuration or setup problems encountered while using Calibration Control can be resolved or isolated using the following steps.


1. **Version:** Confirm that the most current version is installed: Download Ape Software.
2. **Reset:** Try a simple reset, by closing Calibration Control and either holding CTRL and starting Calibration Control or manually deleting the configuration file in the Settings folder: 'general.config' file. The Settings Folder is located at:

C:\Users\Public\Documents\Ape Software\Calibration Control

Note: The database location is defined in the general.config file. Therefore, if the general.config file is deleted, but an 'apecal.mdb' file is still in the Settings Folder, which is the default location for the Data Folder, then the application will automatically connect to that database (i.e., apecal.mdb). If the apecal.mdb file does not exist in the Settings folder, the [Choose Database Format] dialog will appear. Follow the prompts of the dialog to re-connect to the database.

Name	Date modified	Type
Attachments	6/18/2019 3:44 PM	File fold
Emails	6/6/2019 12:27 PM	File fold
Labels	6/6/2019 12:27 PM	File fold
Language	6/21/2019 1:55 PM	File fold
Layouts	6/20/2019 1:28 PM	File fold
Reports	6/6/2019 12:27 PM	File fold
apecal.css	9/5/2011 3:35 PM	Cascadi
apecal_sample.lac...	6/21/2019 1:15 PM	Microsc
apecal_sample.mdb	6/21/2019 1:28 PM	Microsc
connectionlog_v8_...	6/20/2019 3:18 PM	Microsc
connectionlog_v9_...	6/21/2019 1:15 PM	Microsc
Engineering.html	6/20/2019 2:02 PM	Chrome
errorlog_v8_11_9_0...	6/19/2019 3:15 PM	Chrome
errorlog_v8_11_9_0...	6/19/2019 3:15 PM	Chrome
errorlog_v8_11_9_0...	6/19/2019 3:15 PM	Chrome
errorlog_v9_0_0_0....	6/20/2019 1:44 PM	Chrome
general.config	6/21/2019 1:41 PM	CONFIG
index.html	6/20/2019 2:02 PM	Chrome

4. **Isolate Access Database:** Assuming none of the above measures resolve the problem, this step will help identify if there is a problem with the database.
- Open the Data folder from the File menu and then close Calibration Control.
 - Rename the 'apecal.mdb' file to something else like '_apecal.mdb'. This will make CC recreate a new default database in the Data folder.
 - Restart CC and when the dialog below appears, select the [MS Access] button to create a new sample database.
 - If the problem is resolved in the sample database, then zip your '_apecal.mdb' file (your database) and the 'errorlog.html' and email these files to service@apesoftware.com with a description of what you are experiencing.
 - If the problem continues to occur with the sample database, contact Ape Software for assistance.



Problem Steps Recorder

Utilize the built-in Problem Steps Recorder to record the steps you are taking to reproduce the error. Contact Ape Software at service@apesoftware.com with the following information, or as much of it as possible:

- a. Description of what you are doing that leads to the issue you are trying to communicate. Include as many details as you can and include relevant screenshots.
- b. Problem Steps Record File (see above), if possible (file has an *.mht extension)
- c. The 'errorlog.html' file from the Settings folder
- d. The 'general.config' file from the Settings folder
- e. The 'apecal.mdb' file from the Data folder
- f. Zip the above files into a single file/folder if you know how

Crash Reports

The Ape Software Crash Reporter dialog will appear if Calibration Control encounters an error. If this dialog appears, you can easily help us find the cause of this error and we can help you by preventing it from happening again.

Tell Ape About This Crash


If you receive a crash while using Calibration Control, please check the Tell Ape About This Crash checkbox. This sends us the Error Detail and other relevant information. If the checkbox is not checked, we are not receiving any details of any crashes you may encounter. However, if you contact us about a crash, we'll need to locate the Log of All Errors (error log) stored on your computer.

Description

While the Error Detail is one of the most important pieces of the puzzle, you can add your own description of what happened and it can help us put all the pieces together. Try to be as descriptive as possible so we can follow your steps while testing.

Follow up

Finally, your email address is not required to send us details of the crash but it could save both of us some time. If the error is crucial and keeping you from working, your email address could allow us to easily contact you. Even if the error appears to be minor, we may need to contact you if you did not leave a description or if we need more clarification.




Find Lost Records in Equipment Grid

Troubleshooting Missing Records


The two most common things that make it look like some or all of your records are missing in Calibration Control (our **Calibration Management Software**) are (a) un-cleared filters in the filter row or (b) records with hidden status codes.

Clearing the Filter Row

The row at the top of the grid that looks blank is actually a filter row. Learn more about the filter row by checking out the filter row help topic. A *single column* can be cleared by clicking on the clear filter button (funnel with a line through it) within that column. Clear *all of the filters* in the row by clicking on the far left clear filter button.




NOTE: If the phrase 'FILTERED Records' is displayed in the grid header, then not all of the fields have been cleared.



Equipment					
[Equipment] - 126 FILTERED Records					
Equipment ID	Model	Description	Serial	Manu	
SAMPLE-017	CMM-421	CMM	186		
SAMPI F-019	ODM-d80b85	OD Micrometer	386		

Showing Hidden Records

Some of the records may be hidden because they are assigned a status code with a 'Hide these records in the Browse Table' attribute set to true. Learn more about hiding and showing records with certain Status Codes by reading the Hiding Equipment Records by Status Code help topic. Just right-click in the Browse Grid and toggle the Show Hidden Records option.



Lost or Invalid Product Key / Lost All Records


Troubleshoot an Invalid Product Key or Lost Records

In Calibration Control, if it looks like some records are lost or that a product key is absent or no longer valid, the most common scenario is that there are multiple databases on your computer and you are currently connected to a new blank database. Therefore, the solution is to redirect Ape software back to the master database.

How It Happens

Network Database: The database has been moved to a common network location for multiple users, safety, security, and/or ease of backup.

Database Not Available: When launching Calibration Control, the software attempts to connect to the network database, but the network is down, paths are remapped, or permissions are changed. When the following dialog is prompted, you might quickly hit the [Enter] key to get past it, thereby unknowingly creating a new blank database in the default location on your computer's (C:) drive.



How to Fix It

Follow the instructions under the 'Add Additional Users' section of the Creating a Multi-User Environment help topic.

Related Help Topics

- Locating a Database File (MS Access)
- Creating a Multi-User Environment

Troubleshoot Calibration Control Not Working

Why Did Calibration Control Stop Working?

There are a few reasons why Calibration Control would stop working. Use these troubleshooting tips to resolve the problem. If these steps do not work, give us a call to set up a screenshare.

Missing Database Drivers

The most common problem that would cause CC to stop working is that the database drivers are not up to date or were deleted by another program. Sometimes the installation of other versions of MS Office will remove these drivers and they may need to be re-installed. Download the most current MS drivers here, and see if that resolves the issue.

Missing .NET Drivers

If installing the database drivers does not help resolve the issue, try installing the .NET drivers.

Missing Read/Write Privileges

Calibration Control requires users read/write access in a few folder locations that are listed in Minimum Requirements. Without the proper permissions to those folders, users cannot save changes in Calibration Control. The error messages usually indicate if the user is denied the privilege to create, edit, delete in the network folder location where the database resides, (e.g., "Insert Unsuccessful: The database insert failed." or "Data provider failed while executing a provider command. Operation must use an updateable query.").

Cannot Connect to an Unstable Database

If users receive an error message opening Calibration Control that reads, "Data Provider could not be initiated," this points to problems within the database and Calibration Control cannot connect to it. Most of the time, a corrupt MS Access database can be quickly repaired. Refer to our troubleshooting tips for MS Access Database Corruption and How to Repair your MS Access Database. If SQL Server database is corrupt, contact us for troubleshooting help.

Server is Waiting for a Response

The server might require a response relating to the database locking file (*.ldb). See if there are any messages on your server that need to be dismissed.

Get Ape Support

If installing the drivers does not work, send us an email and attach the following files from this folder (C:\Users\Public\Documents\Ape Software\Calibration Control) on your computer:


- connectionlog*.csv
- errorlog*.html
- general.config

If you can, zip up a copy of your database and add it to these files.

Troubleshoot Database Location

Why Am I Seeing the "Choose Database Type" Dialog?

When Calibration Control starts up, it looks in the Settings Folder to find the "general.config" file which tells CC where the database is located. If that file is missing, has been moved, or has been renamed CC prompts the "Choose Database Format" dialog. To navigate to the misplaced database, choose the "BROWSE to Existing Database" option and click [Continue].



In order to browse to the existing database you must first know where the database has been moved to. This is a question that Ape Software does not know, but once the existing database has been located by you then the Calibration Control software will remember the current location.

Last Known Location of Database

When the location of the database is unknown, check the Connection Log which is located in the Settings Folder [C:\Users\Public\Documents\Ape Software\Calibration Control].

Name	Date modified	Type
Attachments	9/10/2019 4:38 PM	File folder
Backup	9/16/2019 12:09 PM	File folder
Emails	9/16/2019 12:09 PM	File folder
Labels	9/16/2019 12:09 PM	File folder
Language	9/24/2019 12:24 PM	File folder
Layouts	9/24/2019 11:30 AM	File folder
Reports	9/16/2019 12:09 PM	File folder
apecal.css	9/5/2011 3:35 PM	Cascading Style S...
apecal.mdb	7/25/2019 12:25 PM	Microsoft Access ...
apecal_sample.laccdb	9/24/2019 12:30 PM	Microsoft Access ...
apecal_sample.mdb	9/24/2019 12:16 PM	Microsoft Access ...
connectionlog_v9_0_0.csv	9/9/2019 3:07 PM	Microsoft Excel C...
Engineering.html	9/24/2019 12:22 PM	Chrome HTML Do...
errorlog_v9_0_0_0.html	8/1/2019 3:22 PM	Chrome HTML Do...
general.config	9/24/2019 12:43 PM	CONFIG File
index.html	9/24/2019 11:55 PM	Chrome HTML Do...

Within the Connection Log Excel document, locate the most recent date and its database path will show a network location of where the database used to be. This may give some insight into where the new database location would be.


To reset Calibration Control to factory default, delete the "general.config" file and restart Calibration Control.

Troubleshooting Label Printing

Label Print Failed

If you receive an error message window upon printing a label, you'll notice there are several error messages. The first will be the *most likely* cause but might not always be the most accurate. The **Code Error** message will show a very short description of the error. Finally, the **System Error** will show the most accurate description of the error, though it may require some research.

1. Take a screenshot of the error message to refer back to.
2. Eliminate the most likely cause by reinstalling the label printer drivers and/or disconnecting/reconnecting your label printer. Search for your specific label printer on Brother's support site to find the proper troubleshooting methods for your device.
3. Take a look at the Label Print Failed error message below and the different solutions to find one that works for you.



Solution 1

Some users have encountered this error after updating from Windows 7 to Windows 10. The most successful solution is downloading and installing these Brother label printer drivers: [b-PAC3CCI Setup](#).

Solution 2

You may see a slight difference in the error message if you are using a Windows 7 computer and you may see there is an error with **vcruntime140.dll**. Downloading and installing the newest version of Microsoft Visual C++ Redistribution package has solved this error as well. You can find the latest supported Visual C++ downloads here: <https://support.microsoft.com/en-us/help/2977003/the-latest-supported-visual-c-downloads>


Get Ape Support

If you've tried all the solutions provided and you are still unable to print a label, email us and include the error message screenshot. We may have encountered the same error before and may have a quick solution for you or we can troubleshoot together with a call/ screenshare.

No Label File Error

Label File Could Not Be Loaded

When attempting to print a label in Calibration Control, the following dialog will display when **(a)** there is a problem with the label drivers, **(b)** printer, or **(c)** the label file is missing.



This situation most often occurs after a change in computer settings, such as **(a)** new computers, **(b)** wiping hard drives and reinstalling the OS, **(c)** cleaning after a virus, and **(d)** installing or uninstalling other software.

The main issue is that the [printer drivers](#) are not reinstalled or not installed properly, as was previously accomplished.

Label Printer Drivers

Download and install the Brother label printer drivers for your model. When using a label printer connected to the network (rather than a computer) or a label printer is connected to another computer, you will need to manually add the label printer to the Windows Devices and Printers.

P-touch Editor Software

If you can see your label printer enabled in the Windows Devices and Printers but are still unable to print labels, try installing the free Brother P-touch Editor Software.

If this still does not resolve the problem, contact Ape Software for assistance.

Troubleshoot SQL Server Connection

For Use with Ape Software Databases

WARNING: Ape Software applications can only connect to databases created by and specifically for the given Ape Software application. Blank databases, databases created by non-Ape professionals, and modified Ape Software databases will not connect or not connect properly. Also, SQL Server databases created by migrating an existing MSAccess database will also not work and often have severe functional issues sometimes not immediately visible.


Universal Data Link (UDL) File

Use a UDL file to create a direct connection to a SQL Server database with as few other variables as possible, like connecting through other software. This helps to simplify and focus the troubleshooting process on the most likely failurepoints (e.g., SQL Server configuration or network permissions).

A UDL file is simply a text file with a .udl extension. Therefore, creating a UDL file is as simple as creating a new text file and changing its extension to UDL.

Create a UDL File

Create a UDL file by **(1)** right-clicking in a folder or desktop, **(2)** selecting 'New', **(3)** selecting 'Text Document', and then **(4)** creating a text file with a .udl extension (e.g., MyTestFile.udl).




Note: If a warning appears related to making the file unusable by changing its extensions, ignore the warning.

Data Link Properties

Double-click on the new UDL file and the following Data Link Properties dialog appears.

Provider Tab


Select the appropriate SQL Server Provider.




Connection

In the Connection tab . . .

1. Enter the fully qualified SQL Server name.
2. Select the authentication method (i.e., Windows or {1st choice} or SQL Server {2nd choice} and enter the appropriate user name and password, if SQL authentication is used.
3. Select the database. By default, this is apecal but must match whatever name is actually used on the SQL Server.
4. Click the [Test Connection] button.



If the above test fails, then work with your organization's Database Administrator (DBA) to resolve the connection before attempting to use the Ape Software application. If using SQL Server authentication, the **User Security** section of the AddApe Database to SQL Server help topic is often useful.



Otherwise, if the above process leads to a 'Test connection succeeded' dialog and the Ape Software application cannot connect to the SQL Server, contact Ape Software to continue troubleshooting the connection.

SMTP Test Procedure


Introduction

If you are having difficulty getting the *Test Email Settings* to work, this test procedure will determine whether or not you have access to your SMTP server and help to identify the problem.

Network firewalls or antivirus software sometimes block email traffic preventing you from accessing the mail server, which is required in order for the Auto Notification utility to function.

This procedure should only be used if you are getting a system error or timeout message.

If you are getting an *authentication error*, the problem is probably due to an invalid username/password or unknown source.



The screenshot shows the "Email Options" dialog box with the following settings:

- Email Server Settings:**
 - SMTP Server: smtp.gmail.com
 - Email: service@apesoftware.com
 - Password: [Redacted]
 - Seconds Between Emails: 0
 - Port (0 = Auto): 0
 - Secure Socket Option: None
 - Disable CSS in Emails:
 - Connect Asynchronously:
 - Use Secure Connection:
- Output Folder:** [Empty field]
- Email File Mask:** email_yyyyMMdd_(RND32).eml
- Enabled:
- Use Legacy Email (pre 9.2)
- Email Settings Test:**
 - Send Test Email to: [Empty field]
 - Or Send to Email Output Folder
 - Send Test Email button


Buttons: Cancel, Save

Procedure

Verify Domain:

To begin, open a Command Prompt Window. Click on the Start menu from your computer's desktop and type "CMD" in the run or search box. In Windows 10, right click on the Window (start) icon and click "Command Prompt".


Then type in “*nslookup*” and insert your SMTP server address, like the example below for Gmail:



```
Command Prompt
Microsoft Windows [Version 10.0.14393]
(c) 2016 Microsoft Corporation. All rights reserved.

C:\Users\Ape Software>nslookup smtp.gmail.com
```

After pressing "Enter" you should see a result similar to the one below. This will verify that your SMTP server address is valid.



```
Command Prompt
Microsoft Windows [Version 10.0.14393]
(c) 2016 Microsoft Corporation. All rights reserved.

C:\Users\Ape Software>nslookup smtp.gmail.com
Server: cdns01.comcast.net
Address: 2001:558:feed::1


Non-authoritative answer:
Name: gmail-smtp-msa.l.google.com
Addresses: 2607:f8b0:400e:c06::6d
           173.194.203.109
           173.194.203.108
Aliases: smtp.gmail.com

C:\Users\Ape Software>
```

Adding TELNET:

Once you have verified the domain, we can try to connect to the SMTP server. This requires the use of TELNET. You will probably need to add TELNET since this function is not available by default.

To add this feature, open up your Control Panel and select “**Programs and Features**”. In Windows 10, click on the Start button then click on the gear icon to get to settings. In the search bar type "turn Windows features on or off". Within the Windows Features dialog, check the box next to "Telnet Client" and click OK to save the change.




Once Windows has installed this feature, click "Close" to complete this step.

Connect to SMTP server:

Now go back to the Command Prompt window and type the following:

```
telnet "SMTP server address" "port"
```



If you are able to connect, you should see a result similar to the one below:


```
Telnet smtp.gmail.com
220 smtp.gmail.com ESMTP c185sm4499581pfb.57 - gsmtp
```

If the result says "Connect Failed", you do not have a connection to your SMTP server using that port. Reasons for this include:

- The email server is not functioning
- The traffic is being blocked by network security

UDL File for Database Connection

Universal Data Link for Database Connection


Use a Universal Data Link (UDL) file when the Ape software Program Files are installed on a central file server and the central database is a SQL Server. The UDL file is then used to define the SQL Server and database name. Finally, the UDL file is consumed by creating a program shortcut with a switch pointing to the UDL file.

Without using a UDL file, an Ape setting file on each of the Client computers (i.e., general.config) requires a property that defines the location of the database server and database. This means the database location is defined, and can be changed, in multiple locations.


Create Universal Data Link (UDL) File

Follow these steps to create a UDL file:


1. Open Windows Explorer to the location where the UDL will be created.
2. Right-click and select New -> Text Document, naming the file something like "apedb.txt".



3. Change the extension of the new text file to ".udl" with an end result of "apedb.udl".
4. Double-click on the new file to open the Data Link Properties dialog. Select the Provider tab and select the "Microsoft OLE DB Provider for SQL Server".



5. Select the Connection tab and enter the (1) Server Name where the Ape Calibration database is located, (2) selectWindows NT Security, and (3) select the database name.




6. Click the [OK] button to save the UDL file settings.


Create Shortcut

Now create a new Shortcut to the "apecal.exe" executable with a switch at the end in the format of -udl:"[path statement ofudl created in previous steps]" and include the quotations. Ensure both path statements are valid for all Client users.

1. While still in Windows Explorer, right-click and select New -> Shortcut.



2. In the Create Shortcut dialog, add the "[executable path]" -udl:"[udl path]" statement like below.



3. Click the [Next] button, give the Shortcut a name like "Calibration Control" and then click the [Finish] button.

Copy Shortcut to Clients

Remember that even though the Program Files for Calibration Control are installed on a File Server, the program itself actually runs on the Client computers. Therefore, each Client computer must have the Minimum Requirements.

1. Ensure the Minimum Requirements are installed on the Client computer.
2. Copy the newly created Shortcut to a Client computer.
3. Test the new shortcut.

If you need help making this feature work, please Contact Us for help.

Startup INI

Single Connection Configuration for All Users

In network environments where the Calibration Control (our **Calibration Management Software**) program files are installed on a single file server and executed from multiple clients, the Startup INI (startup.ini) gives Network Administrators the ability to configure database connection properties in a single location. Each time Ape starts, it looks for the "startup.ini" file in the same folder as the executable. If found, the connection settings in the "startup.ini" file take precedence over connection settings, or absence of settings, at the Client level.


Additionally, all data stored in the INI file (e.g., database name, location, server, password, etc.) is 128 bit encrypted.

Before Using the INI Utility

Before creating the “startup.ini”, install the Ape Software on the computer/server and Create a Multi-User Environment. Note that the Minimum System and Configuration Requirements is a highly relevant and useful document that will help prevent deployment issues and should be followed. In addition, the Startup INI feature may be hidden from your toolbar. To learn how to show this feature, follow the Feature Visibility help topic.


Start the INI Utility

Select the Startup INI menu option of the Utilities tab of the ribbon menu.



Startup INI Editor

Work with your internal IT group to assist with setup if further help is needed. All enabled fields (fields that can be typed in) are required except for the Password fields when the Database Type is set to [MS] Access.

A screenshot of the 'Startup INI Editor' dialog box. The dialog has a title bar with a question mark and a close button. It contains two sections: 'Database Type' and 'Authentication Method'. Under 'Database Type', the 'SQL' checkbox is checked, and the 'Access' checkbox is unchecked. Under 'Authentication Method', the 'SQL' checkbox is checked, and the 'Windows' checkbox is unchecked. Below these are five text input fields labeled 'Server:', 'Database:', 'User:', 'Password:', and 'Confirm:'. There is also a 'Password Visible' checkbox which is currently unchecked. At the bottom, there are three buttons: 'Cancel', 'Delete startup.ini', and 'OK'. Below the buttons, there is instructional text: 'Create/update the startup.ini file in PROGRAM FOLDER.', 'All users who start the EXE in this PROGRAM FOLDER (C:\Program Files (x86)\Ape Software\Calibration Control) will automatically use the setup.ini file if it exists.', and 'You must have read/write privileges to the PROGRAM FOLDER to use this utility.'

Updating Access DB

This help topic covers possible difficulties related to updating an MS Access DB for Calibration Control (CC). Refer to the Update a SQL Server Database help topic for SQL Servers.

Problem Description

Revisions of CC often require / make changes to the structure of the DB. These changes include tables, index, fields, and query creation, deletion, or modification. The speed and ability of CC to update an Access DB is affected by the speed and stability of a network in addition to the size and complexity of the update itself. Therefore, the combination of a slow or unstable network with a large update (e.g., version 9.0) can result in a failure to update. The following are methods to turn such a failure into a success.

Automatic Backups

Before the update process begins, CC automatically creates a backup of the Access DB and places that backup file in the Backup folder below the Database Folder (refer to the Folders help topic). The name format of the backup files is 'apecal-auto-backup-[large number].dbk'. The backup files are renamed copies of the Access DB. This means that the backup files are still valid MS Access DB files, just with a different extension (*.dbk).

It's important to understand the backup process (above) since an update activity may fail / freeze at any point. Therefore, the backup process may need to be re-started. If a re-start is required, it should start with a database file that has not been altered by a previous failed update attempt.

Local Update

The simplest method to bypass a slow or unstable network (e.g., noisy/weak Wi-Fi) is to move/copy CC's Access DB (usually named apecal.mdb) to the local HDD (Hard Disk Drive) of a computer with the new version of CC installed. Update the DB locally and then move the DB back to its network location. Here are a few explicit steps:

Step 1 – Move DB to Local HDD: Move the CC DB to the local HDD (refer to the Locating Database and Move Database help topics). If previous update attempts already failed, use the first backed up copy (*.dbk) of the database (before an update attempt was made) instead of the main database file (*.mdb). Do this by copying and renaming the backed up file to the correct file name (e.g., apecal.mdb) and placing that copy on the local computer's HDD.

Step 2 – Delete Config File: Delete the general.config file (located at C:\Users\Public\Documents\Ape Software\Calibration Control). When the general.config file is deleted, CC provides a dialog that helps you navigate to the copied DB on the local HDD (see step 1, above).

Step 3 – Update DB: Start CC and allow the software to update the local DB.

Step 4 – Return Updated DB: After the update succeeds, close CC and place the DB back in its original network location. Ensure the local DB is moved, deleted, or re-named so you don't accidentally connect to it when re-starting CC.

Step 5 – Connect to Updated DB: Delete the general.config file again and re-start CC, this time navigating to the original network location.

If everything worked as it should, you will be up-and-running. If not, Ape Software can provide a free update service.

Current Access Database Drivers

If you are repeatedly receiving a message saying that you need to set ACE to 'True' in your general.config, then you probably need the download and install database drivers. You can find these drivers on the Minimum Requirements page or just download from here.

Ape Software Free Update Service

If you need Ape to update your DB for you, simply send us an un-altered zipped copy.

Un-Altered DB Copy: Ensure you send us a copy of your DB that has not previously failed the update process. Refer to the Automatic Backups section of this help topic for more info.

Send Method: If your zipped DB is less than 5 MB, you can email a copy to service@apesoftware.com. Otherwise, you can upload your zipped file to a cloud service you are familiar with or we can send you an upload link.

If you have any specific questions, please contact us.

Update a SQL Server Database

Update an Ape Software SQL Server Database

Use these instructions to update the database structure of a Calibration Control database installed on a SQL Server 2012 or higher.


Before Updating

Open **SQL Server Management Studio** and sign in to a SQL Server as a SYSTEM ADMIN because you'll require enough permissions to create and update databases and their objects. Otherwise, your IT Group may need to perform the update.

All of the SQL files needed to update the database reside in the SQL Tools folder under the CC program folder located in your computer at C:\Program Files\Ape Software\Calibration Control.


Before attempting to update an apecal database **protect your database with a backup**. Right-click the database, hover over Tasks, and then select Back Up.

Follow along with our Backup/Restore SQL Server Database help topic for more detailed instructions.



Updating Database

1. Determine the current database structure version by first expanding the database tree and expand the Tables folder. Then, right-click the `dbo.tbIDBVersion_DO_NOT_EDIT` table and Select Top 1000 rows. There should only be one record.



If you do not have a StructureVersion field in the file name (versions 7.2.5 and older), use the following table as a guide to interpret the dbo.tblDBVersion table and how to choose the appropriate sql_update SQL file to begin with:

□ DB Version Interpretation

Indication	Database Version
dbo.tblDBVersion_DO_NOT_EDIT table does not exist	1
Rev = 1.1, 7.0, 7.0.1, etc.	Rev = Version
Rev = DoNotEdit	Use StructureVersion field
StructureVersion = 66	66

2. Open each SQL Update file with SQL Server Management Studio. Once again, ensure you backup your database before proceeding.

Note: If your database name is not the standard "apecal", change the database name before proceeding. The database name may appear in the following lines of the update files:

```
ALTER DATABASE apecal -- ***
SET OFFLINE WITH ROLLBACK IMMEDIATE
ALTER DATABASE apecal --
***SET ONLINE
USE apecal --
***GO
ALTER DATABASE apecal SET PAGE_VERIFY CHECKSUM --
***GO
```

3. Apply the sql_update SQL files in their sort order beginning with the database version indicated by the Db Version table(above).

For example, if the StructureVersion field of your database equals 66, then begin with file **sql_update_22 (db ver 66 to 73).sql** then click the [Execute] button to apply the changes.

Note: Ensure that none of the text in the SQL file is highlighted when the [Execute] button is pressed because this will run only the highlighted text.

4. Once the query is complete (without any red query errors), continue to the next file until you finish with the final SQL file.

5. If you encounter any query errors, take a screenshot of the error(s) and contact us for assistance.

List of Updates

SQL Update Script File

sql_update_01 (pgm ver 1 to 1.1).sql (read as version 1.0 to 1.1)

sql_update_02 (pgm ver 1.1 to 7.0).sql (same as above)

sql_update_03 (pgm ver 7.0 to 7.0.1).sql

sql_update_04 (pgm ver 7.0.1 to 7.2.0).sql

sql_update_05 (pgm ver 7.2.0 to 7.2.2).sql

sql_update_06 (pgm ver 7.2.2 to 7.2.5).sql

sql_update_07 (pgm ver 7.2.5 to 7.2.5.8).sql

sql_update_08 (pgm ver 7.2.5.8 to db ver 15).sql (read as *program* version 7.2.5.8 to database version 15 {new system})

sql_update_09 (db ver 15 to 24).sql (read as database version 15 to 24)

sql_update_10 (db ver 24 to 30).sql

sql_update_11 (db ver 30 to 31).sql

SQL Update Script File

sql_update_12 (db ver 31 to 38).sql

sql_update_13 (db ver 38 to 39).sql

sql_update_14 (db ver 39 to 40).sql

sql_update_15 (db ver 40 to 41).sql

sql_update_16 (db ver 41 to 42).sql

sql_update_17 (db ver 42 to 43).sql

sql_update_18 (db ver 43 to 46).sql

sql_update_19 (db ver 46 to 51).sql

sql_update_20 (db ver 51 to 64).sql

sql_update_21 (db ver 64 to 66).sql

sql_update_22 (db ver 66 to 73).sql

sql_update_23 (db ver 73 to 81).sql

sql_update_24 (db ver 81 to 84).sql

Get Ape Support


Contact us and we can setup a screen share and update your database together or if you prefer, you can send us a backup of your database and we can update your database for you.

Backup and Restore SQL Server Database


Create a Full Database Backup

While there are many ways to back up your database, a full database backup backs up the entire database (including part of the transaction log) which makes it simple to restore. However, as the size of the database increases, full database backups will take longer to complete and require more storage space. For larger databases, consider using a series of differential backups in addition to a full database backup.


1. Using SQL Server Management Studio, log in to the appropriate server with sufficient permissions and, if needed, expand the server tree.
2. Expand the **Databases** folder.
3. Right-click the database that you wish to backup, hover over Tasks, and then click **Back Up...**



4. The Back Up Database dialog will appear and you will see your selected database in the dropdown list. You can select a different database from the dropdown if you choose to do so.




5. By default, the Backup type will be **Full**. You must create a full backup before you can select a differential backup type.
6. Choose **Database** under the Backup component (if not already selected).
7. Under the Destination section, your database will back up to a default SQL Server Backup location, however you can Add/Remove additional location(s) or remove the default file path in order to Add and rename the backup file name (still using the default backup location). Ensure you have read/write permissions to the selected file path.
8. Note: If you have previously created a backup of the same database, go to the **Media Options** page and review the settings. Choose whether this backup will append to the existing database or overwrite it.
9. Recommendation: Go to the **Backup Options** page and review the settings. Under the Compression section, compress the backup to save storage space (database compression is not available in SQL Express).
10. Click OK to start the backup process. Once the backup is successful, a dialog will appear to notify you.




Restore a SQL Server Database Backup

Follow these instructions for restoring a full database backup. Once again, using SQL Server Management Studio, connect to the appropriate server and expand the server tree.

1. Right click the **Databases** folder and select **Restore Backup**
2. Select **Device** as the Source.




3. Click the browse button
4. Then, click on the Add button. The dialog will automatically show available backups in the default backup location. If you previously saved your backup in a different location, copy the path and paste it in the Backup File Location box. Note: Only .BAK file types will be available in this dialog.
5. Select your backup file and click OK.




6. The database name may be automatically loaded if you previously restored this backup. You can change the Database name if needed.

7. At the point, you may see "Ready" at the top of this dialog. This means you can restore your database now or continue and adjust additional settings.



8. The restore process may fail if there are existing connections. To prevent this, go to the **Options** page and check "Close existing connections to destination database" located under the Server connections section.

- Optionally, if you want to overwrite an existing database, select "Overwrite the Existing database". The databasename will need to match the existing database.
- Click OK to start the restore process. Once the restore is successful, a dialog will appear to notify you.




Downgrading a Database

Return a Database to the Previous Version

Sometimes after a new version of Calibration Control is deployed a decision is made to return to a previous version. Unfortunately, many Ape updates include changes to the database which cannot be uninstalled as easily as the software. This is why a manual database backup is strongly encouraged before updating the MS SQL Server database.

Restore an Access Database

- Refer to the Folders used by Ape help topic before continuing.
- Open the Data folder and rename the apecal.mdb file (e.g., higher-version-apecal.mdb).
- In the Backup folder under the Data folder, rename the most recent database backup file from apecal-auto-backup-[long number].dbk to apecal.mdb.
- Install the previous version of Ape or download and install one of the older Ape versions that your product key works with.
- Restart the previous version of Ape to confirm that the database works.



Restore a SQL Server Database

Use the same instructions in the Installing SQL Server help topic – Restore (1st Method), but choose the file that was backed up as part of the update process.